ODYSSEY[®] iQ

microprocessor foot

college park

TECHNOLOGY *for the* HUMAN RACE

DOMESTIC
user guide

benutzerhandbuch • guía del usuario
guide d'utilisation • guida utente

FIGURE 3


FIGURE 4

FIGURE 5


FIGURE 6


FIGURE 7



FIGURE 8

FIGURE 9

FIGURE 10


FIGURE 11

FIGURE 12

LANGUAGES

EN-English	6
DE-Deutsch	25
ES-Español	44
FR-Français	63
IT-Italiano	82

Federal Law restricts that the Odyssey iQ Microprocessor Foot is to be purchased, configured, and fit only by a board-certified prosthetist, licensed by the state in which they practice.

The operations contained in this document are intended for operator use.

INTENDED USE STATEMENT

The Odyssey iQ is a prosthetic foot designed to replace one or more functions of the biologic human foot.

What's Included With Your Device?

Hardware

Odyssey iQ Microprocessor Foot Battery Charger w/ Adapter

These diagrams are to help familiarize you with the unique parts of the Odyssey iQ foot. These parts are referenced in the instructions and used when speaking with a technical service representative.

KEY COMPONENTS (FIGURE 1)

- | | | |
|------------------------------|--|--|
| A. Integrated Pyramid | B. Ankle-Housing Assembly | C. Toe Spring |
| D. Heel Spring | E. Permanent Wedge | F. Rounded Heel Wedge (Optional) |
| G. Battery | H. Power Button / LED Indicator | I. Charging Receptacle / Dust Cover |
| J. Serial Number | • CPI Sock (not shown) | • Foot Shell (not shown) |

TECHNICAL SPECIFICATIONS

STRUCTURAL SPECIFICATIONS

Foot Size	Weight Limit	Build Height	Foot Weight*
23-25 cm	275 lbs / 125 kg	3.5 in / 8.8 cm	1,027 g
26-30 cm		3.6 in / 9.0 cm	

*26cm foot w/shell

GAIT MATCHING® GUIDELINES

The gait match determines the firmness of the foot based on the user's specifications (foot size, patient weight, and impact level).

FIRMNESS CATEGORIES

Refer to the chart below to determine the correct firmness category.

Note: Incorrect category selection may result in poor device function. Contact College Park Technical Service if you have questions about category selection.

WEIGHT LBS	0-140	141-180	181-220	221-275
WEIGHT KG	0-63	64-81	82-100	101-125
SIZE CM	23-30			
LOW IMPACT	1	2	3	4
MODERATE IMPACT	2	3	4	5

DEVICE OPERATION

Battery	Lithium polymer, 3.7V, 980 mAh	Time to full charge	4 hours (time to 50% charge = 2 hours)
Charging Cable	5V, 3A DC (Use only provided power supply to charge the battery)		

WIRELESS

Connection	Bluetooth 5.0/Bluetooth Low Energy (BLE)	Maximum Speed	24 Mbps
Operating Frequency	2.402 - 2.480 GHz	Effective Radiated Power	12.00 dBm

ENVIRONMENTAL USE CONDITIONS

Charging	32F to 104F (0C to +40C)
Operating	14F to 122F (-10C to +50C), Pressure: 70-106 kPa, Humidity: 0%-100% relative humidity
Storage & Transport	-4F to 140F (-20C to +60C) *, Pressure: 70-106 kPa, Humidity: 0%-90% relative humidity

Note: If storing device above or below operating temperature, allow the device to return to within operating temperature range before use. The device must be brought up to the operating temperature, allow the device to sit for 15 minutes.

IP RATING

IP67 | Dust tight and protected against freshwater at a maximum depth of 1 meter for 30 minutes.

THE BATTERY

The Odyssey iQ foot contains an internal lithium polymer battery. This battery supplies 3.7 V 980 mAh. It is advised to charge the battery daily. For most users, the battery will last up to 16 hours of continuous use, depending on the condition of the battery and the frequency of use.

The battery is shipped with a partial charge (up to 30%). We recommend charging the battery to 100% before use.

LED INDICATOR

The Odyssey iQ foot is equipped with a state-of-charge indicator. This indicator tells the user how much life is left in their current battery. The Odyssey iQ foot must be powered on to use this feature.

Press the power button for 1-second to activate the LED indicator light (*Figure 1-H*). The number of blinks indicates the state of the charge. Charge level can also be checked via the Stride Studio app.

FEATURE

POWER BUTTON

Description

Press and hold for 1 second to indicate the battery status.

FOOT LED INDICATIONS – BATTERY

COLOR	INDICATOR	BATTERY STATUS
Solid Green		Fully charged, 91-100%
4 flashes		75 - 90%
3 flashes		50 - 74%
3 flashes		25 - 49%
1 flash		1 - 24%

CRITICALLY LOW BATTERY STATUS

The foot will go into safety mode when in a critically low state. Safety mode means that the foot goes to the set plantarflexion and dorsiflexion resistances and will no longer compensate for terrain or speed changes and won't be able to enter Stance iQ mode.

BATTERY CHARGER

Odyssey iQ feet are supplied with a charger for the lithium polymer battery. The charger is recommended for daily use and will ensure that the battery will receive a full charge and provide maximum running time. There are four power adapter options (US, UK, EU, or AU) to match the needs of different regions. Use only the provided power supply to charge the lithium polymer battery.

CHARGING THE BATTERY

CHARGING THE ODYSSEY IQ FOOT BATTERY

1. Only charge the foot on a flat, level surface and in a cool, dry location away from sources of heat and moisture. Plug the wall adapter into the main AC supply (Figure 2-B).
2. To connect to the front of the ankle, open the dust cover and insert the connector..

Note: The connector is directionally specific (Figure 2-A).

3. The foot will audibly indicate when it begins charging or is unplugged.

Note: Ensure charger port is free of moisture and debris before charging..

Recharge time from empty is approximately 4 hours.

POWER BUTTON

The power button is located on the top of the ankle, near the anterior opening of the foot shell (Figure 1-H). To turn the device on or off, press and hold the power button for 3 seconds. When the device is powered on or off, a green light on the LED indicator will display for about 3 seconds.

FEATURE

POWER BUTTON

Description

Turns foot on/off. Initiates Bluetooth pairing when prompted in app.

POWER BUTTON

ON/OFF

Press and hold the power button for 3 seconds.

APP OVERVIEW

STRIDE STUDIO APP FOR USERS

Users may choose to download the Stride Studio app, but it is not required.

DOWNLOAD STRIDE STUDIO APP

College Park's Stride Studio App is available for both iOS & Android devices. Download Stride Studio from Apple App Store or Google Play Store.

MINIMUM OPERATING SYSTEMS:

- Apple iOS 15
- Android Version 8 (Oreo)

Select Stride Studio icon on the device to launch the app.

LOGGING IN

The landing screen allows you to select whether you are the user or prosthetist.

*Users: There is no login required, Select **"User"** (Figure 3).*

PAIR DEVICE SCREEN

ENABLE BLUETOOTH

Go to Settings on your phone or tablet, then tap Bluetooth.

Launch the Stride Studio app and select "User". The app will prompt you to allow it access to Bluetooth. You must grant it access. You will land on the "Devices" screen.

BLUETOOTH PAIRING

The foot must be paired to the Stride Studio app when it is used for the first time. Select “Connect New Device” (Figure 4), and a list of available devices will be displayed. Select the serial number of the foot from the menu (Figure 5). A Bluetooth pairing request will pop up. Press the “Pair” button (Figure 6). The app will display an “Unauthorized” pop up. Press the power button on the front of the ankle for 1 second to complete the pairing process (Figure 7).

Note: The foot will be identified by the Serial Number, which is located on the left side of the ankle.

Note: Phone or tablet must be closer than 33 feet (10 meters) of the Odyssey iQ foot.

Note: Only one device may be paired to the foot at a time. Prosthetists must disconnect from the foot before the user can pair in the app. Disconnect the foot from the “Devices” screen or by closing the app.

For Bilateral Users: To program easy switching between left and right feet, follow the Linking process in the Setup Wizard.

USER SELECTED MODES

STANCE IQ MODE

Stance iQ mode is a function where the foot detects that the user is standing still and enters high resistance.

If Stance iQ mode is enabled, the foot will lock plantar and dorsiflexion for added support, after 8 seconds of standing still. Once the user begins walking again, it will revert to normal function allowing for plantar and dorsiflexion.

WALKING MODE

In this mode, Stance iQ is disabled and the foot will not enter high resistance.

ANKLE LOCK MODE

Enters the foot into high resistance until deactivated.

CUSTOM SHOE PROFILES

Allows prosthetist and user to program additional shoe profiles. Calibrates the angle for shoes with heel heights between 0-0.80 in (0-2.0 cm).

BILATERAL FOOT LINKING

Set up by the prosthetist, so that prosthetists and bilateral users can easily switch connection to the left or right foot during app use.

USER DASHBOARD

The user dashboard screen is the main hub of the Odyssey iQ foot (*Figures 8-A, 8-B*). It provides access to basic foot information, mode selection, and critical information.

- 1 About - Links to User Guide and social media resources.
- 2 Lifetime Step Counter - Total step count for lifetime of device.
- 3 Mode Selector - See Mode Selector section.
- 4 Reset - Resets the step counter to zero.
- 5 Step Counter - Tracks number of steps between resets.
- 6 Switch Foot - See Switching Between Left and Right Feet section.
- 7 Custom Shoe Profiles - See Custom Shoe Profiles section.
- 8 Battery Status - Indicates battery charge level of Odyssey iQ foot.
- 9 Alerts Settings - Turn non-critical alert sounds on/off.

MODE SELECTOR

The slider bar allows you to toggle between walking mode **(1)** and Stance iQ mode **(2)**, or from Stance iQ mode **(2)** to ankle lock **(3)** (*Figure 9*). An alert will appear each time the slider is moved to Stance iQ mode **(4)** (*Figure 10*).

SWITCHING BETWEEN LEFT AND RIGHT FEET

The foot linking process can be set up so that bilateral users can easily switch between their left and right feet during app use. This icon will only appear if linking has been enabled by your prosthetist (*Figure 8*).

Tap the “Switch Foot” icon to switch between left & right feet.

CUSTOM SHOE PROFILES

This section calibrates the angle for shoes with heel heights between 0-0.80 in (0-2.0 cm). Assign a new shoe profile to optimize variance in heel heights or heel densities (very soft or very hard-heeled shoes). (Figure 11).

To change shoe profiles, begin in a standing position, wearing your preferred shoes. When a profile is selected, it will take effect immediately and remain active until a different profile is selected.

- 1 **Existing Mode** - Programmed modes are displayed here. A checked box indicates the current mode.
- 2 **Add a Shoe** - Click to add a new shoe profile.
- 3 **Edit** - Opens Edit screen.

ADD CUSTOM SHOE PROFILES

To add a new shoe profile, begin in a standing position, wearing your preferred shoes.

Select “Add a Shoe” from the “Shoes” screen (Figure 11). Enter a name for the new shoe, select “Calibrate” and follow the calibration process (Figure 12). When complete, tap the “Save” icon to apply the new shoe profile.

- 1 **Name** - Enter a name for the new shoe profile.
- 2 **Calibrate** - This section calibrates the angle for shoes with heel heights between 0-0.80 in (0-2.0 cm).
Note: Bilateral shoes must be calibrated individually.
- 3 **Save** - Applies settings to shoe profile button.

TROUBLESHOOTING

CAUTION

Caution: The Odyssey iQ foot should never be serviced while connected to the end user. Ensure that the device is disconnected and powered off before any service or maintenance is performed. This device should never be serviced while in use. Never let children use this handle this device unsupervised. Take caution when using this device around pets that may cause damage to the device.

The Odyssey iQ foot features an LED indicator located on the top of the ankle, near the anterior opening of the foot shell. This LED is used to display battery related information.

If the Odyssey iQ is not functioning as expected, try the following:

- Turn the system off, wait several seconds and power back on.
- Make sure the battery has sufficient charge. If the battery is too low, charge the battery.

WARRANTY INSPECTION AND MAINTENANCE INFORMATION

CLEANING INSTRUCTIONS

Clean the Odyssey iQ foot by wiping it with a soft, lint-free cloth. Never apply harsh chemicals.

USE IN WATER

Clean the Odyssey iQ foot by wiping it with a soft, lint-free cloth. Never apply harsh chemicals.

- *The Odyssey iQ foot has been approved for use in freshwater at a maximum depth of 1 meter for 30 minutes.*
- *Use caution when walking on wet surfaces.*
- *The foot should be switched into ankle lock mode when using it in water or near wet surfaces.*
- *Always dry foot thoroughly after exposure to moisture. Rinse and dry if exposed to salt water or chlorinated water.*

COSMETIC FINISHING

College Park does not recommend cosmetically finishing the prosthesis if it will be exposed to moisture.

WARRANTY INSPECTION AND MAINTENENCE INFORMATION

College Park recommends that you schedule your patients for check-ups per the warranty inspection schedule below.

High patient weight or activity level may require more frequent inspections. We recommend you visually inspect the following applicable parts for excessive wear and fatigue at each warranty inspection.

- | | |
|----------------------------------|---------------------|
| ▪ <i>Hydraulic Assembly</i> | ▪ <i>CPI Sock</i> |
| ▪ <i>Composites and Adapters</i> | ▪ <i>Foot Shell</i> |
| ▪ <i>Battery</i> | ▪ <i>Wedge</i> |

WARRANTY INSPECTION SCHEDULE FOR THE ODYSSEY IQ: SIX MONTHS, THEN ANNUALLY.

TECHNICAL ASSISTANCE / EMERGENCY SERVICE 24-7-365

College Park's regular office hours are Monday through Friday, 8:30 am – 5:30 pm (EST). After hours, an emergency Technical Service number is available to contact a College Park representative

LIABILITY

The manufacturer is not liable for damage caused by component combinations that were not authorized by the manufacturer.

CAUTION

College Park products and components are designed and tested according to the applicable official standards or an in-house defined standard when no official standard applies. Compatibility and compliance with these standards are achieved only when College Park products are used with other recommended College Park components. This product has been designed and tested based on single patient usage. This device should NOT be used by multiple patients.

CAUTION

If any problems occur with the use of this product, immediately contact your medical professional. The prosthetist and/or patient should report any serious incident* that has occurred in relation to the device to College Park Industries, Inc., and the competent authority of the Member State in which the prosthetist and/or patient is established.

*“Serious incident” is defined as any incident that directly or indirectly led, may have led, or might lead to any of the following: (a) the death of a patient, user, or other person, (b) the temporary or permanent serious deterioration of a patient’s, user’s, or other person’s state of health, (c) a serious public health threat.

COMPLIANCE

This device has been tested according to ISO 10328 standard to two million load cycles. The expected life duration under normal usage is 2 million steps. Depending on the amputee’s activity this corresponds to a duration of use of two to three years.

ISO 10328 - LABEL

WEIGHT LIMIT (KG)	LABEL TEXT
125	ISO 10328 -P6-125 kg

ISO 10328 - “P” - “m”kg*)

**) Body mass limit not to be exceeded!*

For specific conditions and limitations of use, see intended use section of manufacturer’s written instructions.

SAFETY

Only start up the product in accordance with the information contained in the supplied documents. Failure to follow these technical instructions or use of this product outside the scope of its Limited Warranty may result in injury to the patient or damage to the product.

LEGEND OF SYMBOLS

Note: Possible technical damage.

Info: Basic information regarding this product.

Caution: Possible risk of accident or injury.

Warning: Possible risk of severe accident or injury.

SAFETY INSTRUCTIONS

Info: Use on airplanes

Airlines may not permit the use of this device on their aircraft. Check with the airline before travelling to ensure this device is allowed for use on the plane.

Info: Disposal

These products may not be disposed of with household waste in some jurisdictions. Disposal that is not in accordance with the regulations of your country may have a detrimental impact on health and the environment. Please observe the information provided by the responsible authorities in your country regarding return and collection processes.

Caution: Risk of pinching where the ankle joint bends

Ensure that fingers and other body parts are not in this area when bending the ankle joint.

Caution: Unsupervised use

It is not recommended for children to operate this device without the supervision of an adult. Use extreme caution around small children and household pets.

Caution: Risk of accident while operating a vehicle

A lower extremity amputee's ability to drive a vehicle is determined on a case-by-case basis. Factors include the type of fitting (amputation level, unilateral or bilateral, residual limb conditions, design of the prosthesis) and the amputee's abilities. All persons are required to observe their country's national and state driving laws when operating vehicles. For insurance purposes, drivers should have their driving ability examined and approved by an authorized test center.

Caution: Improper use

Any type of excessive strain, overload or improper use may lead to faulty control or malfunction of the Odyssey iQ foot, resulting in a risk of injury. The Odyssey iQ foot was developed for everyday use and must not be used for unusual activities. These unusual activities include, for example, sports with excessive strain and/or shocks, or extreme sports. Careful handling of the prosthesis and its components not only increases their service life but, above all, ensures your personal safety! Should the prosthesis be subjected to unusual stresses (such as a fall), immediately contact a certified prosthetist and have the prosthesis inspected for any damage.

Caution: Consequences of product deterioration

Wear and tear on system components can lead to malfunction of the Odyssey iQ foot, resulting in a risk of injury. Follow the specified service intervals. The service life of this device is 3 years. Patient concerns about the function should be reported to the prosthetist immediately, including but not limited to noise, sudden loss of function, etc.

Caution: Mechanical overloading

External mechanical influences or loads, such as impacts and vibration, can lead to faulty control or malfunction of the Odyssey iQ foot and result in a risk of injury. The Odyssey iQ foot should not be subjected to mechanical vibrations or impacts.

Caution: Manipulation of system components

Independent changes and/or modifications to system components may lead to faulty control or malfunction of the Odyssey iQ foot, resulting in a risk of injury. No modifications of Odyssey iQ foot except those described in this information document are authorized. The Odyssey iQ foot and damaged components may only be opened or repaired by certified College Park technicians.

Caution: Battery damage

Damage can occur to the battery when dropping, knocking, crushing, vibrating, or puncturing. Avoid damaging lithium batteries and devices. Always inspect for signs of damage, such as hissing, leaking, cracking/ bulging, and smoking before use. If any of these signs are present, immediately place away from flammable materials. Remove the device from service and contact your prosthetist. Do not charge.

**Caution: Penetration of dirt and moisture into charge port**

The penetration of dirt or moisture into the charge port may lead to faulty control or malfunction of the Odyssey iQ foot and result in a risk of injury. Dirt/debris should be removed from the Odyssey iQ charge port. Ensure it is dry before charging.

**Caution: Water and humidity**

The foot can be operated in freshwater but cannot be used in extreme situations such as diving or jumping into water. It is protected against freshwater at a maximum depth of 1 meter for 30 minutes. The foot is not resistant to corrosion and should not come into contact with corrosive materials, saltwater, chlorinated water, or pH extremes..

- Use caution when walking on wet surfaces.
- Under certain circumstances, use in wet environments may increase the risk of slip and fall.
- The foot should be switched into ankle lock mode when using it in water or near wet surfaces.
- After contact with water, ensure that the water has drained out of the system and that the drainage port has not become blocked (Figure 2: drain hole in shell)
- After the foot encounters moisture, wipe it dry using a lint-free cloth.
- Make sure that the charge port is clean and dry before charging.
- If the foot encounters salt or chlorinated water, it should be rinsed immediately and dried off
- If the device malfunctions after exposure to water, contact College Park.

**Caution: Thermal overloading**

Extreme temperature conditions can lead to faulty control or malfunction of the Odyssey iQ foot and result in a risk of injury. Avoid areas outside the specified operating temperature range. The operating temperature range must be between **-10 °C and 50 °C (14 °F and 122 °F)**.

**Caution: Critically low battery status or loss of power**

The foot will go into safety mode when in a critically low state. Safety mode means that the foot goes to the set plantarflexion and dorsiflexion resistances and will no longer compensate for terrain or speed changes and won't be able to enter Stance iQ mode.

**Caution: Magnetic interference**

The Odyssey iQ foot and connected components can malfunction when near high-tension power lines, transmitters, transformers, or other sources of strong electromagnetic radiation (such as security systems for goods in department stores). This can result in a risk of injury.

Caution: Too close to HF communication devices (e.g. mobile phones, Bluetooth devices, WIFI devices)

The device may be susceptible to electromagnetic interference from portable and mobile RF communications devices such as mobile (cellular) telephones or other equipment, even if that other equipment complies with CISPR EMISSION requirements. Portable RF communications equipment (including peripherals such as antenna cables and external antennas) should be used no closer than 30 cm (12 inches) to any part of the Odyssey iQ foot, including cables specified by the manufacturer. Otherwise, degradation of the performance of this equipment could result.

The Odyssey iQ foot is suitable for use in any environment except where immersion in fluids other than freshwater is possible, or where exposure to highly electrical and/or magnetic fields can occur (e.g. electrical transformers, high power radio/TV transmitters, RF surgical equipment, CT and MRI scanners).

Caution: EM disturbances

Do not use the Odyssey iQ foot near active HF SURGICAL EQUIPMENT and the RF shielded room of an EM SYSTEM for magnetic resonance imaging, where the intensity of EM DISTURBANCES is high. High levels of Electro Magnetic disturbances may cause system to stop functioning properly, either not responding to input signal(s) or no movement of the joints.

Caution: Operating the product near active implanted systems

When operating the product, there is a risk of temporary influences of active implantable systems (e.g. pacemakers, defibrillators etc.) because of electromagnetic interference of the product.

- **When operating the product in the immediate vicinity of active implantable systems, ensure that the minimum distances stipulated by the manufacturer of the implant are observed.**
- **Make sure to observe any operating conditions and safety instructions stipulated by the manufacturer of the implant.**

Warning: Using with other equipment

Use of this equipment adjacent to or stacked with other equipment should be avoided because it could result in improper operation. If such use is necessary, this equipment and the other equipment should be observed to verify that they are operating normally.

Warning: Use only specified equipment

Use of accessories, transducers and cables other than those specified or provided by the manufacturer of this equipment could result in increased electromagnetic emissions or decreased electromagnetic immunity of this equipment and result in improper operation.

Warning: Operation of equipment in hospitals

The Odyssey iQ foot was designed for use in residential environments (home, restaurants, etc.) not hospitals or industrial areas. If device is used in environments such as hospitals or industrial areas, user might have to relocate to operate device appropriate so that is away from other HF radio devices.

Warning: Use in an oxygen rich environment

The device is not designed for use in presence of flammable gases or oxygen rich environments.

SYMBOLS

FEATURE	DESCRIPTION	FEATURE	DESCRIPTION
	Direct current	IC: 23763- RSL10SIP	IC ID of BT Module
	Caution		FCC Mark 21 CFR Part 15
	See operating instructions/consult instructions for use		Non-ionizing radiation
	Type BF Applied Part		CE Mark
	Refer to instruction manual/booklet		Device contains electronic components and/or batteries that should not be disposed of in regular waste
IP67	Dust tight and protected against freshwater at a maximum depth of 1 meter for 30 minutes.		On/Off
	Keep away from water/keep dry		Recyclable Lithium Polymer Battery
	Temperature limitation		Lithium
	Medical device manufacturer		Medical Device
	Lot number		Complies with Australia Radio communications requirements
	Serial number		Note: Possible technical damage
	Humidity Rating		Info: Basic information regarding this product
	Pressure Rating		Caution: Possible risk of accident or injury
FCC ID: 2APD9- RSL10SIP	FCC ID of BT Module		Warning: Possible risk of severe accident or injury
		CAUTION	
		WARNING	

COMPLIANCE

STANDARD	EDITION	DESCRIPTION
ISO 10328	Structural testing of lower-limb prostheses-requirements and test methods	2016
ANSI/AAMI 60601-1	Medical Electrical Equipment - Part 1: General Requirements for Basic Safety and Essential Performance	AAMI ES60601-1: 2005 & A1:2012 & A2:2021
IEC 60601-1-2	Medical electrical equipment - Part 1-2: General requirements for basic safety and essential performance - Collateral Standard: Electromagnetic disturbances - Requirements and tests	2014 Ed.4+A1
IEC 60601-1-6	Medical Electrical Equipment – Part 1-6: General Requirements for Basic Safety and Essential Performance – Collateral Standard: Usability	2010 Ed. 3+A1;A2
IEC 62366-1	Medical Devices - Part 1: Application of Usability Engineering to Medical Devices	IEC 62366-1 ED. 1.1 B:2020
IEC 60601-1-11	Medical Electrical Equipment – Part 1-11: General Requirements for Basic Safety and Essential Performance - Collateral Standard: Requirements for Medical Electrical Equipment and Medical Electrical Systems Used in the Home Healthcare Environment	2015 Ed.2+A1
IEC 62304	Medical Device Software - Software Life Cycle Processes	IEC 62304 ED. 1.1 B:2015
FCC Part 15	Radio Frequency	
IEC 62133	Secondary Cells And Batteries Containing Alkaline Or Other Non-Acid Electrolytes - Safety Requirements For Portable Sealed Secondary Lithium Cells, And For Batteries Made From Them, For Use In Portable Applications - Part 2: Lithium Systems	IEC 62133-2 ED. 1.1 B:2021
IEC 60529	Degrees of Protection Provided by Enclosures (IP Code)	2013 Ed.2.2

EMC COMPLIANCE – SPECIFIC MITIGATIONS

The Odyssey iQ foot was tested to the listed standards at the appropriate levels for Home Health Care Equipment below to ensure safety of the product regarding immunity and emissions. All devices maintained their performance during and after the tests were completed.

This device complies with part 15 of the FCC rules. Operation is subject to the following two conditions: (1) This device may not cause harmful interference. (2) This device must accept any interference received, including interference that may cause undesired operation.

This device complies with Industry Canada license-exempt RSS standard(s). Operation is subject to the following two conditions: (1) this device may not cause interference, and (2) this device must accept any interference, including interference that may cause undesired operation of the device.

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PHENOMENON & STANDARD	TEST LEVEL	REMARKS
Radiated Emissions CISPR11 ed5.0 (with A1:2010), CISPR 11 ed6.1 (2015 +A1:2016)	Group 1, Class B	<i>The Odyssey iQ foot uses RF energy only for its internal function. Therefore, its RF emissions are very low and are not likely to cause any interference in nearby electronic equipment.</i>
Electro-Static Discharge Immunity Test IEC 61000-4-2 ed2.0 (2008-12)	Contact ± 8 kV Air ± 2 kV, ± 4 kV, ± 8 kV, ± 15 kV	<i>The Odyssey iQ foot is attached to the patients socket which is design by a certified prosthetist.</i>
Radiated, Radio-Frequency, Electromagnetic Immunity IEC 61000-4-3 ed3.0 (with A1:2007+A2:2010)	10 V/m 80 MHz – 2,7 GHz 80% AM at 1 kHz	<i>Portable and mobile RF communications equipment should be used no closer to any part of the Odyssey iQ foot, including cables, than the recommended in the technical manual. Device does not need to be operated in a shielded environment.</i>
Power Frequency Magnetic Field Immunity Test IEC 61000-4-8 ed2.0 (2009-09)	30 A/m, 50 Hz or 60 Hz	<i>The Icon Knee should not be operated closer than 15 cm of sources of The Odyssey iQ foot should not be operated closer than 15 cm of sources of power frequency magnetic field.</i>

TEST SPECIFICATIONS FOR ENCLOSURE PORT IMMUNITY TO RF WIRELESS COMMUNICATIONS EQUIPMENT

TEST FREQUENCY (MHZ)	BAND (MHZ)	SERVICE	MODULATION	MAXIMUM POWER (W)	DISTANCE (M)	IMMUNITY TEST LEVEL (V/M)
385	380-390	TETRA 400	Pulse modulation 18 Hz	1,8	0,3	27
450	430-470	GMRS 460, FRS 460	FM ± 5 kHz deviation 1 kHz sine	2	0,3	28
710	704-787	LTE Band 13, 17	Pulse modulation 217 Hz	0,2	0,3	9
745						
780						
810	800-960	GSM 800/900, TETRA 800, iDEN 820, CDMA 850, LTE Band 5	Pulse modulation 18 Hz	2	0,3	28
870						
930						
1720	1700-1990	GSM 1800; CDMA 1900; GSM 1900; DECT; LTE Band 1. 3. 4, 25; UMTS	Pulse modulation 217 Hz	2	0,3	28
1845						
1970						
2450	2400-2570	Bluetooth WLAN, 802.11 b/g/n, RFIO 2450, LTE Band 7	Pulse modulation 217 Hz	2	0,3	28
5240	5100-5800	WLAN 802.11 a/n	Pulse modulation 217 Hz	0,2	0,3	9
5500						
5785						

TEST SPECIFICATIONS FOR ENCLOSURE PORT IMMUNITY TO PROXIMITY MAGNETIC FIELDS

TEST FREQUENCY	MODULATION	IMMUNITY TEST LEVEL (A/M)
30 kHz	CW	8
134,2 kHz	Pulse modulation 2,1 kHz	65
13,56 MHz	Pulse modulation 50 kHz	7, 5

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This image shows a single sheet of white paper with horizontal ruling lines. The lines are evenly spaced and run across the width of the page. There are no margins, text, or other markings on the paper.

This image shows a blank sheet of white paper with horizontal ruling lines. The lines are evenly spaced and run across the width of the page. There are no margins, text, or other markings on the paper.

This image shows a single sheet of white paper with horizontal ruling lines. The lines are evenly spaced and run across the width of the page. There are no margins, text, or other markings on the paper.

MEDENVOY SWITZERLAND

Gotthardstrasse 28, 6302 Zug, Switzerland

MADE IN THE USA

©2024 College Park Industries, Inc. All rights reserved.

1170 UG ODIQ 241112

COLLEGE PARK INDUSTRIES, INC

27955 College Park Dr. Warren, MI 48088 USA

EMERGO EUROPE

Westervoortsedijk 60, 6827 AT Arnhem, The Netherlands

Australian Sponsor

EMERGO AUSTRALIA

Level 20, Tower II, Darling Park, 201 Sussex Street,
Sydney, NSW 2000 Australia