#### MULII I ILOU I II IIIIU Q CANL IIIO INOU IIUNO

## **Prolift TLSO MODEL T108 TLSO**

# FOLLOW THESE INSTRUCTIONS FOR WEARING AND CARING FOR YOUR XBACK BRACE

# Steps for putting on and adjusting your brace:



 SLIDE YOUR ARMS THROUGH SHOULDER STRAPS AS IF YOU WERE PUTTING ON A JACKET.



3. CLOSE THE FRONT PANELS LOW AND SNUG ON THE ABDOMEN



 Pull shoulder tension straps to desired level of back support.



7. Adjust shoulder strap Pads to the most comfortable position




2. Position the brace in the center of back



4. Pull the compression straps to the desired level of support and anchor to the front panel.



6. ANCHOR SHOULDER TENSION STRAPS TO THE ABDOMINAL PANEL.



8. TRIM SHOULDER TENSION STRAPS IF TOO LONG BUT MAKE SURE ENDS OF THE STRAPS REACH THE MIDDLE OF THE ABDOMINAL PANEL BEFORE YOU CUT.

#### **Care Instructions:**

Hand wash only in cool water with a mild liquid soap and lay flat to air dry. When drying do not expose to direct heat or sunlight. Do not dry clean. If your brace has a rigid shell, the shell may be cleaned with a mild topical cleaner. Do not use strong detergents or bleaches as both may damage the brace materials.

### **Warranty:**

All XBACK braces are warranted to be free of material or workmanship defects for 5 years. Any defective product will be repaired or replaced at the discretion of XBack Bracing Services, Inc. after review of the defect by XBack Bracing Services, Inc. There are no other implied or expressed warranties.



Simply Innovative Braces™

XBack Bracing Services, Inc. P.O. Box 100 341 West Main Street Bridsboro, PA 19508 USA 610.404.4900 Phone 610.404.4905 Fax info@xbackbrace.com www.xbackbrace.com

