

NEW Flex-Foot® Balance

BY ÖSSUR®

The Flex-Foot Balance is an ideal starter foot for those who have recently undergone surgery or who have little experience on a prosthetic limb. It is ideally suited for those who walk only short distances or those with diabetes or a vascular condition.

The lightweight Balance foot incorporates a full length toe lever and integrated multi-axial design. Optimal ground compliance contributes to a balanced and natural gait. Designed for the slowest walkers, Balance allows household ambulators the advantages of advanced technological design in a composite foot geared specifically to meet their needs for comfort and stability.



Balance is compatible with the Iceross Original or Dermo liners.



Balance is compatible with the Total Knee 1900 (below left) or 2000 (below right).



Experience Our Commitment

T: 800.767.7776 x 3 • F: 800.869.7776
estore.spsco.com • www.spsco.com

KEY BENEFITS of Balance:

Lightweight foot with full length toe lever

Low build height accommodates long residual limbs

Affordable Option for low impact amputees

Easy to finish replaceable foot cover (included with the prosthetic foot)



RATINGS & SPECIFICATIONS

- Low impact levels (inhouse ambulator engaged in limited but routine activities with low impact such as bowling, gardening, shopping, etc.)
- Weight Rating: 79-214 lbs
- Foot Weight: 490 grams
- Adapter: Male Pyramid
- Build Height: 2 3/16"
- Heel Height: 3/8"

OPTIONS

- Sizes 21-28, categories 1 through 3
- Male pyramid
- Foot comes with sandal toe foot cover and a Flex-Foot Sock
- FBP0XYZ Flex-Foot® Balance w/pyramid
 - X: categories 1-3
 - YY: sizes 21-28
 - Z: foot cover L(left) or R(right): Beige - default cover color. Add BR(brown) suffix to the part number

12-MONTH WARRANTY for foot module
2-MONTH for Foot cover

LCODE: L5972 & L5986


PRODUCT FEATURES



IMPACT LEVEL



BALANCE CLEARANCE



CATEGORY SELECTION GUIDE FOR BALANCE

cm/lbs	21	22	23	24	25	26	27	28
79-97	0	0	0	0	-	-	-	-
99-134	1	1	1	1	1	1	-	-
137-176	-	-	2	2	2	2	2	2
179-214	-	-	-	-	3	3	3	3