

COMBO a ToeOFF® Family product

Fabricate a custom KAFO in less than an hour!

Low Profile and functional orthotic solution for management of genu recurvatum combined with footdrop. Combo is fitted individually to each patient. The different components are combined to accommodate leg circumference and length.

FEATURES & BENEFITS:

- Uprights easily shaped and adjusted to optimum height for patient
- Modular components allow you to accommodate variances in thigh and calf circumferences
- Low profile
- Lightweight
- Easy to don and doff
- Easily removed when wearer is relaxing (watching TV, etc.)

5 Components Required to Build a COMBO**

1. Joints with Uprights - universal size
2. Tibia Connector - generally same size as ToeOFF or BlueRocker.
3. Knee Control with Popliteal Interface - determined by knee width at the knee
4. Thigh Cuff - Take circumference of thigh at location where cuff to be worn. Can not be less than 4" to Tibia Plateau
5. Thigh Cuff SoftKIT - same as Thigh Cuff

**** TOTAL Price for all 5 COMBO elements ...\$518.00 (not including cost of ToeOFF or BlueRocker).**

COMBO refers only to the hyperextension attachment to a ToeOFF or BlueRocker. ToeOFF and BlueRocker are sold separately through SPS.

Size	1. Joints w/ Uprights
Universal	28780000

Size	2. Tibia Connector
Small	287840011
Medium	287840012
Large	287840013

Size	Knee Width	3. Knee Control
Small	2 ¼ - 4 ¼"	287820011
Medium	3 ¼ - 5 ¼"	287820012
Large	5 - 6 ¼"	287820013

Size	Thigh Circ.	Medial-Lateral	4. Thigh Cuff Left	4. Thigh Cuff Right
Small	11 ¼ - 15 ¾"	4 ¼ - 6"	287811011	287812011
Medium	15 - 19 ¾"	5 ¼ - 6"	287811012	287812012
Large	19 - 23 ¾"	6 - 7 ½"	287811013	287812013
XLarge	22 ¼ - 27 ½"	6 ¾ - 8 ¼"	287811014	287812014

Size	5, Thigh softKIT
Small	287830011
Medium	287830012
Large	287830013
XLarge	287830014



ToeOFF or Blue Rocker SOLD SEPARATELY
Please contact SPS Customer Service for part number and sizing information on ToeOFF and BlueRocker

ToeOFF or Blue Rocker softKIT SOLD SEPARATELY
Please contact SPS Customer Service for part number and sizing information on ToeOFF and BlueRocker softKIT

Technical Specifications:
Weight Limit: 330 lbs.
Height Limit: 6.5'

Use of COMBO with any product other than ToeOFF or BlueRocker is prohibited.