

Catalyst 9 takes performance and value to a new level

- 60 day trial & 36 mos. warranty
- 9 Activity/Weight Categories to fit broad range of amputee needs
- Rotation with progressive resistance
- Narrow profile
- Exceeds ISPO 22675 test standards
- Excellent value

RANGE OF MOTION

- Max Plantar Flexion: 50°
- Max Dorsiflexion: 23°
- Axial rotation: 14°

Sizes: 22-30cm

Max. Weight: 366lb

Product Weight: 449g for 25cm

Build height: 7 1/8" from floor to top of pyramid dome

Foot Cover: Caucasian, light brown or dark brown

Warranty: 3 yrs for hard components;
6 mos for foot shell & sock

Call SPS Customer Service to order:

Purchase includes foot shell, Spectra sock and 3 adjustment wedges.



Narrow shank profile

Trulife has narrowed the profile of the Seattle Catalyst 9 by 1/2", allowing a wider range of cosmesis shaping options.