

CoreLINE Tri-Panel Knee Immobilizer

- Adjustable three piece design accommodates multiple circumferences
- Lightweight aluminum posterior, medial and lateral stays are malleable to conform to the contours of the knee
- Side panels can be repositioned to optimally position the medial and lateral stays as needed for improved fit
- Posterior stay is sewn in, while lateral and medial stays are removable
- Loop lock closure with plastic D rings

Indications

- Post-operative knee support and immobilization
- Grade 1 & grade 2 collateral ligament sprains
- Patellar subluxations

Product Length	Max Thigh Cir.	Max Calf Cir	Product Number
12"	20.5"	16.5"	CL0131
14"	20.5"	16.5"	CL0132
16"	22"	18"	CL0133
18"	23"	19"	CL0134
20"	24"	20"	CL0135
22"	24"	20"	CL0136
24"	24"	21"	CL0137



LCode.....L1830*

**All coding opinions expressed herein are based upon information provided by the manufacturer. Regardless of the source of the information provided, the final responsibility for correct coding within all established laws, rules, standards and practices is the sole responsibility of the practitioner submitting the claim.*