

CoreLINE Posterior Night Splint

- Adjustable flexion straps provide a gentle stretching of the plantar fascia and Achilles tendon
- Dorsiflexion straps adjust from neutral position (90°) up to 10° dorsi-flexion
- Lightweight, sturdy one piece shell has a cool foam liner, padded calf and padded ankle strap for additional patient comfort



Indications

- Plantar fasciitis
- Achilles tendonitis
- Lower extremity overuse injuries

Size	Men's Shoe Size	Women's Shoe Size	Product Number
Small	Up to 7	Up to 7.5	CL0082
Medium	7.5 - 9	8 - 9.5	CL0084
Large	9.5 - 11.5	10 - 12.5	CL0086

LCode.....L4396*

Rear View

