

CoreLINE Static Positioning AFO

- Designed with wider toe post than other commercially available models to address a greater range of shapes and deformities. Increases patient comfort and prevents edge pressures.
- Kydex plastic shell allows re-forming and increases breakage resistance
- Universal sizing, without straps ensures easy fit and lower inventory
- Available with Fleece or Vel Foam lining



CL0330

Indications

- Prevent Heel Ulcers
- Prevent Contractures
- For recumbent positioning

Size	Calf Circumference	Fleece	Vel Foam
Regular	< 16"	CL0330	CL0331
Large	16" +	CL0335	CL0336

LCODE.....L4396*



CL0331