

CoreLINE Cock-Up Wrist Splint

- Non-removable anatomically designed aluminum dorsal stay maintains positioning with padding for comfort
- Palmer stay can be removed and adjusted/formed for improved fit
- Breathable mesh fabric is lightweight and comfortable for all day use
- Designed to fit below palmar crease for full finger dexterity
- 8" long; Fits wrists up to 9.5" maximum circumference



Indications

- Carpal Tunnel Syndrome
- Support of weak or injured wrists

Size	Wrist Circumference	Left	Right
XSmall	< 5"	CL0270	CL0280
Small	5 - 6.5"	CL0272	CL0282
Medium	6.5 - 7.5"	CL0274	CL0284
Large	7.5 - 8.5"	CL0276	CL0286
XLarge	8.5"+	CL0278	CL0288

LCode.....L3908*