

Element™ with DuraShock Kit

Shock, Rotation and Smooth Rollover in one device

- *Element Foot provides smooth rollover and superior stability*
- *DuraShock offers vertical shock and 60° of axial rotation*

FEATURES & BENEFITS:

- *Maintenance free design*
- *Less than 7" of clearance required*
- *Lightweight and flexible modular design*
- *Rated for patients up to 275 lbs*
- *Interchangeable heel elastomers for fine tuning heel stiffness*

Element Category Selection Chart

Weight (lbs)	22cm	23 cm	24 cm	25 cm	26 cm	27 cm	28 cm	29 cm	30 cm
100-119	22-A11	23-A11	24-A11	25-A32	26-A32				
120-139	22-A11	23-A11	24-A11	25-A32	26-A32	27-A43	28-A43		
140-159	22-A11	23-A11	24-A11	25-A32	26-A32	27-A43	28-A43	29-C64	30-C64
160-179	22-A22	23-A22	24-A22	25-A43	26-A43	27-A43	28-A43	29-C64	30-C64
180-199	22-A22	23-A22	24-A22	25-B43	26-B43	27-B54	28-B54	29-C64	30-C64
200-224	22-A32	23-A32	24-A32	25-B53	26-B53	27-B54	28-B54	29-D75	30-D75
225-249	n/a	n/a	n/a	25-C54	26-C54	27-C65	28-C65	29-D75	30-D75
250-275	n/a	n/a	n/a	25-C64	26-C64	27-C65	28-C65	29-D75	30-D75

Mfg's Suggested LCodes^{††}:

L5981, L5984, L5988

††All coding opinions expressed herein are based upon information provided by the manufacturer. Regardless of the source of the information provided, the final responsibility for correct coding within all established laws, rules, standards and practices is the sole responsibility of the practitioner submitting the claim.

