

EXTREME Liners (AKDT or NEW AKFR)

KEY BENEFITS of the EXTREME LINER:

- **Universal Application** - Indicated for either transfemoral (AK) or transtibial (BK) amputees. Reduce on-hand inventory cost and space requirements through application on either AK or BK users.
- **80% Less Vertical Stretch** - As compared to other ALPS gel liners, virtually eliminating pistoning while minimizing bunching behind the knee during flexion.
- **GripGEL™ provides more control and stability** - ALPS new GripGEL is more tactile than EZGel to gently but firmly support the residual limb and sensitive tissues.
- **Budget Friendly** - Provides outstanding quality, durability and functionality at an extremely affordable price.

INDICATIONS EXTREME LINER:

- Low to highly active transtibial and transfemoral users
- Users requiring firm support of soft tissue with added stability



Extreme Locking Liner (AKDT)


NEW Extreme Cushion Liner (AKFR)

Activity Level: K2, K3 and K4

Sizing: Measure 6 cm from distal end. Available in 6mm or 3mm thickness

Warranty: 12 Months for single; 24 Months for pair

LCodes: AKDT - L5673
AKFT - L 5679



AKFR Cushion Liner 6mm profile

AKFR Cushion Liner 3mm profile

AKDT Locking Liner 6mm profile

AKDT Locking Liner 3mm profile

Size	Cushion 6mm Uniform	Cushion 3mm Uniform	Locking 6mm Uniform	Locking 3mm Uniform
16-19 cm	AKFR16-6	AKFR16-3	AKDT16-6	AKDT16-3
20-23 cm	AKFR20-6	AKFR20-3	AKDT20-6	AKDT20-3
24-25 cm	AKFR24-6	AKFR24-3	AKDT24-6	AKDT24-3
26-27 cm	AKFR26-6	AKFR26-3	AKDT26-6	AKDT26-3
28-31 cm	AKFR28-6	AKFR28-3	AKDT28-6	AKDT28-3
32-37 cm	AKFR32-6	AKFR32-3	AKDT32-6	AKDT32-3
38-43 cm	AKFR38-6	AKFR38-3	AKDT38-6	AKDT38-3
44-53 cm	AKFR44-6	AKFR44-3	AKDT44-6	AKDT44-3

NEW Extreme Sleeve also available. Call SPS for sizing and pricing information.



XX in the part number indicates the size. AKDTXX-6 or AKFRXX-3