

# Shields™ Brace from Hely & Weber

- Kuhl™ perforated neoprene construction provides comfortable compression
- Dynamic patella stabilization provided through entire range of motion
- Independent calf and thigh strap application adjusts for most leg contours
- Low profile 8 1/2" design prevents brace slippage



# Titan™ Wrist - Lacing Orthosis

- Constructed of soft, comfortable orthopedic felt laminated to brush nylon exterior
- Universal right and left sizes limit inventory needs
- Adjustable pull and lacing system handles a wide range of wrist sizes
- Both palmar & dorsal stays are malleable for customization and removable for "step-down" rehabilitation



Call SPS Customer Service at 800.767.7776 for additional information on all Hely & Weber products now available at SPS.