

## Titan™ Wrist - Lacing Orthosis from Hely & Weber

- Constructed of soft, comfortable orthopedic felt laminated to brush nylon exterior
- Universal right and left sizes limit inventory needs
- Adjustable pull and lacing system handles a wide range of wrist sizes
- Both palmar & dorsal stays are malleable for customization and removable for "step-down" rehabilitation



## TITAN™ Thumb Orthosis

- Constructed of soft, comfortable orthopedic felt laminated to brush nylon exterior
- Adjustable lacing pulls easily and fits a wide range of wrist sizes
- Oblique thumb stay is malleable to allow for desired thumb position and rigid enough for excellent thumb stability
- Provides abduction control, angle adjustability and easy reapplication
- Moveable stay pod makes this brace universal and reduces inventory

## Shoulder Immobilizer

- Padded pressure sensitive straps are totally adjustable and comfortable
- Body Strap is easily adjusted or can be removed
- Deep pocket design
- Two adjustment points for each strap allows for better fit



## WEBLY® Ankle Orthosis

- Constructed of heavy gauge ballistic nylon, with non-stretch, non-fraying nylon straps
- Limits inversion/eversion while not restricting plantar and dorsiflexion
- Neoprene tongue and Achilles relief padding aid in player/patient comfort
- Use for chronic ankle sprains and prophylactically
- Low profile brace fits in shoes and work boots



## Three Panel Knee Immobilizer

- Provides rigid support for the knee
- 3/8" pile foam laminate with black flannel interior ensures patient comfort
- Two detachable panels for circumferential adjustments for custom fitting
- Posterior and medial/lateral stays are malleable for custom fitting
- Optional patella strap provides extra knee control
- Available in 8 lengths from 12" - 26" (2" increments)



## Shields™ Brace

- Kuhl™ perforated neoprene construction provides comfortable compression
- Dynamic patella stabilization provided through entire range of motion
- Independent calf and thigh strap application adjusts for most leg contours
- Low profile 8 1/2" design prevents brace slippage



Call SPS Customer Service at 800.767.7776 for additional information on all Hely & Weber products now available at SPS.