


PATENTED SPHERICAL AXIS (BALL & SOCKET) JOINT



FC2 ANKLE

70° ANTERIOR-POSTERIOR ROM
26° ADJUSTABILITY AT ANKLE AXIS
PRE-CONTOURED FOR LOW PROFILE

PATENTED MANUFACTURED IN USA STAINLESS STEEL HEAT TREATED

KEY FEATURES OF FC2 ANKLE

- Plantar/Dorsi Flexion Stop Set-Screws can be set up to 70° and adjusted while the patient is wearing the orthosis
- Reduces fabrication time with SELF-SQUARING Molding Dummy

