

Lerman MINERVA™ III

NEW Lerman Minerva III from Trulife provides comfortable cervical and upper thoracic spine stabilization with added occipital support.

- Occipital plate cradles head and neck
- Easily accessible c-shaped trachea opening
- Ventilated shell
- No tools required to adjust height
- Universal design reduces inventory
- Height (chin to distal edge): 12 - 15"

INDICATIONS:

- C3-T2 Stable Fractures
- Muscular/Ligamentous Injury
- Degenerative Joint Disease
- Cervical Pain/Nerve Impingement

CONTRAINDICATIONS:

- Unstable Fractures
- C1, C2 Fracture

Mfg's Suggested L-code^{**}: L0200

***All coding opinions expressed herein are based upon information provided by the manufacturer. Regardless of the source of the information provided, the final responsibility for correct coding within all established laws, rules, standards and practices is the sole responsibility of the practitioner submitting the claim.*

