

LowPro Summit Locks ONLY \$99

A savings of over 16%

An improved design of the original zero clearance, rotation controlling, proximal lock mechanism, the **LowPro Summit Lock** offers much easier engagement.

CD102 LowPro Summit with bowtie regularly \$119.00

CD102M LowPro Summit with Manta regularly \$129.00



New Bowtie liner attachment shown is easier to roll-on. Manta shape is also available.

New LowPro lock is glued to the inside of the socket, making it much easier to engage and cover.

New tab is easier to engage.

No space is required for the lock distally, maximizing clearance for the knee center. Use the socket adaptor of your choice.

Coyote Design

Available thru



(800) 767-7776 x 3

www.spscocom

Valid for CD102 and CD102M. Offer not available with any other specials or discounts. Available only July 1, 2009 - July 31, 2009. Applies only to CD102 and CD102M.

Patent# 6666894. Other patents pending.