

GET THE MOST FROM YOUR LAMINATIONS



NEW NANO-Res uses nano scale elastomeric elements that function like billions of shock absorbers to help attenuate impact energy. The nano elements create a lamination with unmatched strength and durability. Easy to use, the NANO-Res has a 12-15 min. gel time without bubbles and very good fiber wet out for consistent reliable results.



NEW EAR1 has an optimized molecular weight that lowers its viscosity and surface tension versus competing products. Net result is laminates with superior fiber wet out, equating to high translation of fiber properties. Easy to use, the EAR1 has a 10-15 min. gel time without bubbles to produce consistent reliable results quickly.

T (800) 767.7776 x 3 • F (800) 869.7776

© 2008 SPS. SPS and Experience Our Commitment are trademarks of SPS. All other trademarks are properties of their respective owners. All rights reserved.

SPS

Experience Our Commitment

AONE-1108r0