

NEW RHEO KNEE® BY OSSUR®

THE REDESIGNED RHEO KNEE FEATURES A HOST OF UPGRADES, INCLUDING IMPROVED AESTHETICS, A HIGHER WEIGHT LIMIT, INCREASED TORQUE AND FASTER SWING SPEED.

A truly dynamic knee, the new RHEO KNEE uses Artificial Intelligence to continuously adapt to the user's walking style and environment. The RHEO KNEE recognizes and responds immediately to changes in speed, load and terrain, restoring the user's ability to walk naturally, comfortably and confidently at any speed.

NEW FEATURES

- Higher weight limit of 275 pounds
- Increased torque capabilities
- Streamlined, anatomical shape
- Reduced swing resistance
- Faster swing speed
- Enhanced swing phase programming
- Wireless Bluetooth communication



T (800) 767.7776 x 3 • F (800) 869.7776 • www.spsco.com


Experience Our Commitment

RHEO KNEE®

INDICATIONS

- Community and workplace ambulators with the ability or potential for cadence variation and ramp/stair descent
- Low to moderate impact levels

FEATURES & BENEFITS

- 275 lbs weight limit allows use by heavier and/or more active users
- Anatomical shape simplifies cosmetic finishing
- Smoother, more consistent swing resistance produces faster swing speed and improved user confidence
- 10% lighter than the original RHEO KNEE
- Wireless BlueTooth communication with PDA (compatible with existing RHEO KNEE PDA devices after PDA software update)
- Continuous adaptation and dynamic learning through enhanced artificial intelligence
- Automatically adapts to more user gaits and cadences for added flexibility and user satisfaction (no manual adjustment for toe loading required)
- Proportional stance resistance enables more natural motion and freedom of movement
- Can be used with full height Flex-Foot® prosthetic feet

TECHNICAL SPECIFICATIONS

- Weight limit: 275 lbs
- Knee flexion available: 120 degrees
- Weight: 3.36 lbs
- Build Height: 9 1/4"
- Extension Assist: Internal Torsion Spring
- Battery life: 24 - 48 hrs (depending on use)
- Recharge time: 2 - 4 hrs
- Attachments:
 - Proximal: male pyramid
 - Distal: male pyramid
- Warranty/Service:
 - 36 month limited warranty
 - Includes service at 20 months



T (800) 767.7776 x 3 • F (800) 869.7776 • www.spsco.com

© 2008 SPS. SPS and Experience Our Commitment are trademarks of SPS.
All other trademarks are properties of their respective owners. All rights reserved.

SPS
Experience Our Commitment

ORK2-0909r1