

NEW ALPHA ORIGINAL MP LINERS

Same great attributes, design and materials you've come to know and love in the ALPHA Original Liner at a new fabulous price!

KEY BENEFITS ALPHA ORIGINAL MP LINERS:

- Surrounds limb with skin-friendly, mineral-oil based gel containing Vitamin E
- Reliable performance for a wide range of amputees
- Heat moldable to facilitate customization within practitioner's facilities
- Distinct gel patterns accommodate various TT limb sizes

INDICATIONS:

- Low to high activity transtibial users
- Users needing skin protection against abrasion and breakdown

SIZES:

- 6 Sizes (Small to Extra Large) accommodate proximal circumferences of 9" to 28" - see fitting charts below for proper selection

LCODE^{††}: L-5679 (cushion), L-5673 (locking liner)



Profile view of new Alpha Original MP liner available in grey/green 6mm UNIFORM only



Original MP Liner Sample Part Number
Model: Cushion (ALC) or Locking (ALL)
Style: 5 (Uniform)
Fabric: 0 (Gray/Green)
Thickness: 6 (6mm)
Size: 3 (Small), 0 (Medium), 4 (Medium Plus), 6 (Large),
8 (Large Plus) or 7 (Extra Large)
Single Liner: E (Each)

ALC - 5 0 6 6 - E

ALC	_____	_____	_____	_____	_____	_____
5	_____	_____	_____	_____	_____	_____
0	_____	_____	_____	_____	_____	_____
6	_____	_____	_____	_____	_____	_____
6	_____	_____	_____	_____	_____	_____
E	_____	_____	_____	_____	_____	_____

^{††}All coding opinions expressed herein are based upon information provided by the manufacturer. Regardless of the source of the information provided, the final responsibility for correct coding within all established laws, rules, standards and practices is the sole responsibility of the practitioner submitting the claim.

NEW ALPHA WARRANTY & RETURN POLICY

We are pleased to announce SPS has implemented a new Warranty and Return Policy on all Ohio Willow Wood Alpha Liners purchased through SPS.

Critical elements of the new policy are:

- **Each Alpha Liner will now have a warranty of twelve (12) months from fit date or SPS invoice date.**
(This revision extends the previous Alpha liner warranty by 6 months per liner)
- **Alpha Liners have a ten (10) day trial period from the fit date or SPS invoice date.**

This policy applies to any Alpha Liner fitted or purchased after January 1, 2010, including the new ALPHA Original MP liners showcased above. Bags opened for the purpose of patient trial and fitting may be returned within the ten day trial period.

For more information, contact SPS Customer Service at 800.767.7776 x 3 or your SPS Sales Manager.