

# NEW Plie® 2.0 MPC Knee



Front View



Three-Quarter View



Back View

The new second generation Plie® 2.0 MPC Knee provides a seamless blend of performance and safety, without sacrifice. The synergistic microprocessor and hydraulic technologies deliver a highly responsive prosthetic knee, while facilitating natural and intuitive mobility. Microprocessor control of both swing and stance phases delivers advanced stumble recovery capability. These combined advances result in outstanding user confidence and satisfaction.

Freedom Innovations' game-changing Plie® Set-Up Wizard eases performance optimization and increases fitting efficiency. The new Wizard streamlines programming and set-up of the knee, utilizing automated parameter settings which ensure optimal performance.

Please check with your SPS Sales Manager concerning availability and training.

## Features & Benefits:

- Microprocessor Controlled Swing and Stance
- Finely Tuned Control System delivers near instantaneous adaptation to gait dynamics
- Unique Set-Up Wizard utilizes automated settings to optimize performance
- Optimized Swing Flexion Resistance
- Replaceable Protective Covers

## Technical Specifications:

- Battery Life: Greater than 24 hours depending upon individual activity and storage practices
- Knee Flexion available: 125°
- Product weight: 2.6 lb

## Selection Criteria:

- Transfemoral amputees
- US K3 ambulator (all daily activities)
- Community and workplace ambulators with the ability/potential for varied cadence and ramp/stair descent
- Clearance 9 ¼"
- User weights up to: 275 lb moderate activity  
200 lb high activity

Mfg's Suggested LCodes<sup>††</sup>: L5828, L5845, L5848, L5850, L5856, L7367, L7368 and L5930 (if applicable)

<sup>††</sup>All coding opinions expressed herein are based upon information provided by the manufacturer. Regardless of the source of the information provided, the final responsibility for correct coding within all established laws, rules, standards and practices is the sole responsibility of the practitioner submitting the claim.

## SPS 4 Hour Plie 2.0 MPC Training Courses for 2010

DATE	LOCATION
Apr 13, 2010	Anaheim, CA
May 12, 2010	Dallas, TX
May 14, 2010	Chicago, IL

To Register, go online at <http://www.spsco.com/services/educational.asp> or call Kim Fell at 800.767.7776 x 173