

# DISCOVER A NEW WAVE OF COMFORT



## » Wave Comfort Foot System™

**Tri-Spring Design** - compliments the individual's natural gait by allowing each spring to flex and extend in progressive rhythm while providing a supple feel

**Wave Spring Heel**- provides stable progression and fluid transition from heel-strike to toe-off via the unique heel rocker design

**Low Profile Design**- 4" clearance and lightweight structure (13.2 oz for 26cm) facilitate wide utilization for low to moderate activity users (weighing up to 275 lbs)



» **Call SPS Customer Service  
for more information today**

