


NEW i-LIMB™ Socket Components

Socket Components Kits

Available in Single or Dual Site Kits for 50 Hz or 60 Hz

Each Kit Contains ONE (1) each of the following:

- 50 mm Lamination Ring
- 45 mm Lamination Ring
- 50 Hz (or 60 Hz) Electrode
- 3 Way Cable (11.8") with IDC block
- Coupling piece with large lock ring
- Coaxial plug with small lock ring




50 Hz Single Site Socket Kit (PL091067A)

50 Hz Dual Site Socket Kit (PL091068A)

60 Hz Single Site Socket Kit (PL091069A)

60 Hz Dual Site Socket Kit (PL091070A)



Switch Block*, p/n SA069080B

Allows the terminal device to be powered on and off remotely (3" flying leads)

**may be purchased separately*



Switch Block Power Cable*, p/n SA000219A

Connects the Switch Block (SA069080B with 3" flying leads) to the terminal device via the 4 Way Coaxial Plug (p/n PL091036A)

**may be purchased separately*



Quick Wrist Disconnect, p/n SA000128B

Allows for easy engagement / disengagement between the hand and socket. Provides users with the ability to supinate and pronate their hand manually. NOT Sold separately.



Coupling Piece*, p/n PL091032A

Large Lock Ring* (p/n 091005A) Included. May purchase rings in packets of 10/pkg.

** Coupling piece may be purchased separately*




4 Way Coaxial Plug*, p/n PL091036A

Connects one or two electrodes and the battery to the terminal device.

Plug composed of:
4 Way Coaxial Plug
Plastic Screw and Washer
Small Lock Ring


** Complete coaxial plug may be purchased. Components not sold separately.*



4 Way Coaxial Socket, p/n SA000127A

Connects two electrodes and the battery to the i-LIMB™ Hand.

Connection Socket for PL091036A 4- Way Coaxial Plug; **NOT sold separately.**



2 Way Power Cable, p/n 091004A

Connects batteries to the terminal device via the 4 Way Coaxial Plug, p/n 091036A
Length of cable: 11.8"



3- Way Cable*

Available in three sizes:

11.8"	p/n PL091029A
23.6"	p/n PL091030A
39.4"	p/n PL091031A

Connects electrodes to 4 Way Coaxial Plug p/n (PL091036A)

Includes IDC Block (p/n091007A); *May purchase IDC Block and individual cables separately.*

SPS 4 Hour i-LIMB Training Courses for 2010

DATE

Apr 13, 2010
May 12, 2010
May 14, 2010

LOCATION

Anaheim, CA
Dallas, TX
Chicago, IL

To Register, go online at <http://www.spsco.com/services/educational.asp> or call Kim Fell at 800.767.7776 x 173