

Darco Web™ Ankle Support

Comfortable support for the ankle.

Special Features:

- The open back ankle support uses elastic laces with Velcro® tabs to compress and secure the ankle.
- Low profile design with space age materials comfortably restrict ankle motion.
- Thin enough to wear with athletic shoes or casual dress shoes.
- Constructed of a durable vinyl shell lined with high tech spacer padding to help wick off moisture for comfort.
- Fits either the right or left ankle.

Indications:

- Chronic Weak Ankles
- Post Ankle Trauma
- Post-op Ankle Recovery
- Prophylactic Ankle Protection
- Ankle Support



Latex Free.

| ITEM NO. | DESCRIPTION | COLOR | SIZE | WOMEN | MEN |
|----------|-------------|-------|------|---------|---------|
| WAS0B | Web Ankle | Black | XS | 4 - 6 | 3 - 5 |
| WAS1B | Web Ankle | Black | S | 7 - 9 | 6 - 7 |
| WAS2B | Web Ankle | Black | M | 10 - 11 | 8 - 10 |
| WAS3B | Web Ankle | Black | L | 12 - 13 | 11 - 12 |
| WAS4B | Web Ankle | Black | XL | 14 - 15 | 13 - 14 |

 **SPS**

Tel: 800.767.7776 x 3

Fax: 800.869.7776

estore.spsco.com

www.spsco.com