



Experience Our Commitment

Next Step Suspension Belt

The Next Step Belt . . .

provides the trans-femoral amputee suspension, increased comfort and optimal rotational control. Designed by a Certified Prosthetist, the body of the Next Step Belt features soft textured neoprene against the skin for maximum grip and a nylon outer layer for durability. The unit is bound with stretch binding that assists in "cupping" the hip for greater patient comfort and purchase. Adjustable and easy to don. The straps are made of heavy-duty 2" nylon webbing with hook and loop closure. Available in eight sizes, the unit is attached to the prosthesis by means of chafes(*).



(*) not included with Next Step Belt



Next Step Belt

- Designed to provide the trans-femoral amputee optimal comfort, suspension and rotational control
- Ideal form of primary suspension, ancillary rotational control for lanyard and liner wearers
- Manufactured from latex-free neoprene with soft-grip textured neoprene against the skin for maximum grip with durable nylon outer layer and heavy-duty nylon straps with contact closure

Suggested
L-5698
Billing Code

Description	Catalog #	Description	Catalog #
Next Step Belt 2" x 24"	89-1450	Next Step Belt 2" x 32"	89-1454
Next Step Belt 2" x 26"	89-1451	Next Step Belt 2" x 34"	89-1455
Next Step Belt 2" x 28"	89-1452	Next Step Belt 2" x 36"	89-1456
Next Step Belt 2" x 30"	89-1453	Next Step Belt 2" x 38"	89-1457

MEASURING AND MOUNTING: attach two 2" chafes to your desired location. Measurements are taken anterior to posterior, similar to measuring for a Silesian belt. (Formula for sizing: measure from anteriorly-mounted chafe across opposite hip to posterior chafe and subtract 6" from measurement).