

SPS has expanded its offering of Otto Bock products to include a large number of items with increased weight ratings (noted by **blue font** in the tables below), which are featured in this issue. Please call SPS Customer service for new lower price points and look for your copy of the SPS mini-catalog of the Otto Bock product offering in your mail. If you need additional copies, please contact Kim Martin via kmartin@spsco.com.

DOUBLE ADAPTERS



4R72=75
4R72=60

4R72=45
4R72=32

Part Number	4R72=32	4R72=45	4R72=60	4R72=75
Weight	85 g	95 g	110 g	125 g
Min. Clearance	32 mm	45 mm	60 mm	75 mm
Material	Titanium		Titanium	
Weight Limit	330 lbs	330 lbs	330 lbs	330 lbs

TUBE CLAMP ADAPTERS Ø 30 MM



4R52=5

4R69=5

4R21=5

4R56=1
4R56

Part Number	4R52=5	4R69=5	4R21=5	4R56	4R56=1
Diameter	30 mm	30 mm	30 mm	30 mm	30 mm
Weight	75 g	75 g	130 g	85 g	85 g
Min. Clearance	48 mm	47 mm	48 mm	53 mm	58 mm
Material	Titanium	Aluminum	Stainless Steel	Titanium	
Angle	-	-	-	10°	20°
Weight Limit	300 lbs	220 lbs	275 lbs	220 lbs	220 lbs

PYLON TUBE Ø 30 mm



2R38=5

2R49

Part Number	2R38=5	2R49
Diameter	30 mm	30 mm
Tube Length	440 mm	400 mm
Weight	275 g	240 g
Min. Clearance	45 mm	46 mm
Material	Titanium	Aluminum
Weight Limit	300 lbs	220 lbs

SOCKET ADAPTERS



4R54=5
4R74
4R23

4R77

4R73=D-5
4R73=A-5.1

4R51

Part Number	4R54=5	4R74	4R23	4R77	4R73=A-5.1	4R73=D-5	4R51
Weight	50 g	55 g	85 g	70 g	60 g	60 g	80 g
Min. Clearance	7 mm	12 mm	7 mm	11 mm	10 mm	7 mm	17 mm
Material	Titanium	Aluminum	Stainless Steel	Titanium		Titanium	
Pyramid Position	-	-	-	-	Axially offset 7mm	Diagonally offset 5 mm	-
Weight Limit	300 lbs	220 lbs	275 lbs	330 lbs	330 lbs		330 lbs

SOCKET ATTACHMENT BLOCK



5R1=6-H

Part Number	5R1=6-H*
Wood Diameter	120 mm
Weight	240 g
Min. Clearance	24 mm
Material	Wood
Weight Limit	330 lbs

* Includes lamination dummies

LAMINATION ANCHORS



4R63=5.1
4R63=5

4R42

4R41


4R43

4R89

Part Number	4R63=5*	4R63=5.1*	4R42	4R41	4R43	4R89
Weight	95 g	95 g	130 g	170 g	95 g	180 g
Min. Clearance	11 mm	11 mm	14 mm	20 mm	-	15 mm
Material	Stainless Steel			Stainless Steel		
Weight Limit	275 lbs	275 lbs	330 lbs	275 lbs	275 lbs	275 lbs

* Includes Lamination dummies

SPECIAL ADAPTERS



4R101

4R118

Part Number	4R101	4R118
Weight	205 g	75 g
Min. Clearance	25 mm	10 mm
Material	Aluminum	Aluminum
Weight Limit	220 lbs	275 lbs
Sliding Range	± 11 mm	10 - 25 mm

Tel: 800.767.7776 x 3 • Fax: 800.869.7776
www.spsco.com • estore.spsco.com

© 2011 Southern Prosthetic Supply. SPS and Experience our Commitment are trademarks of SPS. All others are the properties of their respective owners.

FOOT ADAPTERS



2R31=M10-5

2R8=M10

Part Number	2R31=M10-5	2R8=M10
Material	Titanium	Steel
Weight	70 g	130 g
Min. Clearance	7.5 mm	7.5 mm
Bolt Accessories	2D7=M10	2D6=M10
Weight Limit	300 lbs	275 lbs

Use M10 models for 1S49, 1S66, 1S90 (all sizes) and 1S67, size 26

6R6= Foam Covers



- For Knee Joints: 3R21, 3R23, 3R30, 3R32, 3R46, 3R60=KD, 3R60-PRO=KD and 3R106=KD
- Appropriate for modular transtibial and knee disarticulation prostheses
- Drilled with center hole compatible for either 30 mm or 34 mm diameter systems
- Specifications:
 - Universal (left or right)
 - Length = ~55 cm
 - Material: PUR flexible foam
 - Color: Beige

3R6= Foam Cover



- For Knee Joints: 3R15, 3R49, 3R55, 3R21, 3R23, 3R30, 3R32, and 3R46
- Appropriate for modular transfemoral prostheses
- Drilled with stepped center hole
- Pre-shaped
- Specifications:
 - Select left (L) or right (R)
 - Length = ~94 cm
 - Calf Circumference: 36 cm, 40 cm or 44 cm
 - Knee Flexion: 30°
 - Material: PUR flexible foam
 - Color: Beige

Specify side (L/R) and size when ordering; example 3R6=L40

6R8= Foam Covers



- Appropriate for modular transtibial prostheses
- Drilled with center hole for 30 mm or 34 mm diameter systems
- Specifications:
 - Universal (left or right)
 - Specify tube dia. of 30 mm (6R8=30) or 34 mm (6R8=34)
 - Length = ~48 cm
 - Material: PE foam
 - Color: Beige

Specify tube diameter when ordering; example 6R8=30

3S26= Cosmetic Foam Cover (AK)



- For Knee Joints: C-Leg® (3C88-1 and 3C98-1) Compact (3C86 and 3C96)
- Appropriate for modular transfemoral prostheses
- Drilled with stepped center hole
- Pre-shaped
- Specifications:
 - Select left (L) or right (R)
 - Length = ~94 cm
 - Calf Circumference: 44 cm
 - Knee Flexion: 30°
 - Material: PUR flexible foam
 - Color: Beige

Specify side (L/R) when ordering; example 3S26=R

3R24= Foam Covers



- For Knee Joints: 3R17, 3R20, 3R22, 3R33, 3R34, 3R36, 3R40, 3R95 and 3R95=1
- Appropriate for modular transfemoral prostheses
- Drilled with stepped center hole
- Pre-shaped
- Specifications:
 - Specify left (L) or right (R)
 - Length = ~94 cm
 - Calf Circumference: 36 cm, 40 cm or 44 cm
 - Knee Flexion: 30°
 - Material: PUR flexible foam
 - Color: Beige

Specify side (L/R) and calf circumference when ordering; example 3R24=L36

6R18= Cosmetic Foam Cover (PE) pre-shaped



- Appropriate for modular transtibial prostheses
- Drilled with 30 mm or 34 mm hole diameters
- Pre-shaped
- Specifications:
 - Tube Dia: 30mm (6R18=30) or 34 mm (6R18=34)
 - Universal (left or right)
 - Length = ~48 cm
 - Material: PE foam
 - Color: Beige

Specify tube diameter when ordering; example 6R18=34

99B14= Perlon Cosmetic Stockings



- Intended for outer cosmetic finish of modular knee disarticulation and transfemoral prostheses

Size	0	1	2	3	4
Length	~59 cm	~61.5 cm	~66 cm	~72 cm	~83.5 cm
Ankle Radius	~9 cm	~9 cm	~9 cm	~9.25 cm	~9.5 cm
Foot Length	~18cm	~19 cm	~20 cm	~20.5 cm	~21 cm
Collar Radius	~15 cm	~16 cm	~17 cm	~18 cm	~19 cm
Colors	tan (-) or brazil (B)				

99B116= Soft Touch Protective Outer Surface Cover



- Special Surface coating adds water and dirt resistance
- Extremely stretchable
- Fit transtibial prosthesis
- Specify size and color when ordering; example 99B116=214

Size	Small (2)	Med (4)	Large (6)
Calf Cir.	30-33 cm	34-39 cm	40-45 cm
Ankle Cir.	19-21 cm	22-25 cm	26-29 cm
Foot Length	≥ 21cm	≥ 23 cm	≥25 cm
Colors	0, 2, 4, 6, 8 10, 12, 14, 16 or 18		

99B16= Perlon Knee Stockings



- Intended for outer cosmetic finish of modular transtibial prostheses
- Banded collar radius is ~11.5 cm for all sizes

Size	1	2	3
Length	~38.5 cm	~40.5 cm	~44.5 cm
Ankle Radius	~9 cm	~9.2 cm	~9.25 cm
Foot Length	~17.5 cm	~18.5 cm	~19 cm
Colors	tan (-) or brazil (B)		