

\$100 Price Reduction** on Flex-Foot® Talux®

KEY BENEFITS of Talux:

- Fluid, natural walking motion**
- Multi-Axial Range of Motion**
- Achilles strap** enhances forward propulsion

IMPACT LEVEL



OPTIONS

- Available in sizes 23-30, Categories 1-8
- Choice of male pyramid or 30mm carbon tube


36-MONTH WARRANTY



**Limited Time offer available through March 31, 2010

***All coding opinions expressed herein are based upon information provided by the manufacturer. Regardless of the source of the information provided, the final responsibility for correct coding within all established laws, rules, standards and practices is the sole responsibility of the practitioner submitting the claim.*

TALUX



CATEGORY SELECTION GUIDE FOR TALUX - LOW ACTIVITY

cm/lbs	23	24	25	26	27	28	29	30
100-115	1	1						
116-130	1	1						
131-150	2	2	2	2				
151-170	3	3	3	3	3	3		
171-195	4	4	4	4	4	4	4	4
196-220	5	5	5	5	5	5	5	5
221-255			6	6	6	6	6	6
256-285			7	7	7	7	7	7
286-325						8	8	8

CATEGORY SELECTION GUIDE FOR TALUX - MODERATE ACTIVITY

cm/lbs	23	24	25	26	27	28	29	30
100-115	1	1						
116-130	2	2	2	2				
131-150	3	3	3	3	3	3		
151-170	4	4	4	4	4	4	4	4
171-195	5	5	5	5	5	4	5	5
196-220			6	6	6	6	6	6
221-255			7	7	7	7	7	7
256-285						8	8	8