

different

ARC 2.0

ADJUSTABLE POSITIONING

Post-op requirements vary, therefore the aluminum support can be positioned in a variety of ways including:



0° TO FULL INTERNAL ROTATION



UP TO 70° EXTERNAL ROTATION




15° – 45° ABDUCTION

2.0 Material

New Bledsoe 2.0 Material, with dual layer fabric, keeps patients cool.

Its moisture management system captures moisture and moves it away from the skin to another layer of fabric for quick evaporation, drying the skin faster than typical breathable materials.



2.0 MATERIAL CLOSE UP



Shoulder BRACES

ARC 2.0

Size	Part Number
U	AE050400

ARC 2.0 WITH PILLOW

Size	Part Number
U	AE050500

ORIGINAL ARC WITH PILLOW

Size	Standard	Cool
U	AE050300	AE055300

UNIVERSAL ARM SLING

Size	Part Number
U	SA201000

SLING AND SWATHE

Size	Part Number
U	SA204000

CLINIC SHOULDER IMMOBILIZER

Envelope Length	Size	Part Numbers
11"	XS	SA200501
13.5"	S	SA200503
15.5"	M	SA200505
17.5"	L	SA200507
20"	XL	SA200509

CLINIC ARM SLING

Envelope Length	Size	Part Numbers
11"	XS	SA200001
13.5"	S	SA200003
15.5"	M	SA200005
17.5"	L	SA200007
20"	XL	SA200009



BLED SOE BRACE SYSTEMS
2601 Pinewood Drive
Grand Prairie, TX 75051

UNITED STATES

Toll Free (888) 253-3763
Tel (972) 647-0884
Fax (972) 660-5495

INTERNATIONAL

Tel (972) 522-7484
Fax (972) 660-5495

www.bledsoebrace.com
www.soebledsoe.com

FEATURED ON THE COVER - ARC 2.0

Key: 1. Padded neck pad 2. Deep sling 3. Moisture wicking material 4. Cryotherapy cutout 5. Quick-pull tabs
6. Breathable material 7. Pistol grip 8. Adjustable positioning 9. Malleable waistband 10. Quick-release buckles



ShoulderBraces



it's soe bledsoe

ARC 2.0

ABDUCTION ROTATIONAL CONTROL

Shoulder braces should not move and migrating pillows can cause the shoulder joint to shift, resulting in unnecessary pain during rehabilitation. Our signature aluminum waistband is moldable to each patient's unique torso shape and prevents anterior brace migration that is common with most pillows. Our newest member of the ARC family, the ARC 2.0, now features 2.0 Material for maximum breathability and patient comfort.

- » Universal sling design folds to fit every patient with one brace—right or left, football player or gymnast
- » Unique pistol grip adjusts with quick-pull tabs and keeps the arm from migrating forward out of the sling
- » Optional underarm strap relieves pressure on the neck for larger patients and patients positioned in external rotation

INDICATIONS:
For non-surgical or post-operative support and immobilization of the injured shoulder when varying degrees of internal rotation, external rotation and abduction are required. Positions include gunslinger, neutral plane, and statue of liberty.



PISTOL GRIP



QUICK-PULL TABS



ONE HAND BUCKLE



COLD THERAPY CUTOUT

ARC 2.0

WITH PILLOW

The ARC 2.0 with Pillow with our signature ARC sling and new 2.0 Material in a traditional pillow design.

- » Universal sling design folds to fit every patient with one brace—right or left, from football player to gymnast
- » Unique pistol grip adjusts with quick-pull tabs and keeps the elbow seated in the sling and prevents the hand from migrating forward out of the sling
- » Optional underarm strap relieves pressure on the neck for larger patients and reduces internal rotation migration
- » Traditional pillow design with 15° Abduction Pillow

INDICATIONS:
For non-surgical or post-operative support and immobilization of the injured shoulder where glenohumeral positioning is key.

ORIGINAL ARC

WITH PILLOW

Our signature ARC sling, with the comfort and durability you have come to expect from Bledsoe, in an economical design.

- » Breathable foam sling design folds to fit every patient with one brace—right or left, football player or gymnast
- » Unique pistol grip keeps the arm from migrating out of the sling
- » Traditional pillow design with 15° Abduction Pillow

AVAILABLE:
Standard Sling
Cool Sling is perforated for added breathability

INDICATIONS:
For non-surgical or post-operative support and immobilization of the injured shoulder where glenohumeral positioning is key.



Clinic Shoulder Immobilizer

- » Immobilizing strap secures shoulder in place
- » Oversized envelope accommodates large casts and dressings
- » Adjustable padded strap with easy-touch closure on front and back
- » Thumb loop keeps the elbow seated in the sling

INDICATIONS:
For non-surgical or post-operative support and immobilization of the shoulder, mild sprains and strains of the shoulder, or for fractures of the arm.



Clinic Arm Sling

- » Oversized envelope accommodates large casts and dressings
- » Adjustable padded strap with easy touch closure on front and back
- » Thumb loop keeps the elbow seated in the sling

INDICATIONS:
For non-surgical or post-operative support and immobilization of the shoulder, mild sprains and strains of the shoulder, or for fractures of the arm.



Universal Arm Sling

- » No trimming necessary, simply fold to proper length for most patients
- » Large envelope accommodates large casts and dressings
- » Adjustable slide buckle closure

INDICATIONS:
For non-surgical or post-operative support and immobilization of the shoulder, mild sprains and strains of the shoulder, or for fractures of the arm.



Sling and Swathe

- » Lightweight, soft foam sling supports the arm
- » Swathe immobilizes the arm and shoulder
- » Adjustable slide buckle straps

INDICATIONS:
For non-surgical or post-operative support and immobilization of the shoulder, mild sprains and strains of the shoulder, or for fractures of the arm.

