

ODYSSEY iQ®

microprocessor foot
lightweight, low-profile design
simple setup



Work smarter, not harder. The new Odyssey iQ is a hydraulically dampened microprocessor foot that is lightweight and incredibly low-profile. With an IP67 rating, long-lasting battery, and fast responses, it accommodates a wide range of users and activities. It also features simple setup and a unique method of spring-assisted dorsiflexion. Users have the option to lock the hydraulics for maximum stability in quiet standing or engage the Stance iQ mode. The Odyssey iQ was created to be user-friendly for patients and clinicians—providing access to valuable functions in College Park's Stride Studio app.

K3
iQ

TECH SPECS

FUNCTION

Microprocessor

BUILD HGT

23-25 cm | 3.5 in
26-30 cm | 3.6 in

MOUNTING

Endo (IP)

WEIGHT LIMIT

275 lbs

L-CODE#

L5973

PDAC APPROVED!

WATER CAPABILITY



FRESH WATER PROOF

SIZES

23-30 cm

WEIGHT*

1027 g

IMPACT LVL



WARRANTY

3 years

(3 year extended warranty available)

ENVIROSHELL†



standard



*26 cm model, including Enviroshell. ‡The listing of codes with these products should not be construed as a guarantee for coverage or payment. Ultimate responsibility for the coding of services/products rests with the individual practitioner.
† Caucasian, Tan, Brown, or Jet Black.

ODYSSEY iQ PART NUMBER

OD iQ SIDE SHELL COLOR* SIZE FIRMNESS CATEGORY

L OR R **C, T, B, OR J** **EN** **23-30** **D** **1-5**

*(C)aucasian, (T)an, (B)rown, or (J)et Black. Shell width is standard.
Example part number: OD iQ LCEN 24 D3

ODYSSEY iQ FIRMNESS CATEGORIES

PATIENT WEIGHT	0-140lbs	141-180	181-220	221-275
SIZE	23-30cm			
LOW IMPACT	1	2	3	4
MODERATE IMPACT	2	3	4	5

INCLUDED

Endo foot includes:
sandal toe foot shell,
CPI sock, wedge kit,
charging kit, and
remote charging kit.

ENVIROSHELL PART NUMBER


OD iQ S SIDE SHELL COLOR** SIZE

L OR R **C, T, B, OR J** **23-30**

** (C)aucasian, (T)an, (B)rown, or (J)et Black.





Enviroshell included with every foot order









HARDWARE AND TOOLS

	PART NUMBER	DESCRIPTION
	CP SSET	Super Set Screws (8pc) M8 x 14mm
	TFH 800	Foot Horn, Stainless Steel
	CP HK 4	Hex Key, 4mm


REPLACEMENT PARTS

	PART NUMBER	DESCRIPTION
	HZ HWK 2125	Horizon Wedge Kit, 21-25cm
	HZ HWK 2630	Horizon Wedge Kit, 26-30cm
	GA GEL	Cyanoacrylate Adhesive, 409 Gel (CA glue)


ACCESSORIES

	PART NUMBER	DESCRIPTION
	OD iQ WAR3	Odyssey iQ Microprocessor Foot 3 Year Extended Warranty
	OD iQ DC	Odyssey iQ Dust Cover
	OD iQ CHG US	Odyssey iQ Charging Kit, US
	OD iQ CHG G	Odyssey iQ Charging Kit, Global
	IC CB US	Icon Knee Charging Blade, US
	IC CB G	Icon Knee Charging Blade Kit, Global
	OD iQ RCK	Remote Charger Kit

SHELLTREAD

	PART NUMBER (LEFT/RIGHT)	DESCRIPTION	ADDITIONAL INFORMATION
	BE TLN x	BE TRN x	Shelltread (left or right), Neutral (indicate size)
			Increases durability and prevents premature wear of foot shell.

SOCKS

	PART NUMBER	DESCRIPTION	ADDITIONAL INFORMATION
	OD iQ SS 2325	CPI Sock, 23-25cm	Highly durable sock that acts as a noise barrier and protects the foot from debris.
	OD iQ SS 2628	CPI Sock, 26-28cm	
	OD iQ SS 2930	CPI Sock, 29-30cm	
	OD iQ SSK 2325	CPI Sock, 6-pack, 23-25cm	
	OD iQ SSK 2628	CPI Sock, 6-pack, 26-28cm	
	OD iQ SSK 2930	CPI Sock, 6-pack, 29-30cm	