



## SPS TF ÖSSUR LEG ORDER FORM PACKAGE

Complete the TF Össur Leg Order Form and email to [customerservice@spsco.com](mailto:customerservice@spsco.com).  
SPS will then create and send you a quote detailing all items in the Össur Leg.

[spsco.com](http://spsco.com) | 800.767.7776

### BILLING & SHIPPING INFORMATION

SPS Account #:	PO#	
<u>Ship To:</u>		
Company:	Contact:	
Address:	Phone:	Fax:
City/State/Zip:	Email:	

Completion of this order form with the most accurate and up-to-date information, including all patient information, measurements, and construction information enables us to provide the highest quality product for you and your patient.

### REQUIRED PATIENT INFORMATION

Patient ID:	Height:	ft.	in	Impact Level:	High	Med	Low	Amputation Level:	TF	KD	
	Weight:		lbs.	K-Level:	K1	K2	K3	K4	Left	Right	Bilat Order









### NOTES:

\*This form is available in electronic form you can fill out electronically and email or print.

# SPS TF ÖSSUR LEG ORDER FORM PACKAGE

## PRODUCTS

### 1: LINERS (check one then fill in desired sizes and options)


Locking Options		Seal-In Options					
	TF Locking					  	LINER SIZE
	LINER SIZE						Standard Conical

### 2: SOCKETS (check options)

No Socket	Connect TF Socket		
Limb circumference (mm):	Requested size:	Medium	Large
Limb length (mm):	Requested size:	Short	Standard

### 3: SUSPENSION & SOCKET ADAPTER

LOCKING OPTIONS					
Icelock 621 Ratchet & Extra Pin 	Icelock 651 Smooth & Extra Pin 	Icelock 600XM 	Icelock 211 Clutch & Extra Pin Icelock 214 4-Hole Clutch & Extra Pin 	Icelock 234 Lanyard 	
Works With					
Icelock 672 Male Pyramid Ti  Icelock 673 Male Pyramid Al	Icelock 674 4-Hole Adapter Al  Icelock 675 4-Hole Adapter Ti				
EXPULSION OPTIONS			UNITY ELEVATED VACUUM OPTIONS		
Icelock 544 Expulsion Kit 	Icelock 552 Expulsion Valve 	Unity Valve TF 	Icelock 544 Unity Plate 		
			Works With		
		Icelock 544 Socket Adapter 	3-Prong Socket Adapter 4-Prong Socket Adapter		

 Waterproof rated Össur products provide permanent protection against harmful ingress of fresh, salt or chlorinated water, even after submersion.

# SPS TF ÖSSUR LEG ORDER FORM PACKAGE

## PRODUCTS

### 4: ENDOSKELETAL COMPONENTS (check all that apply)



#### 4-HOLE ADAPTERS

- 4-Hole Male
- 4-Hole Female
- 4-Hole Male w/ Rotation
- 4-Hole Female w/ Rotation
- 4-Hole Male w/ 8mm Ax Offset
- 4-Hole Female Pyramid 20/20 Offset
- 4-Hole Female Pyramid 20/20 Offset, 10 degree, Right
- 4-Hole Female Pyramid 20/20 Offset, 10 degree, Left



#### INSERTS FOR PRONG

- Male Insert for Prong
- Female Insert for Prong
- LP Female Inserts for Prong



#### DOUBLE ADAPTERS

- Female Double Adapter 32mm
- Female Double Adapter 45mm
- Female Double Adapter 60mm
- Female Double Adapter 75mm



#### PYLONS

- Height Adjustable Standard
- Height Adjustable Long
- Pylon Al, 220lb Limit
- Female Short      Female Kit Short
- Pylon Ti Waterproof, 365lb Limit
- Female Long      Female Kit Long



#### ROTATION ADAPTER

- Knee Rotation Adapter (not included in package price)



#### TUBE CLAMPS

- Male Tube Clamp
- Female Tube Clamp
- 4-Hole Tube Clamp



#### SINGLE ADAPTERS

- Male Single Adapter Short
- Male Single Adapter Long
- Female Single Adapter Short
- Female Single Adapter Long

Other:

Össur may select different components to fit your dimensions if necessary. If you do not want Össur to make any changes, check this box:

### 5: KNEE (check one)

\*Locking knee cables will be installed on the anterior/lateral section of the check/definitive socket unless otherwise requested in the Notes section.



Locking Knee



Balance Knee OFM1 w/ IKF Adapter w/ Male Pyramid w/ Loop Adapter



Balance Knee OFM2



Total Knee 1900 A835300 Pyramid Top Extension Assist option



Total Knee 2000 A835300 Pyramid Top Extension Assist option



OHP3 Knee w/ IKF Adapter w/ Male Pyramid w/ Loop Adapter



Total Knee 2100 A835300 Pyramid Top Extension Assist option



OH7 Knee<sup>†</sup> w/ IKF Adapter<sup>†</sup> w/ Male Pyramid<sup>†</sup> w/ Loop Adapter<sup>†</sup>



Mauch Knee Mauch Knee Plus



Rheo Knee



Rheo Knee XC






















<sup>†</sup> Discontinued 31 Dec 2022


# SPS TF ÖSSUR LEG ORDER FORM PACKAGE

## PRODUCTS

**7: FOOT: WITH OR WITHOUT UNITY** (check one and fill out the options below)

No Foot **FOOT OPTIONS:** Foot Size: \_\_\_\_\_ cm. Category (1-9) \_\_\_\_\_ Left Right Foot Shell: Beige Brown

		Suggested L-Codes†			Suggested L-Codes†														
	Balance Foot S  with Torsion  with Unity	5972, 5986 5972, 5984, 5988 5781*		K2 Sensation with DP Flexion with Unity	5972, 5986 5972, 5968 5781*														
	Balance J with Unity	5976 5781*		Assure with Unity	5981 5781*														
	Pro-Flex ST 	5981, 5986		Vari-Flex with Unity	5981‡ 5781*														
	Talux**	5981‡, 5986‡		Pro-Flex LP Pro-Flex LP Torsion with Unity	5981, 5986 5981, 5984, 5988 5781*														
	Pro-Flex Pivot with Unity	See Össur reimbursement guide		Proprio Foot with Unity	5973‡ 5781*														
	Pro-Flex LP Align	5981, 5986, 5990		ReFlex Rotate with Unity	5987‡, 5984‡ 5781*														
	Pro-Flex XC  Pro-Flex XC Torsion with Unity	5981‡ 5987‡, 5984‡ 5781*		Torsion Shock Adapter with Male Pyramid Adapter with Female Pyramid Adapter	5984, 5988														
	ReFlex Shock with Unity	5987‡ 5781*																	
	Pro-Flex Modular 	5980‡, 5986‡	Desired build height (mm): <small>With foot cover on. No return or exchange after foot is cut.</small> <table border="1"> <thead> <tr> <th rowspan="2">Foot Size</th> <th colspan="2">Build Height</th> </tr> <tr> <th>Min</th> <th>Max</th> </tr> </thead> <tbody> <tr> <td>22-24cm</td> <td>157mm</td> <td>313mm</td> </tr> <tr> <td>25-27cm</td> <td>168mm</td> <td>327mm</td> </tr> <tr> <td>28-30cm</td> <td>176mm</td> <td>340mm</td> </tr> </tbody> </table>			Foot Size	Build Height		Min	Max	22-24cm	157mm	313mm	25-27cm	168mm	327mm	28-30cm	176mm	340mm
Foot Size	Build Height																		
	Min	Max																	
22-24cm	157mm	313mm																	
25-27cm	168mm	327mm																	
28-30cm	176mm	340mm																	

 Waterproof rated Össur products provide permanent protection against harmful ingress of fresh, salt or chlorinated water, even after submersion.

\*Unity code added to other codes, e.g. Proprio Foot with Unity: 5973, 5781.

\*\*Discontinued 31 Dec 2022

† The responsibility for accurate coding lies solely with the healthcare provider, Össur assumes no responsibility or liability for the provider's coding decisions.

‡ PDAC approved.

CLEAR FORM

4 of 4