

Patient Aids

- Skin Care 532-536
- Left Foot Gas Pedal 536
- Shoe Horns 536
- Xerosox Prosthetic Protectors 537
- Disposable Shorts 537
- Donning Devices 538-539
- Canes, Crutches and Accessories 540-545
- Walkers and Accessories 545-550
- Commodes 550

ALPS - Silicone Skin Lotion

- 100% silicone skin lotion helps to smooth dry, sensitive skin
- Hypoallergenic, with no perfumes or dyes
- Leaves no oily feeling and can be used sparingly to achieve relief for dry, chapped skin conditions

Description	Product No.
4 oz Bottle	SL108
32 oz Bottle	SL108/32

**SL108****Cleanser -16 oz bottle (not pictured)**

- Non-alkaline, soap-free cleanser which can be used for either hand or machine washing of prosthetic liners and sleeves
- Patients with sensitive skin may also find that this soothing cleansing agent is ideal for washing skin areas that come into contact with the prosthesis

Product No. **PD595**
**PA604****Lubricant (16 oz bottle)**

- Designed for donning silicone prosthetic liners
- Eliminates the sticky feel and allows the patient to roll up the liner or sheath easily
- Active Ingredient: Silicone
- 70% Isopropyl Alcohol

Product No. **PA604**
**PO840****Ointment (4 oz tube)**

- Vitamin A and D formulation that soothes, protects and helps to promote healing and reduce the pain of skin rash or superficial skin injuries

Product No. **PO840**
**ATP008****Antiperspirant**

- Designed for prosthetic patients to help manage the continuing problem of perspiration
- 4 oz spray

Product No. **ATP008**
Fitting Lotion

- 4 oz bottle

Product No. **FL801**
**FL801**

Iceross® Clean & Simple Set 1

- Includes one bottle of soap, moisturizing lotion and lubricant spray

Product No. **S-010001**

Ossur - Iceross® Clean & Simple Soap

- Iceross Clean & Simple soap is a mild liquid soap for daily hygienic treatment of the skin and the silicone liner
- Free from fragrances and dyes, it is designed for sensitive skin and contains natural ingredients that nourish and soften the skin
- It has a pH of 5.5 and is hypo allergenic (250 ml/8.8 oz)

Product No. **S-210250**

Iceross® Clean & Simple Moisturizing Lotion

- A non-greasy restorative lotion that nourishes and softens the skin
- Free from fragrances and dyes
- It has a pH of 5.5 and is hypo allergenic
- 250 ml/8.8 oz

Product No. **S-410250**

Iceross® Clean & Simple Lubricant Spray

- Aids the donning and doffing of Iceross Seal-In X5 Liners and Iceross Original Liner without fabric cover
- 250 ml/8.8 oz

Product No. **S-610250**

Liquid Body Powder

- Effective, safe, non-allergenic
- Controls chaffing, perspiration and odor while soothing areas of friction
- 8 oz bottle

Product No. **57751**

SPS - Liquid Silk Body Powder

- Also known as Dry-Lite™
- 18 oz bottle

Product No. **LS100**

**57751****LS100**

Cornstarch Baby Powder

- Corn Starch, 15 oz

Product No. **199386**



Talc Baby Powder

- 14 oz

Product No. **51393**



B.F.I. First Aid Powder

- First aid to help reduce the risk of skin infection in minor cuts, scrapes and burns
- Antiseptic dusting powder; soothing, protective, absorbent
- Active Ingredients: Bismuth-Formic-Iodide, 16%
- 1.25 oz

Product No. **86214**



Cold Spray

- Relieves hot conditions
- Lowers surface temperature
- 6 oz can

Product No. **C7454**



It Stays Adhesive

- Gentle roll-on body adhesive that washes off with water
- 2 oz bottles; 12 bottles per case

Product No.

6



SPS - Odor Eliminator

- All natural product that kills odor causing bacteria; does not mask the odor, but destroys the bacteria (safely, and without chemicals) that cause offensive foot, body or prosthesis odor
- The regeneration of bacteria can take weeks, which is why Odor Eliminator lasts, even after showering; after using Odor Eliminator for a few months, the patient will find the need for re-application less frequent
- 1/2 oz (15 g) jar

Product No.

OE100



Patient Aids

Diabetic Basics - Antimicrobial Instant Foot Sanitizer

- Colorless, water-based spray product sanitizes, moisturizes with skin softeners and emollients, and helps prevent infection, killing 99.9% of germs
- An option for amputee use
- For limb hygiene
- 8 oz spray bottle

Product No.

BAS-29001



Moleskin

Description	Size	Product No.
12" x 5 yards	White	58369
12" x 5 yards	Flesh	58364
* contains latex		




AmpuBalm

- 3³/₄ oz jar; 24 jars per case
- Protects abrasions from further chafing or injury

Product No.

582001



Adaptskin

- Protection for amputees who experience skin breakdown, heat rash, abrasions, chafing, dry skin and other residual limb irritations
- Hypoallergenic, penetrates to form a waterproof barrier
- Compatible with mineral oil, silicone and urethane liners
- Adaptskin 50™ is Medium viscosity for more general applications and daily use
- Adaptskin 90™ is a thicker, more concentrated product for high abrasion areas or severe skin damage

Description	Product No.
Adaptskin 50	AS5004J
Adaptskin 90	AS9004J

**2nd Skin™ Dressing**

- Increases healing speed
- Prevents wounds from drying out
- Reduces risk of scarring
- Latex free

Description	Product No.
3" circles, 48 per jar, non-sterile	10-603
2 pads, 3" x 6 1/2", non-sterile	47-209
4 pads, 2" x 3", sterile moist burn pads	47-019

**Left Foot Gas Pedal**

- Designed for people with limited or no use of their right foot
- Quick release gas pedal can be installed in any vehicle with an automatic transmission
- No adjustments necessary after installation
- Pedal stop will not allow right foot to inadvertently apply throttle

Product No. **D7M**

**Pro Metal Shoe Horn**

- Dozen per box

Product No. **583120**



Xerosox Prosthetic Protectors

- Designed as a cast/bandage/prosthetic protector to create a watertight vacuum seal
- Eliminates air from the cast or bandage by simply squeezing the built-in hand pump
- Removal is easy: Gently lift the cuff to let the air back in
- Non-skid grid provides secure footing

Xerosox Prosthetic Leg Protector			
Description	Circumference	Length	Product No.
Small	11-15.5"	32"	PL-14
Medium	15-17"	35"	PL-16
Large	16.5" and up	37"	PL-18
Xerosox Pump Leg Protector			
Description	Circumference	Length	Product No.
Small Half Leg	10-13"	21"	HL-13
Large Half Leg	13" and up	23.5"	HL-15
X-Small Full Leg	7.5-11"	19"	FL-12
Small Full Leg	14-16.5"	29"	FL-14
Medium Full Leg	16.5-21"	33"	FL-16
Large Full Leg	21" & up	37"	FL-18
Xerosox Pump Arm Protector			
Description	Circumference	Length	Product No.
Small Half Arm	7.5-9.5"	16.5"	HA-13
Large Half Arm	9.5" and up	18.5"	HA-15
X-Small Full Arm	6-7"	15"	FA-12
Small Full Arm	7-9"	21"	FA-14
Medium Full Arm	9-10"	26"	FA-16
Large Full Arm	10" & up	28.5"	FA-18



Disposable Boxer Shorts

- 25 per bag

Size	Product No.
Medium	960400
Large	960401
X-Large	960402
XX-Large	960404



Disposable Exam Shorts

- Available in multi-layer material for comfort and fluid resistance
- Elastic waist provides a snug fit
- Color: Blue
- Latex free
- 30 per case

1 case

NON27209L



Zip Whiz Easy-Proth[®]

- Prosthetic donning sleeve; Double length design prevents friction to the skin caused by excessive force

Size	Color	Proximal Circumference	Product No.
X-Small	Yellow	Up to 15 ¹ / ₂ "	5-EU-300-100XS
Small	Red	15 ¹ / ₂ " to 19 ¹ / ₂ "	5-EU-300-100S
Medium	Blue	19 ¹ / ₂ " to 23 ¹ / ₂ "	5-EU-300-100M
Large	Green	23 ¹ / ₂ " to 27 ¹ / ₂ "	5-EU-300-100L
X-Large	Lt. Blue	27 ¹ / ₂ " to 31 ¹ / ₂ "	5-EU-300-100XL

**Knit-Rite - Power Tube Prosthetic Pull Sock[™]**

- A low friction donning sock that is simple, effective and economical
- Open on both ends and fashioned as a straight, non-tapered tube

Size	Max. Circumference	Product No.
X-Small	16"	3M-PT08-XS
Small	20"	3M-PT10-SM
Medium	24"	3M-PT12-MD
Large	28"	3M-PT14-LG
X-Large	32"	3M-PT16-XL
XX-Large	36"	3M-PT18-2X

**Power Pull Donning Sock**

- The single piece integral design provides a more direct and effective pull transfer and eliminates stress points for a more durable product; 100% Nylon


Size	Max. Circumference	Product No.
Size A	less than 8"	3MPP04AS
Size B	8" - 12"	3MPP06BS
X-Small	12" - 16"	3MPP08XS
Small	16" - 20"	3MPP10SM
Medium	20" - 24"	3MPP12MD
Large	24" - 28"	3MPP14LG
X-Large	28" - 32"	3MPP16XL
XX-Large	32" - 36"	3MPP182X



Ultra Low Friction Sheath™ (ULFS)

- For push-in donning problems - apply over socks or skin, leaving extra length in outer layer

Width	Product No.			
	14" Single Layer with Hole	14" Double Layer	21" Single Layer with Hole	21" Double Layer
8"	2HU11408	2HUL1408	—	—
10"	2HU11410	2HUL1410	2HU12110	2HUL2110
12"	2HU11412	2HUL1412	2HU12112	2HUL2112
14"	2HU11414	2HUL1414	2HU12114	2HUL2114
16"	2HU11416	2HUL1416	2HU12116	2HUL2116
18"	—	—	2HU12118	2HUL2118

**Easy-Slide Donning Sleeve**

- Ultimate donning aid for open toe compression stockings
- Reduces physical strain for both caregiver and user
- Increase in patient's independence and therapy compliance
- Stockings retain therapeutic graduated compression longer
- Easy to use; durable

Size	Color	Shoe Size		Product No.
		Womens	Mens	
X-Small	Pink	<3	<2	E01002
Small	Light Blue	3-6	2-4 ¹ / ₂	E01102
Medium	White	6 ¹ / ₂ -9 ¹ / ₂	5-8	E01202
Large	Dark Blue	10- 13	8 ¹ / ₂ -11 ¹ / ₂	E01302
X-Large	Purple	≥13 ¹ / ₂	≥12	E01352

Easy-Slide Caran

- Same material as Easy-Slide, but designed for use with closed toe compression stockings
- Ingenious velcro release

Size	Color	Circ. of Widest Part of Foot	Product No.
Small	Light Blue	<8"	E01752
Medium	White	8"-10"	E01753
Large	Dark Blue	>10"	E01754

Easy-Slide Arm Donning Sleeve

Size	Product No.
Medium	E01402
Large	E01425

Easy-Off Stocking Removal Aid

Size	Product No.
Small	E01802
Medium	E01902
Large	E02002



Offset Handle Models

- Adjustable 31-39" (79-99 cm)
- Gray vinyl grip, silver shaft

Product No. **5085**

Curved Handle Models

- Adjustable 31-39" (79-99 cm)
- Gray vinyl grip, silver shaft

Product No. **5086**

Unipoise Single Point Cane

- 7/8" Offset Handle; Adjustable from 31"-40"
- 300 lb weight limit

Product No. **G05355**

Low Profile Quad Cane

- Anodized 7/8" (2.2 cm) silver aluminum shaft
- Low profile steel base, chrome finish, narrow or wide
- Offset handle, molded gray vinyl hypalon foam hand grip
- Spring button height adjustment
- Reinforced rubber tips
- Height adjustment 29"-38" (74-96 cm)

Description	Product No.
Narrow Base, 6" x 8" (15 x 20 cm) - Each	G05345S
Wide Base, 8" x 12" (20 x 30 cm)	5096

Medline - Offset Handle Aluminum Cane

- Contoured handle with foam handgrip comes with a wrist strap to provide comfort and a secure grip
- Handle height adjusts in one-inch increment from 29"-38" and is secured with a locking ring for added safety
- Offset handle aluminum cane with wrist strap
- Black
- 250 lb weight limit

Product No. **MDS86420**



5085

5086



G05355



G05345S



MDS86420

Offset Handle Cane - Bariatric

- All-steel construction
- Solid brass, double push button adjustment for extra stability
- Handle height adjusts in one-inch increment from 29"-38" and is secured with a locking ring for added safety
- PVC handgrip
- Double-tubing insert for added strength and safety
- Black
- 500 lb weight limit

Product No.

MDS86420XW



MDS86420XW

Large Base Quad Cane

- Double-plated base with extremely strong welds
- Solid brass, push button adjustment for extra stability
- Handle height adjusts in one-inch increment from 29"-38" and is secured with a locking ring for added safety
- All-steel construction
- Black shaft
- 500 lb weight limit
- 8"x12" black footprint

Product No.

MDS86228XW



MDS86228XW

Aluminum Large Base Quad Cane

- Dual push-button adjustment mechanism with locking ring allows for quick, secure fitting and to be ambidextrous
- Height adjusts from 30" to 39"
- The non-skid rubber tips with metal inserts ensure a secure grip and long lasting performance
- No assembly required
- Aluminum construction
- 250 lb weight limit
- Black

Product No.

MDS86228



MDS86228

Aluminum Small Base Quad Cane (not pictured)

- Dual push-button adjustment mechanism with locking ring allows for quick, secure fitting and to be ambidextrous
- Height adjusts from 30" to 39"
- The non-skid rubber tips with metal inserts ensure a secure grip and long lasting performance
- Aluminum construction
- No assembly required
- 250 lb weight limit
- Black

Product No.

MDS86222W

Forearm Crutches

- 1" (2.5 cm) anodized aluminum tubing
- Vinyl coated cuffs with 180° pivot and nylon pivot block
- Gray handgrips on welded and coated steel handles
- Spring button adjustments at both legs and forearm cuffs
- Metal reinforced rubber tips

Description	Product No.
Adult 5'2" to 6'2"	5090
Junior 4'2" to 5'2"	5090-J



5090

Medline - Aluminum Forearm Crutches, Adult

- Built tough for long-lasting performance
- Recommended for long-term disability
- Contoured vinyl handgrips provide comfort and long wear
- 250 lbs weight capacity

Description	Height	Product No.
Tall	5'10" to 6'6"	MDS805160
Average	5'0" to 6'2"	MDS805161



MDS805161

Quick-Fit Crutches

- Bolt-adjusted handle and push-button adjusted foot pieces with fitting scale allow for accurate fitting
- Aluminum; Latex-free rubber components
- 300 lb weight capacity

Description	Height	Product No.
Tall	5'10" to 6'6"	G50314-8
Adult	5'2" to 5'10"	G51314-8
Youth	4'6" to 5'2"	G52314-8
Child	3'9" to 4'7"	G53314-8



G51314-8

Red Dot Aluminum Crutches

- Clip-lock adjusted handle and push button adjusted height allows for 10 second adjustments
- 1" adjustment increments on both handle and height provide a precise fit which reduces muscle fatigue
- Latex-free rubber components
- 300 lb weight capacity

Product No.	G91-214-8
-------------	------------------




G91-214-8

Heavy Duty Tall Adult Crutches

- Oversized wing nut and push-button for quick height adjustment
- Metal ring helps prevent excessive wear and ensures maximum durability
- Aluminum
- Latex-free rubber components
- 650 lb weight capacity

Description	Product No.
Tall Adult	MDS80334XW
Adult	MDS80335XW



Bariatric Crutches

- Heavy-duty steel construction provides durability
- Push-button adjusted foot pieces with fitting scale allow for accurate, fast fitting
- 1" adjustment increments on both handle and height provide a precise fit which reduces muscle fatigue
- Supports up to 550 lbs

Description	Height	Product No.
Adult	5'1" to 5'9"	G61314B
Adult Tall	5'10" to 6'6"	G60314B



Wood Adjustable Crutches

- Made from strong weather resistant laminated hardwood
- Preshaped bows allow for quick and easy adjustment
- Includes tip, grip and underarm cushions

Description	Height	Product No.
Adult	5'4" to 6'5"	5073-S
Medium	4'1" to 5'7"	5074-S
Junior	4'5" to 5'	5075-S
Child	3'8" to 4'4"	5076-S



Super Crutch Tips

- Cushion base absorbs shock and provides extra comfort
- For sure footed traction the extra wide vacuum seal rings grip when weight is applied and release when weight is removed
- Fits: $\frac{7}{8}$ " (2.2 cm) shafts
- Base Diameter: $2\frac{1}{16}$ " (5.2 cm); Overall Height: $3\frac{1}{16}$ " (7.8 cm)
- Color: Gray
- Sold by the pair

Product No.

G00842



Standard Crutch Tips

- 100% natural rubber
- Metal disc molded into tip prevents foot piece cut through
- Fits: $\frac{7}{8}$ " (1.9 cm) shafts
- Base Diameter: $1\frac{3}{4}$ " (4.4 cm); Overall Height: $2\frac{1}{2}$ " (6.3 cm)
- Color: Gray
- Sold by the pair

Product No.

G00502



Close Style Crutch Handgrips

- Provides sturdy support that alleviate blistered hands, cramped fingers, wrist tension and pressure on sensitive nerves
- For crutches with removable handles
- Neoprene added for strength
- Color: Tan

Product No.

5072



Open Style Crutch Handgrips

- Provides sturdy support that alleviate blistered hands, cramped fingers, wrist tension and pressure on sensitive nerves
- For fixed handle crutches, special glue packed with each pair for added security
- Neoprene added for strength
- Color: Tan

Product No.

5071



Crutch Pads (not pictured)

- High density reinforced rubber with molded-in fabric keeps the cushion from splitting or tearing
- Extra thick at the top for maximum comfort
- Natural rubber compound with corrugated cell construction provides a snug, tensioned fit to the arm piece
- Ribbed design adds additional comfort
- Fits all adult and medium crutches
- Color: Tan

Product No. **5070**

Crutch Tips

- Shock absorbing concave base

Size	Color	Product No.
Small	Tan	6099-T
Medium	Tan	6100-T
Large	Tan	6101-T



AliMed® Aluminum Four-Wheel Rollator

- Lightweight and versatile (13 lbs)
- Easily folds for quick storage
- Large nylon pouch under seat to store necessary items in transit
- Locking loop brakes
- Weight capacity: 300 lbs
- Adjustable handle height (33"-38")
- Dimensions: 23" W x 23" L
- Seat: 13" D x 14" W x 22" H

Product No. **77016**



AliMed® One-Button Foldaway Walker

- Polished aluminum frame
- Pinch-proof release button
- Folds for easy storage and transportation
- High cross brace design allows walker use over commodes
- Cushioned handgrips
- Weight capacity: 350 lbs
- Height adjusts in 1" increments
- Dimensions: 31" - 38" H x 18" D (3 1/2" folded)

Product No. **74261**



Medline - Strider™ Prescription Walker

- Handles adjust to an infinite number of positions to provide stability and posture
- Extra wide base provides stability, safety and balance
- Large 5" wheel foot pieces will easily roll over irregular surfaces
- Options include platform attachment, plus horizontal or vertical handles which may be used in any combination or position
- All adjustments accomplished without tools
- Can be used as a trailing walker or in the traditional forward position



Regular

- Approximate Patient Height: 4'8" to 6'4" (142-193 cm)
- Color: Black
- Weight: 9 lbs (3.6 kg)
- Height Adjustment: 28-40" (71-102 cm)
- Base Width: 27"

Product No.

G07780

Children's Folding Walker

- Height-adjustable "grow with you" design
- Adjustable from 20.5"-24.5"

Product No.

G07750

Extra Wide Folding Walker

- Two button folding allows independent operation of each side so walker fits through narrow spaces
- Easy to use folding mechanism is designed for operation with palm, finger or arm
- Lightweight, polished 1" aluminum tubing with dual steel cross bars supports up to 400 lbs
- Contoured vinyl handgrips

Description	Height	Product No.
Adult	5'4" to 6'2"	G07767
Youth	4'7" to 5'3"	G07768



Easy Care Folding Walker

- Dual folding mechanism is easy to operate with palms or fingers
- Ideal for those with limited dexterity
- Lightweight aluminum frame with steel cross brace provides durability
- Supports up to 300 lbs

Description	Max Height	Product No.
Adult	6.2"	G30755P
Youth	5.3"	G30756P

**Two Button Folding Adult Walker**

- Height adjustable 32.5" to 39.5"
- Each side operated independently to provide easy movement through narrow spaces and greater stability when standing up
- Contoured textured latex free PVC handgrips provide additional comfort
- 5" wheels
- Width inside handgrips: 18"; width inside base legs: 21"; depth opened at base: 16"
- Supports up to 300 lbs

Description	Product No.
Walker	MDS86410W54
Replacement Parts	
Extension Leg	MDS86615E
Glide Caps	MDS86615C
5" Wheel (2/pack)	G07722G

**Guardian Adult Walker**

- Easy Care Walkers; dual folding mechanism is easy to operate with palms or fingers
- Lightweight aluminum frame with steel cross brace for durability
- Ideal for those with limited dexterity
- Folds to only 4" for easy transport and storage
- Weight limit: 300 lbs
- 5" wheels
- Latex free

Size	Max Height	Product No.
Adult	6.2"	G30757W
Youth	5.3"	G30758W



Extra Wide Two Button Walker

- Features easy-to-use push-button mechanisms that may be operating by the fingers, palm or side of hand
- Polished aluminum frame
- Folds for easy storage and transportation
- Dual 5" wheels
- Weight limit: 500 lbs
- Height adjust from 30" to 43¹/₂" in one-inch increments
- Dimensions: Width inside handgrips: 20¹/₂"; width inside base legs: 28"; depth opened at base: 19"

Description	Product No.
Walker	MDS86410XWW
Replacement Parts	
Extension Leg	MDS86615E
Glide Caps	MDS86615C
8 Hole Footpiece set, 5" Wheel	G07722-8B

**Adult Sidestepper Cane/Hemi-Walker**

- Lightweight frame is easy to lift and maneuver
- Dual grip offers extra assistance
- Provides support and stability to those with limited or no dexterity in one arm or hand
- Wide base provides more stability than a quad cane

Product No.	G07770
-------------	---------------

**Adult Hemi Walker (not pictured)**

- Lightweight aluminum frame with deluxe rubber tips for excellent traction
- Adjustable from 30" to 34"
- Folds to 3" for convenient storage
- Depth opened at base: 14"
- Weight limit: 250 lbs

Product No.	MDS86620
-------------	-----------------

Adult Bariatric Folding Walker (not pictured)

- Reinforced aluminum frame with dual steel cross brace
- Extra wide frames with up to 22" between handgrips to accommodate larger users
- Weight limit 650 lbs

Description	Max Height	Product No.
Walker	6.6"	G30754B
Replacement Parts		
8 Hole Footpiece set, 5" Wheel		G07722-8B

Walker Attachment Platform

- Comfort padded gray forearm cuff
- Universally adjustable latex free hand grip & forearm cuff pads
- Attachment requires two cross braces on the button in order to work; it will not work on 1-button or youth walker
- 300 lb weight capacity

Product No.

MDS86615P

Patient Aids

Adult Crutch Attachment Platform

- Platform attaches to the top of the crutch
- Acts as a padded armrest with handgrip
- Useful for those who cannot straighten arm to grip the hand-hold of the crutch

Product No.

G07706**Leg Extension for Walker**

- Includes one pair of glide caps
- Extension legs 14.5" with 5 holes
- Used on MDS86410 series

Product No.

MDS86615E**Glide Caps**

- Includes one pair of glide caps
- Used on MDS86410W54 and MDS86410XWW Two Button Walkers

Product No.

MDS86615C

Footpiece Extension Set

- Front wheel foot pieces allow you to easily maneuver over surfaces which provide difficulty for walker users
- It also provides a more natural walking stride
- Available in a choice of 3" and 5" inline wheels on 360 degree casters which all you to change directions with ease

Description	Product No.
3"	G07731G
5"	G07722G
Swivel Set	G07820-8

**Extra Wide Commode**

- This bariatric commode has a back rail that removes quickly and easily for convenient storage and use
- Features a durable, plastic seat with lid plus armrests for additional comfort
- Includes a 10-quart capacity commode bucket with seat cover and splash shield
- Height adjusts from 19.5" to 23.5" in one inch increments
- Non-skid, non-marring rubber tips provide extra safety
- 28" base width, 20.5 overall base depth; 24" width between arms; 14" seat width; 15" overall seat depth; 20.5" seat height
- Weight limit: 650 lbs

Product No. **MDS89664XW**

**3-in-1 Steel Commode (not pictured)**

Product No. **MDS89664H**

Universal Fit Commode Bucket

- Constructed of high-impact, white polyurethane
- Universal; fits all Medline commodes
- Complete with handle and lid
- Diameter of bucket is 11"; 12 quart/2 gallon capacity

Product No. **MDS80306B**

