

Sock Measurement Instructions

- Sizing and Measurement Custom Socks 312-313

Prosthetic Socks

- Comfort Products 314-318
- Knit-Rite 318-328
- SPS 328

Specialty Prosthetic Socks

- Liner-Liner Prosthetic Sock 329
- Hip Disarticulation Prosthetic Sock 329
- AK and BK Post-Op Prosthetic Socks 330

Casting/Fitting Socks

- Comfort Products™ 330-331
- Knit-Rite® 331
- RX Textiles 332
- SPS 333

Filler/Spacer Socks

- Knit-Rite® 334
- RX Textiles 335

Gel Socks

- Comfort Products 336
- Silipos® 336-337
- SPS 337

Sheaths

- ALPS 338
- Comfort Products 338
- Knit-Rite 339-340
- Silipos 341
- Truform® 342

Shrinkers

- Bell-Horn® 342
- Comfort Products 342-343
- Juzo® 343-344
- Knit-Rite 345-349
- Rx Textiles 350
- Supreme Medical 350
- Truform 351-352

Sock Measurement Instructions

- To provide a proper fitting prosthetic sock for your patient, it is important that a correct measurement be made
- There are three basic measurements that are necessary for a proper fitting prosthetic sock:
 - length
 - top width
 - toe width

A - B Measurement (top width)


- Should be taken from top left corner of sock lying flat to top right corner and correspond to $\frac{1}{2}$ the circumference measurement at the point where the top of the sock should be on the leg or arm

C - D Measurement (length)

- Should be taken from top of sock lying flat to bottom of toe and should correspond to the length from distal to where top of sock should be on the leg or arm

E - F Measurement (toe width)

- Should be taken from left side of sock lying flat, perpendicular to toe fashioning, to right side of sock and should correspond to $\frac{1}{2}$ the circumference 2" up from the distal end of the limb



Seam of sock is exaggerated for this illustration

Note: Remember – a proper measurement means a correct fitting prosthetic sock for your patient.

Guidelines for Sizing and Measuring Custom Socks

Prosthetic Sock Sizes

Nearly 90% of the amputees who use socks can be correctly fitted with a stock size sock because the knitted construction has a stretch which allows a fitting range. Most prosthetists can judge the needed size by working with the limb and the cast. Nevertheless, some measurement guidelines can be helpful. First, measure the limb circumference at the top of the socket; second, measure the toe limb circumference two (2) inches from the distal end; third, measure from the distal end to the top of the socket. The length of the sock should be 2 to 4 inches more than this measurement depending on the type of fitting and suspension preference, the fabric (Wool and A Plus will shrink more than cotton) and the amount of tautness the patient is likely to achieve in donning.

Socks are sized in flat measurements; i.e., half the circumference measurement (less 5% to give a snug fit) for top and toe – then matching those figures to the closest comparable size measurements in all the charts – see previous page. Custom socks need to be ordered only when patient's measurements do not correspond with any of the coded stock sizes or when a sock needs an addition or adjustment, i.e., split for hip joint, suspension system hole, elastic at waist or leg, fly, etc.

Hip Joint/Pelvic Belt Prosthetic Socks

An order for socks with a split for hip joint should state the size of the sock and whether or not a tie or just the edges split is desired. Slits are three (3) inches deep unless otherwise specified.

Ischial Containment/High Lateral Brim

To eliminate excessive rollover, medial cuts for sockets with high lateral wall can be made when requested. Specify the depth of the cut desired.

Continued on next page

Guidelines for Sizing and Measuring Custom Socks - Continued

Irregular Shapes

For prosthetic socks which are not evenly graduated from the top to the toe, take circumference measurements at the irregular points, lengths from the distal end to each of those points and determine total length. Sending a drawing is very helpful. State fabric and ply choice.

Fitting Socks and Re-Orders

Please state customer's name on all orders for custom socks. A "fitting" sock will normally be sent. It is suggested that only one sock be ordered until the right size and material can be established, then the quality can be specified. Records for each customer are kept and reorders can duplicate the original or alterations can be made, as needed, usually without another pattern.

1. Symes and Chopart

Socks have a bulbous end. Standard sizes available. If custom size is needed, prepare a profile drawing of the limb, giving circumference measurements at the top, ankle, and ball. Note on drawing if limb is symmetrical or asymmetrical. If asymmetrical, record a circumference measurement of the partial foot extension.

2. Knee Disarticulation (box toe)

State the length desired, the top circumference at the perineum and the toe circumference one (1) inch from the end of the limb. Specify box toe.

3. Large Knee and Thigh

Send a drawing of the limb, front to back view. Specify length of sock with allowance for turn down. State length from the end to the beginning of the bulge below the knee and from the end of the largest part of the knee. State circumferences at top, knee, below condyles and two (2) inches from the end of the limb.

4. Flexed Knee

Prepare a profile drawing. Specify the length of the sock. Specify circumference at top, above knee, around knee, and at the end of the limb. If flexion angle is desired, please specify.

5. Tailoring

When socks need to be very irregular, an outline pattern with measurements every two inches can be sent or an exact pattern can be duplicated. A single flat seam will be used in construction with all the edges being finished. Tailoring of the knitted material can produce socks that are impossible to fashion on the knitting machines, especially when rounding is desired. If knitted socks or tailored socks are preferred, specify by using the "K" or "T" code.

6. Above Elbow Arm Socks

Prepare a drawing of the front view. Note the total length and the length from the axilla to the distal end of the limb. Take circumference measurements one (1) inch from the distal end and at the axilla. Specify a flat measurement from the upper arm to the point on the chest where the sock should come.

7. Hip Disarticulation and Hemi-Pelvectomy

Prepare a drawing as illustrated. Specify length from groin to waist. Hip disarticulation socks are made four (4) inches longer than this measurement to allow turnover, unless otherwise specified. Give circumference measurements for waist, hips and sound leg (take at angle). If fly or colostomy hole is desired, specify which leg has been amputated. Hip socks are also made with leg extensions and/or elastic at the waist and leg.

WHEN ORDERING please provide the following information for special and irregular prosthetic socks

- | | |
|--------------------|--------------------------|
| 1. Patient's name | 5. Toe width, flat |
| 2. Quantity | 6. Material and weight |
| 3. Length | 7. Shipping instructions |
| 4. Top width, flat | |

Irregular or special size: send complete information and sketch of stump, including length and circumference measurements.

Comfort Products - Regal Acrylic Prosthetic Sock

- Made of Acrylic and Lycra®
- Soft stretchable sock providing the patient with the latest advancement in fiber technology and knitting methods
- With the combination of fibers it gives a degree of stretch and conformance for the residual limb not possible with other fibers while wicking away moisture to ensure greater comfort
- Can be machine washed and dried in the dryer on delicate cycles
- Band color indicates ply: 2=beige, 3=yellow, 4=pink, 5=green, 6=gray, 7=gold
- Replace ? with ply
- Add 3 to the end of product number for sock with distal hole



| Regal™ Acrylic Sock (3, 5 and 6 Ply are Standard – 2, 4 and 7 Ply are Special Order) | | | | | | | | |
|--|-------------------|--------|-------|-------|---------|------|---------|-------------------|
| Size | Fits Current Size | Length | | Top | | Toe | | Product No. |
| | | in. | cm | in. | cm | in. | cm | |
| Narrow X-Short | 00, 0, A0, A1 | 8-10 | 20-25 | 4-6 | 10-15 | 3-3½ | 7.6-9 | PSRNAXS? |
| Narrow Short | 00, 0, A0, A1 | 10-12 | 25-30 | 4-6 | 10-15 | 3-3½ | 7.6-9 | PSRNASH? |
| Narrow Regular | 00, 0, A0, A1 | 14-16 | 35-40 | 4-6 | 10-15 | 3-3½ | 7.6-9 | PSRNARG? |
| Narrow Long | 00, 0, A0, A1 | 18-20 | 45-50 | 4-6 | 10-15 | 3-3½ | 7.6-9 | PSRNALG? |
| Narrow X-Long | 00, 0, A0, A1 | 22-24 | 55-60 | 4-6 | 10-15 | 3-3½ | 7.6-9 | PSRNAXL? |
| Medium X-Short | 1, 2 | 8-10 | 20-25 | 7-8 | 17.5-20 | 4-5 | 10-13 | PSRMDXS? |
| Medium Short | 1, 2 | 10-12 | 25-30 | 7-8 | 17.5-20 | 4-5 | 10-13 | PSRMDSH? |
| Medium Regular | 1, 2 | 14-16 | 35-40 | 7-8 | 17.5-20 | 4-5 | 10-13 | PSRMDRG? |
| Medium Long | 1, 2 | 18-20 | 45-50 | 7-8 | 17.5-20 | 4-5 | 10-13 | PSRMDLG? |
| Medium X-Long | 1, 2 | 22-24 | 55-60 | 7-8 | 17.5-20 | 4-5 | 10-13 | PSRMDXL? |
| Medium XX-Long | 1, 2 | 26-28 | 65-70 | 7-8 | 17.5-20 | 4-5 | 10-13 | PSRMDXXL? |
| Wide X-Short | 3, 4 | 8-10 | 20-25 | 9-10 | 23-25.5 | 5½-6 | 13.5-15 | PSRWDXS? |
| Wide Short | 3, 4 | 10-12 | 25-30 | 9-10 | 23-25.5 | 5½-6 | 13.5-15 | PSRWDSH? |
| Wide Regular | 3, 4 | 14-16 | 35-40 | 9-10 | 23-25.5 | 5½-6 | 13.5-15 | PSRWDRG? |
| Wide Long | 3, 4 | 18-20 | 45-50 | 9-10 | 23-25.5 | 5½-6 | 13.5-15 | PSRWDLG? |
| Wide X-Long | 3, 4 | 22-24 | 55-60 | 9-10 | 23-25.5 | 5½-6 | 13.5-15 | PSRWDXL? |
| Wide XX-Long | 3, 4 | 26-28 | 65-70 | 9-10 | 23-25.5 | 5½-6 | 13.5-15 | PSRWDXXL? |
| Extra Wide Short | 5, 6 | 10-12 | 25-30 | 11-12 | 28-30.5 | 7-8 | 17.5-20 | PSRXWD SH? |
| Extra Wide Regular | 5, 6 | 14-16 | 35-40 | 11-12 | 28-30.5 | 7-8 | 17.5-20 | PSRXWDRG? |
| Extra Wide Long | 5, 6 | 18-20 | 45-50 | 11-12 | 28-30.5 | 7-8 | 17.5-20 | PSRXWDLG? |
| Extra Wide X-Long | 5, 6 | 22-24 | 55-60 | 11-12 | 28-30.5 | 7-8 | 17.5-20 | PSRXWDXL? |

Comfort AG Sock

- High performance alternative to the Regal sock
- CoolMax fibers plaited with X-Static silver
- Available in Standard 3, 5 and 6 ply, Special Order 2, 4, and 7 ply
- Antimicrobial protection
- Wicks moisture away from skin
- Replace ? with Ply
- Add 3 to the end of product number for sock with distal hole



| Comfort AG Prosthetic Sock (3, 5 and 6 Ply are Standard - 2, 4 and 7 Ply are Special Order) | | | | | | | | |
|---|-------------------|--------|-------|-------|---------|------|---------|--------------------|
| Size | Fits Current Size | Length | | Top | | Toe | | Product No. |
| | | in. | cm | in. | cm | in. | cm | |
| Narrow X-Short | 00, 0, A0, A1 | 8-10 | 20-25 | 4-6 | 10-15 | 3-3½ | 7.6-9 | PSRAGNAXS? |
| Narrow Short | 00, 0, A0, A1 | 10-12 | 25-30 | 4-6 | 10-15 | 3-3½ | 7.6-9 | PSRAGNASH? |
| Narrow Regular | 00, 0, A0, A1 | 14-16 | 35-40 | 4-6 | 10-15 | 3-3½ | 7.6-9 | PSRAGNARG? |
| Narrow Long | 00, 0, A0, A1 | 18-20 | 45-50 | 4-6 | 10-15 | 3-3½ | 7.6-9 | PSRAGNALG? |
| Narrow X-Long | 00, 0, A0, A1 | 22-24 | 55-60 | 4-6 | 10-15 | 3-3½ | 7.6-9 | PSRAGNAXL? |
| Medium X-Short | 1, 2 | 8-10 | 20-25 | 7-8 | 17.5-20 | 4-5 | 10-13 | PSRAGMDXS? |
| Medium Short | 1, 2 | 10-12 | 25-30 | 7-8 | 17.5-20 | 4-5 | 10-13 | PSRAGMDSH? |
| Medium Regular | 1, 2 | 14-16 | 35-40 | 7-8 | 17.5-20 | 4-5 | 10-13 | PSRAGMDRG? |
| Medium Long | 1, 2 | 18-20 | 45-50 | 7-8 | 17.5-20 | 4-5 | 10-13 | PSRAGMDLG? |
| Medium X-Long | 1, 2 | 22-24 | 55-60 | 7-8 | 17.5-20 | 4-5 | 10-13 | PSRAGMDXL? |
| Medium XX-Long | 1, 2 | 26-28 | 65-70 | 7-8 | 17.5-20 | 4-5 | 10-13 | PSRAGMDXXL? |
| Wide X-Short | 3, 4 | 8-10 | 20-25 | 9-10 | 23-25.5 | 5½-6 | 13.5-15 | PSRAGWDXS? |
| Wide Short | 3, 4 | 10-12 | 25-30 | 9-10 | 23-25.5 | 5½-6 | 13.5-15 | PSRAGWDSH? |
| Wide Regular | 3, 4 | 14-16 | 35-40 | 9-10 | 23-25.5 | 5½-6 | 13.5-15 | PSRAGWDRG? |
| Wide Long | 3, 4 | 18-20 | 45-50 | 9-10 | 23-25.5 | 5½-6 | 13.5-15 | PSRAGWDLG? |
| Wide X-Long | 3, 4 | 22-24 | 55-60 | 9-10 | 23-25.5 | 5½-6 | 13.5-15 | PSRAGWDXL? |
| Wide XX-Long | 3, 4 | 26-28 | 65-70 | 9-10 | 23-25.5 | 5½-6 | 13.5-15 | PSRAGWDXXL? |
| Extra Wide Short | 5, 6 | 10-12 | 25-30 | 11-12 | 28-30.5 | 7-8 | 17.5-20 | PSRAGXWDSH? |
| Extra Wide Regular | 5, 6 | 14-16 | 35-40 | 11-12 | 28-30.5 | 7-8 | 17.5-20 | PSRAGXWDRG? |
| Extra Wide Long | 5, 6 | 18-20 | 45-50 | 11-12 | 28-30.5 | 7-8 | 17.5-20 | PSRAGXWDLG? |
| Extra Wide X-Long | 5, 6 | 22-24 | 55-60 | 11-12 | 28-30.5 | 7-8 | 17.5-20 | PSRAGXWDXL? |

Ply-Two™ Super Soft Sock

- Constructed of Acrylic and Lycra for strength and durability, while remaining soft and supple
- Nylon reinforced rounded toe eliminates the distal end pocket of conventional prosthetic socks
- Non-allergenic and double-fleeced
- Flat seamed top with a special non-ravel interlocked stitch
- Add 3 to the end of product number for a vulcanized 3S hole; add 3K to the end of product number for a knitted 3S hole (a vulcanized hole will be included unless a knitted hole is requested)



| Size | Fits Std. Sizes | Top | Toe | X-Short 8-10" | Short 10-14" | Regular 14-20" | Long 20-26" |
|--------|-----------------|---|------------------------------------|------------------|-----------------|-------------------|----------------|
| Narrow | Arm, 00, 0, 1 | 4 ¹ / ₂ "-7 ¹ / ₂ " | 2 ³ / ₄ "-5" | PTNXS | PTNS | PTNR | PTNL |
| Medium | 2, 3, 4, 5, 6 | 7 ¹ / ₄ "-12" | 4"-7 ¹ / ₂ " | PTMXS | PTMS | PTMR | PTML |

**Easy Care Prosthetic Sock**

- Made of 52% Wool, 38% Acrylic and 10% Nylon
- Machine washable and dryable on the delicate cycle
- Choose from either 3 ply, 5 ply or 6 ply
- Replace ? with desired length (10"-32", 2" increments)
- *Note: Sizes 2-6, where length is over 18", toe width changes
- Add 3 to end of product number for sock with distal hole

| Size | Top Width | Toe Width | Lengths | Product No. | | |
|------|-----------|---|---------|----------------|----------------|----------------|
| | | | | 3 Ply | 5 Ply | 6 Ply |
| 0 | 6" | 3 ¹ / ₂ " | 10"-32" | PSSN?03 | PSSN?05 | PSSN?06 |
| 1 | 7" | 4" | 10"-32" | PSSN?13 | PSSN?15 | PSSN?16 |
| 2 | 8" | 5" (10-18" length), 4" (20-32" length) | 10"-32" | PSSN?23 | PSSN?25 | PSSN?26 |
| 3 | 9" | 6" (10-18" length), 4 ¹ / ₂ " (20-32" length) | 10"-32" | PSSN?33 | PSSN?35 | PSSN?36 |
| 4 | 10" | 6 ¹ / ₂ " (10-18" length), 5" (20-32" length) | 10"-32" | PSSN?43 | PSSN?45 | PSSN?46 |
| 5 | 11" | 7" (10-18" length), 5 ¹ / ₂ " (20-32" length) | 10"-32" | PSSN?53 | PSSN?55 | PSSN?56 |
| 6 | 12" | 8" (10-18" length), 6" (20-32" length) | 10"-32" | PSSN?63 | PSSN?65 | PSSN?66 |

BK Prosthetic Tube Sock

- Made using a unique air float terry construction throughout entire length of the sock offering comfort and durability
- Low profile seam above the toe are for additional comfort
- Reversible, and can be worn with either the smooth side or the cushioned side against the skin
- Add 3 to end of product number for sock with distal hole
- 3 ply thickness

| Size | Max Size | Length | Product No. |
|-------|----------|--------|--------------|
| A1, 2 | 8"-10" | 12" | CPT12 |
| A1, 2 | 10"-14" | 16" | CPT16 |
| A1, 2 | 16"-20" | 22" | CPT22 |
| 2 | 20"-26" | 28" | CPT28 |



SoftWick Stretch Prosthetic Half Sock

- The ideal solution for managing distal volume changes with a liner
- Distal hole is knit into the sock to ensure a proper fit without unraveling
- CoolMax® and Lycra construction stretches to accommodate a multitude of shapes
- Four sizes ensures a proper fit for almost any application
- The color band knit into each sock allows easy size identification

| Size | Band Color | Length | Top Width | Toe Width | Product No. | |
|------|------------|--------|-----------|-----------|---------------|--------------|
| | | | | | with Hole | without Hole |
| 1 | Black | 3" | 3" | 2½" | SW213K | SW21 |
| 2 | Gray | 3" | 3½" | 3" | SW223K | SW22 |
| 3 | Royal Blue | 6" | 4" | 2¾" | SW233K | SW23 |
| 4 | Turquoise | 6" | 5½" | 3½" | SW243K | SW24 |



Prosthetic
Socks

Bio-Wick™ One Ply Sock

- Microsafe® Anti-Microbial provides continuous protection against odor-causing bacteria
- CoolMax transports moisture away from skin, keeping interior dryer and more comfortable
- Completely seamless; no stitching in toe area
- Exceptional stretch conforms to virtually any shape including Symes
- Made of 45% Cotton, 45% CoolMax®, 8% Acetate, and 2% Lycra
- Non-run construction allows for cutting to length if necessary
- Available with a vulcanized 3S hole

| Size | Fits Sock Sizes | | Top | Toe | Length | Product No. |
|--------------|-----------------|---------------|---|------------------------------------|--------|-------------|
| | Medium/Short | Medium/Reg. | | | | |
| Medium/Short | Arm, 00, 0, 1 | Arm, 00, 0, 1 | 3 ³ / ₄ "-7 ¹ / ₂ " | 3 ³ / ₄ "-6" | 14" | BWMS |
| Medium/Reg. | Arm, 00, 0, 1 | Arm, 00, 0, 1 | 3 ³ / ₄ "-7 ¹ / ₂ " | 3 ³ / ₄ "-6" | 20" | BWMR |
| Medium/Long | Arm, 00, 0, 1 | Arm, 00, 0, 1 | 3 ³ / ₄ "-7 ¹ / ₂ " | 3 ³ / ₄ "-6" | 28" | BWML |
| Wide/Short | 2, 3, 4, 5 | 2, 3, 4, 5 | 5"-11" | 4"-8" | 14" | BWWS |
| Wide/Reg. | 2, 3, 4, 5 | 2, 3, 4, 5 | 5"-11" | 4"-8" | 20" | BWWR |
| Wide/Long | 2, 3, 4, 5 | 2, 3, 4, 5 | 5"-11" | 4"-8" | 28" | BWWL |

**Knit-Rite - Dero® Hole-In-Toe Prosthetic Sock**

- Rotation control prosthetic sock, for use with roll-on suction suspension liners
- Minimizes or eliminates unwanted rotation and vertical movement between the outer socket and silicone suspension liner, gel interface or silicone insert
- Incorporates friction enhancing polymer strips into the prosthetic sock to restrict movement
- Provides silicone suction suspension in transfemoral cases where it may be contraindicated by allowing easier donning and volume management capability
- Polymer strips add virtually no thickness and are vertically oriented in standard A/K and B/K patterns to allow ply combinations, while maintaining the circumferential stretch inherent in the Soft-Sock



| Size | Sock | | | Product No. | | |
|------------------|------|---------|-----------------|-----------------|-----------------|--|
| | Size | Length | 1 Ply | 3 Ply | 5 Ply | |
| AK Regular-Short | 1-2 | 10"-14" | 1DA1RGSH | 1DA3RGSH | 1DA5RGSH | |
| BK Regular-Short | 1-2 | 10"-14" | 1DB1RGSH | 1DB3RGSH | 1DB5RGSH | |
| AK Wide-Medium | 3+ | 16"-20" | 1DA1WDMD | 1DA3WDMD | 1DA5WDMD | |

Proximal Soft-Sock®

- For managing proximal volume changes
- Tapered and open-ended soft-sock
- Knit with polyester/lycra
- Sold each

| Size | Sock Size | Length | Lightweight | Product No. | |
|----------------|-----------|--------|-----------------|-----------------|-----------------|
| | | | | 3 Ply | 5 Ply |
| Regular Short | 1, 2 | 10.5" | 1PS1RGSH | 1PS3RGSH | 1PS5RGSH |
| Regular Medium | 1, 2 | 14" | 1PS1RGMD | 1PS3RGMD | 1PS5RGMD |
| Wide Medium | 3+ | 14" | 1PS1WDMD | 1PS3WDMD | 1PS5WDMD |
| X-Wide Medium | 4+ | 14" | 1PS1XWMD | 1PS3XWMD | 1PS5XWMD |

**Soft-Sock with X-Static**

- Full fashion knit sock made of Polyester/X-Static/Lycra in easy-to-fit sizes
- The fiber wicks moisture away from skin
- The skin remains drier and more comfortable
- Soft-Sock feels soft and cuddly worn next to the skin and insulates against the cold

| Size | Sock Size | Length | Lightweight | Product No. | | |
|-----------------|-----------|--------|-----------------|-----------------|---------------------|--------------------|
| | | | | 1 Ply | 3 Ply Yellow Top | 5 Ply Green Top |
| Narrow Short | B, 0, 1 | 10-14 | 1SX1NASH | 1SX3NASH | 1SX5NASH | |
| Regular X-Short | 1, 2 | 8-12 | 1SX1RGXS | 1SX3RGXS | 1SX5RGXS | |
| Regular Short | 1, 2 | 10-14 | 1SX1RGSH | 1SX3RGSH | 1SX5RGSH | |
| Regular Medium | 1, 2 | 16-20 | 1SX1RGMD | 1SX3RGMD | 1SX5RGMD | |
| Wide Medium | 3+ | 16-20 | 1SX1WDMD | 1SX3WDMD | 1SX5WDMD | |

**Shorty Soft-Sock**

- Designed for use with elevated vacuum and total surface bearing sockets (TSB)
- The short length reduces volume mainly in the distal one third of the socket and remain below the suction sleeve
- Color keyed for ease of patient instruction
- Distal hole; knit with easy care Coolmax
- Sold each

| Toe Stripe | Top | Toe | Length | Product No. |
|------------|---------------------------------|---------------------------------|---------------------------------|-----------------|
| Red | 4 ³ / ₄ " | 2 ³ / ₄ " | 7" | 1SPSS181 |
| Red | 5 ³ / ₄ " | 3 ¹ / ₈ " | 5 ⁵ / ₈ " | 1SPSS191 |
| Black | 3 ¹ / ₈ " | 2 ⁵ / ₈ " | 3 ⁵ / ₈ " | 1SPSS182 |
| Black | 4 ¹ / ₂ " | 4" | 3 ¹ / ₄ " | 1SPSS192 |



SmartKnit® Partial Foot Sock

- For all levels of partial foot amputation
- Completely seamless, fits like a second skin
- 4-way stretch for smooth, comfortable fit
- Washable CoolMax/Lycra, wicks moisture from the skin
- Bacteria-fighting, anti-microbial fabric, non binding Halo top
- One size fits all, sold individually
- 3¹/₂" (9 cm) width x 14¹/₂" (37 cm) length

| Description | Product No. |
|-------------|--------------|
| White | 84107 |
| Black | 84707 |

**A-Plus™ (Acrylic/Lycra) Sock**

- The prosthetic sock made with Acrylic Plus Lycra Spandex for a garment featuring softness plus stretch, moisture wicking plus conforming and easy care plus long-wear
- Available in 2, 3, and 4 ply
- Because the Acrylic/Lycra Spandex yarn snugs the knit close together, a 2 ply (green stitching) approximates the thickness of a 3 ply wool sock; 3 ply (blue stitching) approximates 5 ply wool and the 4 ply (red stitching) approximates 6 ply wool
- All A-Plus socks will measure up to 10% less in width so that the size of the A-Plus sock is the same size that patient would wear in a wool sock*
- The length measurement is stretched in fleecing and will relax to stated measurement when washed
- Available in Lengths of 8"-24" in 2" increments
- Replace ? with Ply and Replace ?? with Length



| A-Plus Standard Socks | | | | |
|------------------------------|---------------------------------|---------------------------------|-----------|-----------------|
| Size | Length | Toe Width | Top Width | Product No. |
| 0 | 6" | 3 ¹ / ₂ " | 8"-24" | 1KL?00?? |
| 1 | 7" | 4" | 8"-24" | 1KL?01?? |
| 2 | 8" | 5" | 8"-24" | 1KL?02?? |
| 3 | 9" | 6" | 8"-20" | 1KL?03?? |
| 3 | 9" | 5 ¹ / ₂ " | 22"-24" | 1KL?03?? |
| 4 | 10" | 7" | 8"-20" | 1KL?04?? |
| 4 | 10" | 6" | 22"-24" | 1KL?04?? |
| Arm Socks | | | | |
| A | 4" | 3" | 6"-20" | 1KL?0A?? |
| B | 5 ¹ / ₂ " | 3 ¹ / ₂ " | 6"-20" | 1KL?0B?? |

*Note: Due to the stretch of A-Plus Socks, widths are up to 10% less.

A-Plus Hole-In-Toe Sock

- Made with a distal hole specifically for use with pin suspension systems
- The knit construction adds little or no additional distal thickness and eliminates the need for cutting the sock which helps prevent raveling
- The range of thicknesses allows precise fit adjustments with all pin systems
- Available in 2, 3, and 4 ply
- Available in Lengths of 8"-24" in 2" increments
- Replace ? with Ply and Replace ?? with Length



| A-Plus with Hole | | | | |
|-------------------------|------------------|------------------|---------------|--------------------|
| Size | Top Width | Toe Width | Length | Product No. |
| 0 | 6" | 3 1/2" | 8"-24" | 1KL?S0?? |
| 1 | 7" | 4" | 8"-24" | 1KL?S1?? |
| 2 | 8" | 5" | 8"-24" | 1KL?S2?? |
| 3 | 9" | 6" | 8"-20" | 1KL?S3?? |
| 3 | 9" | 6" | 22"-24" | 1KL?S3?? |
| 4 | 10" | 7" | 8"-20" | 1KL?S4?? |
| 4 | 10" | 6" | 22"-24" | 1KL?S4?? |

| A-Plus Hole-In-Toe Slim Taper PTB | | | | |
|--|------------------|------------------|---------------|--------------------|
| Size | Top Width | Toe Width | Length | Product No. |
| 0 | 6" | 4" | 10"-16" | 1TH?04?? |
| 0 | 6" | 5" | 10"-16" | 1TH?05?? |
| 1 | 7" | 5" | 10"-16" | 1TH?15?? |
| 2 | 8" | 6" | 10"-16" | 1TH?26?? |

| A-Plus Hole-In-Toe Wide Taper PTB | | | | |
|--|------------------|------------------|---------------|--------------------|
| Size | Top Width | Toe Width | Length | Product No. |
| 0 | 6" | 3" | 10"-16" | 1TH?03?? |
| 1 | 7" | 3 1/2" | 10"-16" | 1TH?13?? |
| 2 | 8" | 4" | 10"-16" | 1TH?24?? |

| A-Plus Slim Taper PTB | | | | |
|------------------------------|------------------|------------------|---------------|--------------------|
| Size | Top Width | Toe Width | Length | Product No. |
| 0 | 6" | 4" | 10"-16" | 1TL?04?? |
| 0 | 6" | 5" | 10"-16" | 1TL?05?? |
| 1 | 7" | 5" | 10"-16" | 1TL?15?? |
| 2 | 8" | 6" | 10"-16" | 1TL?26?? |

| A-Plus Wide Taper PTB | | | | |
|------------------------------|------------------|------------------|---------------|--------------------|
| Size | Top Width | Toe Width | Length | Product No. |
| 0 | 6" | 3" | 10"-16" | 1TL?03?? |
| 1 | 7" | 3 1/2" | 10"-16" | 1TL?13?? |
| 2 | 8" | 4" | 10"-16" | 1TL?24?? |

Soft Sock-Cuddly 'N Soft

- Full fashion knit sock made of CoolMax/Lycra spandex in easy-to-fit sizes
- The fiber wicks moisture away from skin
- The skin remains drier and more comfortable
- Soft-Sock feels soft and cuddly worn next to the skin and insulates against the cold



| Size | Fits Current Size | Length | Product No. | | | |
|-----------------|-------------------|---------|-----------------|-----------------|-----------------|-----------------|
| | | | 1 Ply | 3 Ply | 5 Ply | 6 Ply |
| Child X-Short | A | 6" | 1SP1CHXS | — | — | — |
| Child Short | A | 8" | 1SP1CHSH | 1SP3CHSH | 1SP5CHSH | 1SP6CHSH |
| Child Medium | A | 12" | 1SP1CH12 | 1SP3CH12 | 1SP5CH12 | — |
| Narrow Short | B, 0, 1 | 10"-14" | 1SP1NASH | 1SP3NASH | 1SP5NASH | 1SP6NASH |
| Narrow Medium | B, 0, 1 | 16"-20" | 1SP1NAMD | 1SP3NAMD | 1SP5NAMD | 1SP6NAMD |
| Narrow Long | B, 0, 1 | 22"-28" | 1SP1NALG | 1SP3NALG | 1SP5NALG | 1SP6NALG |
| Regular X-Short | 1, 2 | 8"-12" | 1SP1RGXS | 1SP3RGXS | 1SP5RGXS | 1SP6RGXS |
| Regular Short | 1, 2 | 10"-14" | 1SP1RGSH | 1SP3RGSH | 1SP5RGSH | 1SP6RGSH |
| Regular Medium | 1, 2 | 16"-20" | 1SP1RGMD | 1SP3RGMD | 1SP5RGMD | 1SP6RGMD |
| Regular Long | 1, 2 | 22"-28" | 1SP1RGLG | 1SP3RGLG | 1SP5RGLG | 1SP6RGLG |
| Wide Short | 3+ | 10"-14" | 1SP1WDSH | 1SP3WDSH | 1SP5WDSH | 1SP6WDSH |
| Wide Medium | 3+ | 16"-20" | 1SP1WDMD | 1SP3WDMD | 1SP5WDMD | 1SP6WDMD |
| X-Wide Short | 4+ | 10"-14" | 1SP1XWSH | 1SP3XWSH | 1SP5XWSH | — |
| X-Wide Medium | 4+ | 16"-20" | 1SP1XWMD | 1SP3XWMD | 1SP5XWMD | — |

| Soft Sock Hole-In-Toe | | | | | | |
|-----------------------|-------------------|---------|-----------------|-----------------|-----------------|-----------------|
| Size | Fits Current Size | Length | Product No. | | | |
| | | | 1 Ply | 3 Ply | 5 Ply | 6 Ply |
| Child Short | A | 8" | 1SPSCHSH | 1SP3SCSH | 1SP5SCSH | 1SP6SCSH |
| Narrow Short | B, 0, 1 | 10"-14" | 1SPSNASH | 1SP3SNSH | 1SP5SNSH | 1SP6SNSH |
| Narrow Medium | B, 0, 1 | 16"-20" | 1SPSNAMD | 1SP3SNMD | 1SP5SNMD | 1SP6SNMD |
| Regular X-Short | 1, 2 | 8"-12" | 1SPSRGXS | 1SP3SRXS | 1SP5SRXS | 1SP6SRXS |
| Regular Short | 1, 2 | 10"-14" | 1SPSRGSH | 1SP3SRSH | 1SP5SRSH | 1SP6SRSH |
| Regular Medium | 1, 2 | 16"-20" | 1SPSRGMD | 1SP3SRMD | 1SP5SRMD | 1SP6SRMD |
| Regular Long | 1, 2 | 22"-28" | 1SPSRGLG | 1SP3SRLG | 1SP5SRLG | 1SP6SRLG |
| Wide Short | 3+ | 10"-14" | 1SPSWDSH | 1SP3SWSH | 1SP5SWSH | 1SP6SWSH |
| Wide Medium | 3+ | 16"-20" | 1SPSWDMD | 1SP3SWMD | 1SP5SWMD | 1SP6SWMD |

Sterling Cotton Sock

- Available in 1, 2, 3, or 5 ply weight
- Machine or hand wash
- 1 and 2 ply not available in sizes 5 and 6
- Arm sock (AS & AL) available in 3 & 5 ply only
- Replace ? with Ply and Replace ?? with Length (2" increments)

| Size | Top Width | Toe Width | Length | Product No. |
|------|-----------|-----------|---------|------------------|
| 00 | 5" | 3" | 10"-26" | 1CS?D0??? |
| 0 | 6" | 3½" | 10"-26" | 1CS?10??? |
| 1 | 7" | 4" | 10"-26" | 1CS?11??? |
| 2 | 8" | 5" | 10"-26" | 1CS?22??? |
| 2 | 8" | 4" | 10"-26" | 1CS?02??? |
| 2 | 8" | 4½" | 10"-26" | 1CS?12??? |
| 3 | 9" | 6" | 10"-18" | 1CS?03??? |
| 3 | 9" | 4½" | 20"-26" | 1CS?03??? |
| 4 | 10" | 6½" | 10"-18" | 1CS?04??? |
| 4 | 10" | 5" | 20"-26" | 1CS?04??? |
| 5 | 11" | 7" | 10"-18" | 1CS?05??? |
| 5 | 11" | 5½" | 20"-26" | 1CS?05??? |
| 6 | 12" | 8" | 10"-18" | 1CS?06??? |
| 6 | 12" | 6" | 20"-26" | 1CS?06??? |
| AS | 4" | 3" | 6"-22" | 1CS?AS??? |
| AL | 5½" | 3½" | 6"-22" | 1CS?AL??? |



Sterling Stretch Hole-In-Toe

- Adds Lycra Spandex to the acrylic to allow the sock to conform to the shape of the limb better; Fleecing adds softness to the inside
- Wash and wear stretch; conforms to shape of limb
- Hemmed top and smooth widenings; reinforced toe with distal hole

| Size | Fits Current Size | Length (in) | Top (in) | Toe (in) | Product No. | | |
|-------------------|-------------------|-------------|----------|----------|-----------------|-----------------|-----------------|
| | | | | | 3 Ply | 5 Ply | 6 Ply |
| Narrow Short | 00, 0, AS, AL | 10-12 | 4-6 | 3-3½ | 1HS3NASH | 1HS5NASH | 1HS6NASH |
| Narrow Regular | 00, 0, AS, AL | 14-16 | 4-6 | 3-3½ | 1HS3NARG | 1HS5NARG | 1HS6NARG |
| Narrow Long | 00, 0, AS, AL | 18-20 | 4-6 | 3-3½ | 1HS3NALG | 1HS5NALG | 1HS6NALG |
| Medium Short | 1, 2 | 10-12 | 7-8 | 4-5 | 1HS3MDSH | 1HS5MDSH | 1HS6MDSH |
| Medium Regular | 1, 2 | 14-16 | 7-8 | 4-5 | 1HS3MDRG | 1HS5MDRG | 1HS6MDRG |
| Medium Long | 1, 2 | 18-20 | 7-8 | 4-5 | 1HS3MDLG | 1HS5MDLG | 1HS6MDLG |
| Medium X-Long | 1, 2 | 22-24 | 7-8 | 4-5 | 1HS3MDXL | 1HS5MDXL | 1HS6MDXL |
| Medium XX-Long | 1, 2 | 26-28 | 7-8 | 4-5 | 1HS3MDXX | 1HS5MDXX | 1HS6MDXX |
| Wide Short | 3, 4 | 10-12 | 9-10 | 6-7 | 1HS3WDSH | 1HS5WDSH | 1HS6WDSH |
| Wide Medium | 3, 4 | 14-16 | 9-10 | 6-7 | 1HS3WDMD | 1HS5WDMD | 1HS6WDMD |
| Wide Long | 3, 4 | 18-20 | 9-10 | 4½-5 | 1HS3WDLG | 1HS5WDLG | 1HS6WDLG |
| Wide X-Long | 3, 4 | 22-24 | 9-10 | 4½-5 | 1HS3WDXL | 1HS5WDXL | 1HS6WDXL |
| Wide XX-Long | 3, 4 | 26-28 | 9-10 | 4½-5 | 1HS3WDXX | 1HS5WDXX | 1HS6WDXX |
| Extra Wide Short | 5, 6 | 12-14 | 11-12 | 7-8 | 1HS3EWSH | 1HS5EWSH | 1HS6EWSH |
| Extra Wide Medium | 5, 6 | 16-18 | 11-12 | 7-8 | 1HS3EWMD | 1HS5EWMD | 1HS6EWMD |

Sterling Stretch Prosthetic Sock

- Adds Lycra spandex to the acrylic to allow the sock to conform to the shape of the limb better
- Fleecing adds softness to the inside next to the skin
- Available in 3, 5, or 6 ply weight
- Wash and wear stretch
- Conforms to shape of limb
- Transports moisture from skin
- Reinforced toe, hemmed top and smooth widenings



| Size | Fits Current Size | Length (in) | Top (in) | Toe (in) | Product No. | | |
|-------------------|-------------------|-------------|----------|----------|-----------------|-----------------|-----------------|
| | | | | | 3 Ply | 5 Ply | 6 Ply |
| Narrow X-Short | 00, 0, AS, AL | 8-10 | 4-6 | 3-3½ | 1SS3NAXS | 1SS5NAXS | 1SS6NAXS |
| Narrow Short | 00, 0, AS, AL | 10-12 | 4-6 | 3-3½ | 1SS3NASH | 1SS5NASH | 1SS6NASH |
| Narrow Regular | 00, 0, AS, AL | 14-16 | 4-6 | 3-3½ | 1SS3NARG | 1SS5NARG | 1SS6NARG |
| Narrow Long | 00, 0, AS, AL | 18-20 | 4-6 | 3-3½ | 1SS3NALG | 1SS5NALG | 1SS6NALG |
| Medium X-Short | 1, 2 | 8-10 | 7-8 | 4-5 | 1SS3MDXS | 1SS5MDXS | 1SS6MDXS |
| Medium Short | 1, 2 | 10-12 | 7-8 | 4-5 | 1SS3MDSH | 1SS5MDSH | 1SS6MDSH |
| Medium Regular | 1, 2 | 14-16 | 7-8 | 4-5 | 1SS3MDRG | 1SS5MDRG | 1SS6MDRG |
| Medium Long | 1, 2 | 18-20 | 7-8 | 4-5 | 1SS3MDLG | 1SS5MDLG | 1SS6MDLG |
| Medium X-Long | 1, 2 | 22-24 | 7-8 | 4-5 | 1SS3MDXL | 1SS5MDXL | 1SS6MDXL |
| Medium XX-Long | 1, 2 | 26-28 | 7-8 | 4-5 | 1SS3MDXX | 1SS5MDXX | 1SS6MDXX |
| Wide Short | 3, 4 | 10-12 | 9-10 | 6-7 | 1SS3WDSH | 1SS5WDSH | 1SS6WDSH |
| Wide Medium | 3, 4 | 14-16 | 9-10 | 6-7 | 1SS3WDMD | 1SS5WDMD | 1SS6WDMD |
| Wide Long | 3, 4 | 18-20 | 9-10 | 4½-5 | 1SS3WDLG | 1SS5WDLG | 1SS6WDLG |
| Wide X-Long | 3, 4 | 22-24 | 9-10 | 4½-5 | 1SS3WDXL | 1SS5WDXL | 1SS6WDXL |
| Wide XX-Long | 3, 4 | 26-28 | 9-10 | 4½-5 | 1SS3WDXX | 1SS5WDXX | 1SS6WDXX |
| Extra Wide Short | 5, 6 | 12-14 | 11-12 | 7-8 | 1SS3EWSH | 1SS5EWSH | 1SS6EWSH |
| Extra Wide Medium | 5, 6 | 16-18 | 11-12 | 7-8 | 1SS3EWMD | 1SS5EWMD | 1SS6EWMD |

Sterling Sock with Hole-In-Toe

- Available in 1, 3 or 5 ply
- 1 ply available only in size 1 and 2 (size 2 only in 8" top - 5" toe)
- Replace ? with Ply and Replace ?? with Length (2" increments)

| Sterling Socks with Hole | | | | | |
|--------------------------|-----------|-----------|---------|-----------------|-----------------|
| Size | Top Width | Toe Width | Length | Product No. | |
| | | | | 1 Ply | 3 or 5 ply |
| 0 | 6" | 3½" | 10"-12" | - | 1HS?10?? |
| 1 | 7" | 4" | 10"-14" | 1HS?11?? | 1HS?11?? |
| 2 | 8" | 4" | 10"-16" | - | 1HS?02?? |
| 2 | 8" | 4½" | 10"-16" | - | 1HS?12?? |
| 2 | 8" | 5" | 10"-16" | 1HS?22?? | 1HS?22?? |
| 3 | 9" | 6" | 12"-16" | - | 1HS?03?? |

Sterling Wool Sock

- 100% wool; 50% domestic virgin and 50% Australian virgin wool
- High degree of absorption
- Available 1, 2, 3, 5-8 ply (arm socks available in 3 or 5 ply only)
- 4 ply is not available in any size or style
- Machine wash 1-2 ply warm on gentle cycle
- Hand wash 3, 5-8 ply at warm temperature
- 1, 2, 7, and 8 ply in lengths 10"-26" only
- Replace ? with Ply and Replace ?? with Length (2" increments)



| Sterling Wool Socks | | | | | | |
|----------------------------------|-----------|-----------|---------|------------|---------------------------|--------------|
| Size | Top Width | Toe Width | Length | 1 or 2 ply | Product No. 3 or 5 ply | 6,7 or 8 ply |
| 00 | 5" | 3" | 10"-28" | 1MS?D0?? | 1WS?D0?? | 1WS?D0?? |
| 0 | 6" | 3 1/2" | 10"-28" | 1MS?10?? | 1WS?10?? | 1WS?10?? |
| 1 | 7" | 4" | 10"-28" | 1MS?11?? | 1WS?11?? | 1WS?11?? |
| 2 | 8" | 4" | 10"-28" | 1MS?02?? | 1WS?02?? | 1WS?02?? |
| 2 | 8" | 4 1/2" | 10"-28" | 1MS?12?? | 1WS?12?? | 1WS?12?? |
| 2 | 8" | 5" | 10"-28" | 1MS?22?? | 1WS?22?? | 1WS?22?? |
| 3 | 9" | 6" | 10"-18" | 1MS?03?? | 1WS?03?? | 1WS?03?? |
| 3 | 9" | 4 1/2" | 20"-28" | 1MS?03?? | 1WS?03?? | 1WS?03?? |
| 4 | 10" | 6 1/2" | 10"-18" | 1MS?04?? | 1WS?04?? | 1WS?04?? |
| 4 | 10" | 5" | 20"-28" | 1MS?04?? | 1WS?04?? | 1WS?04?? |
| 5 | 11" | 7" | 10"-18" | 1MS?05?? | 1WS?05?? | - |
| 5 | 11" | 5 1/2" | 20"-28" | 1MS?05?? | 1WS?05?? | - |
| 6 | 12" | 8" | 10"-18" | 1MS?06?? | 1WS?06?? | - |
| 6 | 12" | 6" | 20"-28" | 1MS?06?? | 1WS?06?? | - |
| Arm Socks 3 or 5 Ply Only | | | | | | |
| AS | 4" | 3" | 6"-22" | - | 1WS?AS?? | - |
| AL | 5 1/2" | 3 1/2" | 6"-22" | - | 1WS?AL?? | - |

Sterling Wool PTB Sock

- Replace ? with Ply and Replace ?? with Length (2" increments)
- Available in 3 ,5 or 6 ply

| PTB Size Sterling Wool Socks | | | |
|------------------------------|-----------|---------|-------------|
| Top Width | Toe Width | Length | Product No. |
| 6" | 3" | 10"-16" | 1WS?00?? |
| 6" | 3 1/2" | 10"-16" | 1WS?10?? |
| 6" | 4" | 10"-16" | 1WS?20?? |
| 6" | 5" | 10"-16" | 1WS?30?? |
| 7" | 3 1/2" | 10"-16" | 1WS?01?? |
| 7" | 4" | 10"-16" | 1WS?11?? |
| 7" | 4 1/2" | 10"-16" | 1WS?21?? |
| 7" | 5" | 10"-16" | 1WS?31?? |
| 8" | 4" | 10"-16" | 1WS?02?? |
| 8" | 4 1/2" | 10"-16" | 1WS?12?? |
| 8" | 5" | 10"-16" | 1WS?22?? |

Wool Super Sock

- Soft, resilient, absorbent, thick and protective
- 100% virgin wool - uniquely soft
- Machine wash and dry
- Available in 3, 5, or 6 ply (6 ply not available in size 6)
- Replace ? with Ply and Replace ?? with Length (2" increments)



| Size | Top Width | Toe Width | Length | Product No. |
|------|-----------|-----------|---------|-----------------|
| 0 | 6" | 3 1/2" | 8"-30" | 1K1?00?? |
| 1 | 7" | 4" | 8"-30" | 1K1?01?? |
| 2 | 8" | 5" | 8"-24" | 1K1?02?? |
| 2 | 8" | 4" | 26"-30" | 1K1?02?? |
| 3 | 9" | 6" | 8"-20" | 1K1?03?? |
| 3 | 9" | 5 1/2" | 22"-30" | 1K1?03?? |
| 4 | 10" | 7" | 8"-20" | 1K1?04?? |
| 4 | 10" | 6" | 22"-30" | 1K1?04?? |
| 5 | 11" | 7" | 8"-18" | 1K1?05?? |
| 5 | 11" | 5 1/2" | 20"-30" | 1K1?05?? |
| 6 | 12" | 8" | 8"-18" | 1K1?06?? |
| 6 | 12" | 6" | 20"-30" | 1K1?06?? |
| A | 4" | 3" | 5"-22" | 1K1?0A?? |
| B | 5 1/2" | 3 1/2" | 5"-22" | 1K1?0B?? |

Super Sock with Hole

- Made with a distal hole especially for use with suction suspension mechanisms
- Construction eliminates the need for cutting the sock
- Helps prevent raveling with little or no additional distal thickness
- Available in 3 or 5 ply
- Replace ? with Ply and Replace ?? with Length (2" increments)



| Super Socks with Hole | | | | |
|-----------------------|-----------|-----------|---------|-----------------|
| Size | Top Width | Toe Width | Length | Product No. |
| 0 | 6" | 3 1/2" | 10"-12" | 1K1?S0?? |
| 1 | 7" | 4" | 10"-14" | 1K1?S1?? |
| 2 | 8" | 5" | 10"-16" | 1K1?S2?? |
| 3 | 9" | 6" | 12"-16" | 1K1?S3?? |

X-WOOL™ Prosthetic Sock

- Virgin Wool - plaited inside with X-Static/Lycra corespun yarn, Lycra allows garment to contour without stretching out, accommodates flexion and allows fewer sizes
- Absorbent virgin wool is thick and cushioning
- Silver is naturally antimicrobial, enabling X-Static to safely inhibit the growth of odor causing bacteria and fungi
- Fleeced to add softness; machine washable - resists shrinkage
- Available in 3, 5 and 6 ply



| Size | Fits Sock Size | Length | 3 Ply | Product No. | 6 Ply Gray Top |
|--------------------|-------------------|---------|-----------------|--------------------|-------------------|
| | | | Yellow Top | 5 Ply Green Top | |
| Narrow X-Short | 00, 0, AS, AL | 8"-10" | 1WX3NAXS | 1WX5NAXS | — |
| Narrow Short | 00, 0, AS, AL | 10"-12" | 1WX3NASH | 1WX5NASH | 1WX6NASH |
| Narrow Regular | 00, 0, AS, AL | 14"-16" | 1WX3NARG | 1WX5NARG | 1WX6NARG |
| Narrow Long | 00, 0, AS, AL | 18"-20" | 1WX3NALG | 1WX5NALG | 1WX6NALG |
| Medium X-Short | 1, 2 | 8"-10" | 1WX3MDXS | 1WX5MDXS | — |
| Medium Short | 1, 2 | 10"-12" | 1WX3MDSH | 1WX5MDSH | 1WX6MDSH |
| Medium Regular | 1, 2 | 14"-16" | 1WX3MDRG | 1WX5MDRG | 1WX6MDRG |
| Medium Long | 1, 2 | 18"-20" | 1WX3MDLG | 1WX5MDLG | 1WX6MDLG |
| Medium X-Long | 1, 2 | 22"-24" | 1WX3MDXL | 1WX5MDXL | — |
| Medium XX-Long | 1, 2 | 26"-28" | 1WX3MDXX | 1WX5MDXX | — |
| Wide Short | 3, 4 | 10"-12" | 1WX3WDSH | 1WX5WDSH | 1WX6WDSH |
| Wide Medium | 3, 4 | 14"-16" | 1WX3WDMD | 1WX5WDMD | 1WX6WDMD |
| Wide Long | 3, 4 | 18"-20" | 1WX3WDLG | 1WX5WDLG | 1WX6WDLG |
| Wide X-Long | 3, 4 | 22"-24" | 1WX3WDXL | 1WX5WDXL | — |
| Wide XX-Long | 3, 4 | 26"-28" | 1WX3WDXX | 1WX5WDXX | — |
| X-Wide Short | 5, 6 | 12"-14" | 1WX3EWSH | 1WX5EWSH | — |
| X-Wide Medium | 5, 6 | 16"-18" | 1WX3EWMD | 1WX5EWMD | — |
| Hole-In-Toe | | | | | |
| Narrow X-Short | 00, 0, AS, AL | 8"-10" | 1WX3SNXS | 1WX5SNXS | — |
| Narrow Short | 00, 0, AS, AL | 10"-12" | 1WX3SNSH | 1WX5SNSH | — |
| Narrow Regular | 00, 0, AS, AL | 14"-16" | 1WX3SNRG | 1WX5SNRG | — |
| Narrow Long | 00, 0, AS, AL | 18"-20" | 1WX3SNLG | 1WX5SNLG | — |
| Medium Short | 1, 2 | 10"-12" | 1WX3SMSH | 1WX5SMSH | — |
| Medium Regular | 1, 2 | 14"-16" | 1WX3SMRG | 1WX5SMRG | — |
| Medium Long | 1, 2 | 18"-20" | 1WX3SMLG | 1WX5SMLG | — |
| Wide Short | 3, 4 | 10"-12" | 1WX3SWSH | 1WX5SWSH | — |
| Wide Medium | 3, 4 | 14"-16" | 1WX3SWMD | 1WX5SWMD | — |

Standard Symes Sock – Super Sock and Soft Sock

- Available in 3 or 5 ply
- Lengths available: 18", 20", 22", 24"
- Replace ?? with length

| Super Sock | | |
|---------------------------|--------------------|-----------------|
| <u>Top x Ankle x Ball</u> | <u>Product No.</u> | |
| | <u>3 Ply</u> | <u>5 Ply</u> |
| 6" x 3" x 3 1/2" | 1S13A0?? | 1S15A0?? |
| 6" x 3" x 4" | 1S13B0?? | 1S15B0?? |
| 6" x 3 1/2" x 4" | 1S13C0?? | 1S15C0?? |
| 6" x 3 1/2" x 4 1/2" | 1S13D0?? | 1S15D0?? |
| 7" x 3" x 3 1/2" | 1S13A1?? | 1S15A1?? |
| 7" x 3" x 4" | 1S13B1?? | 1S15B1?? |
| 7" x 3 1/2" x 4" | 1S13C1?? | 1S15C1?? |
| 7" x 3 1/2" x 4 1/2" | 1S13D1?? | 1S15D1?? |
| Soft Sock | | |
| 6" x 3" x 3 1/2" | 1SP3A0?? | 1SP5A0?? |
| 7" x 4" x 4 1/2" | 1SP3E1?? | 1SP5E1?? |
| 7" x 5" x 5 1/2" | 1SP3H1?? | 1SP5H1?? |



SPS - Cool Blue Sock

- Made from CoolMax Lycra® and Spandex fabric the Cool Blue Sock provides excellent protection
- Wicks moisture keeping skin dry and comfortable
- Stretches and recovers to conform to the limb and it's fine knit offers comfort and longevity
- The CoolMax Lycra and Spandex fabric retards odor in the sock
- Comes with or without a hole and is available in several easy-to-fit sizes in your choice of 1, 3, 5 or 6 ply
- Replace ? with ply



| Size | Length | Toe Width | Top Width | <u>Product No.</u> | |
|------|---------|-----------|-----------|--------------------|----------------|
| | | | | without Hole | with Hole |
| B, 0 | 8"-12" | 2 1/2" | 4" | CB?NAXS | CB?SNXS |
| B, 0 | 10"-14" | 2 1/2" | 4" | CB?NASH | CB?SNSH |
| B, 0 | 16"-20" | 2 1/2" | 4" | CB?NARG | CB?SNRG |
| B, 0 | 22"-26" | 2 1/2" | 4" | CB?NALG | CB?SNLG |
| 1, 2 | 8"-12" | 3 1/2" | 6"-6 1/2" | CB?MDXS | CB?SMXS |
| 1, 2 | 10"-14" | 3 1/2" | 6"-6 1/2" | CB?MDSH | CB?SMSH |
| 1, 2 | 16"-20" | 3 1/2" | 6"-6 1/2" | CB?MDRG | CB?SMRG |
| 1, 2 | 22"-26" | 3 1/2" | 6"-6 1/2" | CB?MDLG | CB?SMLG |
| 3, 4 | 10"-14" | 4 1/2"-5" | 8"-8 1/2" | CB?WDSH | CB?SWSH |
| 3, 4 | 16"-20" | 4 1/2"-5" | 8"-8 1/2" | CB?WDRG | CB?SWRG |
| 3, 4 | 22"-26" | 4 1/2"-5" | 8"-8 1/2" | CB?WDLG | CB?SWLG |
| 5, 6 | 10"-14" | 6 1/2" | 10" | CB?XWSH | CB?SXSH |
| 5, 6 | 16"-20" | 6 1/2" | 10" | CB?XWRG | CB?SXRG |

Knit-Rite - Liner-Liner Prosthetic Sock

- Suspension liner skin interface relieves skin, shear irritations and improves comfort with liners
- Washable, thin material facilitates grip; high stretch for improved fit
- X-Static® Silver fibers fight bacteria and odor; transport heat and moisture
- 5.8% Lycra Spandex, 94.2% Polyester/X-Static (5%)

| Size | Sock Size | Length | Product No. |
|-----------------|-----------|---------|-----------------|
| Narrow X-Short | B, 0 | 8"-12" | 1SXLNAXS |
| Narrow Short | B, 0 | 10"-14" | 1SXLNASH |
| Regular X-Short | 1, 2 | 8"-12" | 1SXLRGXS |
| Regular Short | 1, 2 | 10"-14" | 1SXLRGSH |
| Regular Medium | 1, 2 | 16"-20" | 1SXLRGMD |
| Wide Short | 3 | 10"-14" | 1SXLWDSH |
| Wide Medium | 3 | 16"-20" | 1SXLWDMD |



Hip Disarticulation Prosthetic Sock

- 2-ply and 3-ply acrylic/Lycra
- 97% Orlon Acrylic; 3% Lycra spandex; Latex-free
- Modified box toe; fleeced inside
- Sold individually (4 socks to a box)
- On soundside, has 3" leg extension

*Hip Disarticulation
Prosthetic Socks*



| Flat Width | Length | Product No. | | | | |
|-------------|-------------|-----------------|-----------------|---------------------|---------------------|-----------------|
| | | 3-Ply Soft Sock | 5-Ply Soft Sock | 2-Ply Acrylic/Lycra | 3-Ply Acrylic/Lycra | 3-Ply X-Static |
| 8" (20 cm) | 12" (30 cm) | — | — | — | — | 1HX30812 |
| 10" (25 cm) | 12" (30 cm) | — | — | — | — | 1HX31012 |
| 12" (30 cm) | 14" (36 cm) | — | — | — | — | 1HX31214 |
| 12" (30 cm) | 16" (40 cm) | 1HP31216 | 1HP51216 | 1HL21216 | 1HL31216 | — |
| 13" (33 cm) | 14" (36 cm) | — | — | — | — | 1HX31314 |
| 13" (33 cm) | 16" (40 cm) | 1HP31316 | 1HP51316 | 1HL21316 | 1HL31316 | — |
| 14" (36 cm) | 14" (36 cm) | — | — | — | — | 1HX31414 |
| 14" (36 cm) | 16" (40 cm) | 1HP31416 | 1HP51416 | 1HL21416 | 1HL31416 | — |
| 15" (38 cm) | 16" (40 cm) | 1HP31516 | 1HP51516 | 1HL21516 | 1HL31516 | 1HX31516 |
| 16" (40 cm) | 16" (40 cm) | 1HP31616 | 1HP51616 | 1HL21616 | 1HL31616 | 1HX31616 |
| 17" (44 cm) | 16" (40 cm) | 1HP31716 | 1HP51716 | 1HL21716 | 1HL31716 | 1HX31716 |
| 18" (46 cm) | 16" (40 cm) | 1HP31816 | 1HP51816 | 1HL21816 | 1HL31816 | — |
| 18" (46 cm) | 18" (46 cm) | — | — | — | — | 1HX31818 |
| 19" (48 cm) | 16" (40 cm) | 1HP31916 | 1HP51916 | 1HL21916 | 1HL31916 | — |
| 19" (48 cm) | 18" (46 cm) | — | — | — | — | 1HX31918 |

A/K and B/K Post-Op Socks

- 3-ply weight
- 97% Orlon Acrylic; 3% Lycra spandex; Latex-free
- Modified box toe; fleeced inside
- Sold individually (4 socks to a box)

| Length | Top Width | Toe Width | <u>Product No.</u> | |
|--------|-----------|-----------|--------------------|-----------------|
| | | | Pre-Sterilized | Unsterilized |
| 14" | 5" | 3 1/2" | — | 1PL3BS14 |
| 18" | 5" | 4" | 1PL3S418 | 1PL3BS18 |
| 18" | 6" | 5" | 1PL3S518 | 1PL30518 |
| 18" | 8" | 5" | 1PL3S618 | 1PL32618 |
| 18" | 9" | 7" | 1PL3S718 | 1PL33718 |

**Comfort - Cast Sock**

- Comfort Cast socks are constructed from a combed natural cotton with a stretch stitched toe seam
- Available in BK or AK styles and either lightweight (1 ply) or heavyweight (2 ply)
- Add 3 to end of product number for sock with distal hole



| AK Heavyweight | | | |
|-----------------------|----------|----------------|----------------|
| Description | Pack Qty | Round Toe | Pointed Toe |
| 20" | 12 | CSHRA20 | CSHPA20 |
| 24" | 12 | CSHRA24 | CSHPA24 |
| BK Lightweight | | | |
| 20" | 24 | CSLR20 | CSLP20 |
| 24" | 24 | CSLR24 | CSLP24 |
| 30" | 24 | CSLR30 | CSLP30 |
| 36" | 24 | CSLR36 | CSLP36 |
| BK Heavyweight | | | |
| 20" | 12 | CSHR20 | CSHP20 |
| 24" | 12 | CSHR24 | CSHP24 |
| 30" | 12 | CSHR30 | CSHP30 |
| 36" | 12 | CSHR36 | CSHP36 |

Ultra Fit Filler Sock

- The stretchable casting sock that conforms to irregular shapes including Symes and Chopart
- Nylon reinforced toe increases durability
- Full one ply thickness; 18" length
- Can be cut without running or unraveling
- Cotton and Lycra construction
- One size fits 90% of requirements for length and width
- 6 per package



| Description | Max Circ. | Length | Product No. | |
|-------------|-----------|--------|---------------|----------------|
| | | | Without hole | With hole |
| 1/2 Ply BK | 22" | 18" | CUFL | CUFL3 |
| 1 Ply BK | 22" | 18" | CUFH | CUFH3 |
| 1 Ply AK | 28" | 18" | CUFHAK | CUFHAK3 |
| Arm | 12" | 18" | CUFA | CUFA3 |
| Pediatric | 12" | 18" | CUFP | CUFP3 |

Knit-Rite - Accucast™ Cast Socks

- A smooth and 4-way stretchy cotton/Lycra provides a light compression and high conformity for more accurate impressions
- Smooth surface provides a smoother, more accurate model and reduced modification time
- Seamed toe
- Sold per dozen

| Width x Length | Product No. |
|--|-----------------|
| BK-2 ³ / ₄ " (7 cm) x 24" (60 cm) | 6SCST224 |
| AK-6 ¹ / ₄ " (16 cm) x 24" (60 cm) | 6CSR624 |



Cast Socks

- Flat, tubular knit, soft, carded cotton yarn, seamed toe
- Sold per dozen

| Description | Width x Length | Stretch Range | Product No. |
|----------------|---------------------------------------|---------------|-----------------|
| BK Lightweight | 3 ¹ / ₂ " x 24" | 7"-24" | 6S7924LT |
| BK Heavyweight | 3 ¹ / ₂ " x 24" | 7"-20" | 6S7924HV |
| AK Lightweight | 5 ¹ / ₂ " x 22" | 11"-28" | 6S7918LT |
| AK Heavyweight | 6" x 22" | 12"-30" | 6S7918HV |

RX Textiles - Tufftoe

- A Cotton/Lycra (95% Cotton/5% Lycra) elastic fitting sock specifically designed for the locking device pin on all silicone sleeves
- The innovative super tough 100% nylon reinforced toe and the addition of a tightly sewn circular hole at the distal end make this sock virtually indestructible
- A special thread and gel further strengthens the opening providing exceptional strength and durability
- Sold in packs of 6



| Description | Length | Width | Product No. |
|-------------|--------|--------|----------------|
| 1 Ply BK | 11" | 3 1/2" | 3000160 |
| 2 Ply BK | 11" | 4" | 3000260 |
| 4 Ply BK | 11" | 4 1/2" | 3000360 |
| 1 Ply AK | 14" | 3 1/2" | 3000460 |
| 2 Ply AK | 14" | 4" | 3000560 |
| 4 Ply AK | 14" | 4 1/2" | 3000660 |

Elastic Fitting Sock (EFS)

- Made of cotton/Lycra spandex/nylon has the added benefit of exceptional elasticity for improved fitting qualities, range of motion and reduced sock migration
- Machine wash; tumble dry
- 70% cotton, 25% polypropylene, 5% Lycra




| Description | Length | Width | Product No. |
|---------------|------------------------------|---|----------------|
| Small/Medium | 12" relaxed 29" stretched | 4" relaxed 12" stretch at top 10" at the bottom | 3001500 |
| Large/X-Large | 15" relaxed 36" stretched | 4" relaxed 15" stretch at top 10" at the bottom | 3001600 |

Fitting Socks

- 70% wool, 30% polyester
- Great for initial fitting
- Works well as a supplement to stump socks
- Sold per dozen

| Description | Product No. |
|------------------|----------------|
| A. One Ply | 2011201 |
| B. Heavy Fitting | 2011001 |



SPS - Cast Socks


- 100% Cotton
- 12 per package

| Description | Product No. |
|--------------------|----------------|
| Thin, Regular 27" | 2010501 |
| Heavy, Regular 27" | 2010101 |

**CATCAM Casting Pants/Sock**

- 50% cotton/50% tubigrip
- Sold individually

| Description | Product No. |
|-------------|----------------|
| Regular | 2011801 |
| Large | 2011901 |

**CATCAM Casting Pants/Sock**

- 100% cotton
- Sold individually

Product No. **2011401**

Nyserts

| Description | Pack Qty | Product No. |
|-------------|----------|---------------|
| AK | 25 | 247253 |
| AK | 100 | 247251 |
| BK | 25 | 247252 |
| BK | 100 | 247250 |



Casting/Fitting
Socks

Knit-Rite - Spacer Socks

- Quality 1 ply, light or heavy weight fitting socks feature smooth distal end and fabric
- Toe is rounded and reinforced for good fit and long wear
- The blended cotton polyester fabric is soft, comfortable and absorbent
- 6 per package
- "Light" is a looser weave, "Heavy" is a tighter weave

| Weight | Fits Sizes | Length | Product No. |
|--------------------|------------|-------------|-------------------|
| Narrow Light | #B, #0 | 15" (38 cm) | 6S-8035-LT |
| Narrow Heavy | #B, #0 | 15" (38 cm) | 6S-8035-HV |
| Medium Light | #1, #2 | 18" (45 cm) | 6S-8040-LT |
| Medium Heavy | #1, #2 | 18" (45 cm) | 6S-8040-HV |
| Wide Light | #3, #4 | 18" (45 cm) | 6S-8045-LT |
| Wide Heavy | #3, #4 | 18" (45 cm) | 6S-8045-HV |
| X-Wide Heavy | #4, #5 | 24" (60 cm) | 6S-8050-HV |
| XX-Wide Heavy | #5, #6 | 18" (45 cm) | 6S-8055-HV |
| XX-Wide Heavy Long | #5, #6 | 24" (60 cm) | 6S-8060-HV |

***Hole-In-Toe Spacer Socks***

- Same as Spacer Socks but with reinforced hole
- 6 per package


| Weight | Fits Sizes | Length | Product No. |
|--------------|------------|-------------|-------------------|
| Narrow Heavy | #B, #0 | 15" (38 cm) | 6S-S035-HV |
| Medium Heavy | #1, #2 | 18" (45 cm) | 6S-S040-HV |
| Wide Heavy | #3, #4 | 18" (45 cm) | 6S-S045-HV |

Stretch Spacer Sock

- Same features as Spacer Socks with Lycra added for improved fitting qualities, ROM and reduced migration
- Medium single ply thickness
- Elasticity allows comfortable fit for more than 50% greater range than the non-elastic version
- 6 per package



| Cotton/Lycra with Turn Down Cuff, Cotton/Lycra Toe | | | | | |
|---|------------|--------------------|-------------------|-------------------|-------------|
| Description | Fits Size | Unstretched Length | Regular | Hole-In-Toe | Product No. |
| Narrow Short | #A, #B | 12" (30 cm) | 6S-8S35-SH | 6S-SS35-SH | |
| Medium Short | #0, #1 | 12" (30 cm) | 6S-8S40-SH | 6S-SS40-SH | |
| Wide Short | #2, #3, #4 | 12" (30 cm) | 6S-8S45-WS | 6S-SS45-WS | |
| Cotton/Lycra with Nylon Toe, Straight Top | | | | | |
| Narrow | #A, #B | 12" (30 cm) | 6S-8S35-MD | 6S-SS35-MD | |
| Medium | #0, #1 | 15" (38 cm) | 6S-8S40-MD | 6S-SS40-MD | |
| Wide Short | #2, #3, #4 | 12" (30 cm) | 6S-8S45-SH | 6S-SS45-SH | |
| Wide | #2, #3, #4 | 18" (38 cm) | 6S-8S45-MD | 6S-SS45-MD | |



Graduated Length Volume Filler Socks


- A ready-made solution for common distal volume change
- Made by attaching toe section of a standard A-Plus prosthetic sock to a stretchy sheath-type proximal section
- Overall length is 18"; can be cut shorter
- Sold each

| Size | Width | Toe Length | Product No. | | | |
|------|--------------|--------------|-----------------|-----------------|-------------------|-------------------|
| | | | 2-Ply Standard | 3-Ply Standard | 2-Ply Hole-In-Toe | 3-Ply Hole-In-Toe |
| 0 | 3.5" (9 cm) | 3" (7.5 cm) | 1CL20003 | 1CL30003 | 1CL2S003 | 1CL3S003 |
| 0 | 3.5" (9 cm) | 4" (10 cm) | 1CL20004 | 1CL30004 | 1CL2S004 | 1CL3S004 |
| 0 | 3.5" (9 cm) | 5" (12.5 cm) | 1CL20005 | 1CL30005 | 1CL2S005 | 1CL3S005 |
| 0 | 3.5" (9 cm) | 6" (15 cm) | 1CL20006 | 1CL30006 | 1CL2S006 | 1CL3S006 |
| 1 | 4" (10 cm) | 3" (7.5 cm) | 1CL20103 | 1CL30103 | 1CL2S103 | 1CL3S103 |
| 1 | 4" (10 cm) | 4" (10 cm) | 1CL20104 | 1CL30104 | 1CL2S104 | 1CL3S104 |
| 1 | 4" (10 cm) | 5" (12.5 cm) | 1CL20105 | 1CL30105 | 1CL2S105 | 1CL3S105 |
| 1 | 4" (10 cm) | 6" (15 cm) | 1CL20106 | 1CL30106 | 1CL2S106 | 1CL3S106 |
| 2 | 5" (12.5 cm) | 3" (7.5 cm) | 1CL20203 | 1CL30203 | 1CL2S203 | 1CL3S203 |
| 2 | 5" (12.5 cm) | 4" (10 cm) | 1CL20204 | 1CL30204 | 1CL2S204 | 1CL3S204 |
| 2 | 5" (12.5 cm) | 5" (12.5 cm) | 1CL20205 | 1CL30205 | 1CL2S205 | 1CL3S205 |
| 2 | 5" (12.5 cm) | 6" (15 cm) | 1CL20206 | 1CL30206 | 1CL2S206 | 1CL3S206 |

SPS - Filler Socks

- 1 ply; rounded toe
- 50/50 Cotton-Polyester
- Nylon reinforced toe
- 12 per package

| Description | Product No. | |
|-------------|----------------|-------------------|
| | Without hole | With distal hole |
| 3.5" x 19" | 1513-HV | WH-1513-HV |
| 4" x 12" | 1519-HV | O1519-HV |
| 5" x 12" | 1521-HV | O1521-HV |
| 5" x 15" | 1515-HV | O1515-HV |
| 5" x 19" | 1512-HV | O1512-HV |
| 6.5" x 24" | 1511-HV | O1511-HV |
| 9" x 15" | 1525-HV | WH-1525-HV |
| 9" x 23" | 1516-HV | O1516-HV |



Filler/Spacer
Socks

Comfort Products - Angel Prosthetic Sock

- Hybrid textile and gel sock; combines X-Static Silver and CoolMax for natural anti-microbial protection while wicking moisture away from the skin
- Allows the wearer to stay dry while eliminating odor causing bacteria
- Fleeced interior; Exceptional stretch in both length and width-easily conforms to irregular shapes

| Sock Length | Gel Length | Product No. | |
|-------------|------------|-----------------|------------------|
| | | Regular | with 3S Hole |
| 12 | 9 | ANGELY | ANGELY3 |
| 12 | 12 | ANGELY12 | ANGELY123 |
| 14 | 9 | ANGELY14 | ANGELY143 |

**ComforGel™ Liner**

- Is a unique design that incorporates a special gel formulation that is enclosed between an inner and outer covering of cotton stretch fabric
- The cotton stretch fabric allows the liner to accommodate a 20 inch circumference from its relaxed size of 4 inches
- Toe is nylon reinforced for added strength
- Conforms to irregular shapes and contours
- 3 ply lightweight yellow stitching, 5 ply heavyweight green stitching



| Ply | Liner Type | Sock Length | Product No. | |
|-----|------------|-------------|-----------------|------------------|
| | | | Regular | with 3S Hole |
| 3 | BK | 10 | CGELY10 | CGELY103 |
| 3 | BK | 12 | CGELY12 | CGELY3 |
| 3 | AK | 12 | CGELAY | CGELAY3 |
| 3 | BK | 14 | CGELY14 | CGELY143 |
| 5 | BK | 10 | CGELG10 | CGELG103 |
| 5 | BK | 12 | CGELG12 | CGELG123 |
| 5 | AK | 12 | CGELAG | CGELY3 |
| 5 | BK | 14 | CGELG14 | CGELG143 |
| 5 | AK | 14 | CGELAG14 | CGELAG143 |

Silipos - Partial Foot Sock

- Designed to provide comfort and protection to a partial foot amputee
- Specially designed two-ply sock (CoolMax® material) is laminated around the sock's "toe" with the Silipos patented 3 mm vitamin enriched gel



| Description | Product No. |
|-------------|--------------|
| Regular | 10805 |
| X-Long | 10815 |

Symes Gel Liner™

- The Symes Gel Liner is a CoolMax® sock that is lined with the Silipos medical grade mineral oil gel to minimize skin breakdown
- Hypoallergenic
- One size fits most
- Length: 16½" (42 cm), Gel Length: 15" (38 cm)
- Distal Circumference: 10"-14½" (25-37 cm)
- Proximal Circumference: 15"-20" (38-51 cm)
- 5 ply, 3 mm gel - vitamin enriched

Product No.

1302

**SPS - Socket Gel Liner™**

- BK gel liner laminated with patented tri-block polymer gel
- Ensures comfortable fit

Single Socket Gel Liner:

- 3 ply, 3 mm gel
- CoolMax exterior material
- Tapered shape
- Gel against the skin

Double Socket Gel Liner:

- 5 ply, 3 mm gel
- Two layers of CoolMax (gel in between layers)
- Tapered shape



| Description | Length | Toe Circ. | Top Circ. | Product No. | |
|----------------|--------|-----------|-----------|--------------|--------------|
| | | | | Single | Double |
| Short/Narrow | 12" | 9"-12" | 15"-18" | 12725 | 92735 |
| Short/Regular | 12" | 12"-14" | 15"-18" | 12755 | 92745 |
| Medium/Regular | 16" | 12"-14" | 15"-18" | 12765 | 92755 |

Ultra Cushion Gel Liner™

- 3 ply or 5 ply, 3 mm gel
- Two layers of cotton, nylon, and Lycra, coated with Silipos gel
- Best combination of protection, quality and value
- Nylon reinforced toe for added durability



| Size | Description | Length of Liner | Distal Circ. | Proximal Circ. | Product No. | |
|------|----------------|-----------------|--------------|----------------|--------------|--------------|
| | | | | | 3 Ply | 5 Ply |
| I | Short Regular | 12½" | 7"-11" | 8"-13" | 92705 | 92706 |
| II | Medium/Regular | 15"-17" | 8"-12½" | 8"-14" | 92715 | — |
| II | Medium/Regular | 17"-19" | 8"-12½" | 8"-14" | — | 92716 |
| III | Medium/Wide | 17"-19" | 9"-14" | 9"-16" | 92725 | — |
| III | Medium/Wide | 18"-21" | 9"-14" | 9"-16" | — | 92726 |

ALPS - Skin Reliever

- Skin reliever eliminates shear forces, friction, and prevents abrasions on the skin of the stump
- Replaces pelite and socks
- Reduces shearing and abrasions; softens and soothes skin
- Tapered from 3 mm - 1.5 mm
- Warranty: 1 month

| Size | Product No. |
|---|------------------|
| 10-19 cm | ENCP-16-3 |
| 20-27 cm | ENCP-20-3 |
| 28-31 cm | ENCP-28-3 |
| Measurements taken 6 cm from the distal end | |

**Comfort Products - Comfort Sheath**

- Toe: reinforced for greater wear resistance
- Body: unique open mesh design which stretches to conform to shape while allowing for needed ventilation
- Top: turned welt construction providing added strength where needed
- Add to 3 to end of product number for sock with distal hole

| Description | Length | Product No. |
|----------------|---------|-------------|
| BK Extra Short | 6"-8" | SHXS |
| BK Short | 10"-14" | SHS |
| BK Medium | 16"-20" | SHM |
| BK Long | 22"-26" | SHL |
| AK Medium | 16"-20" | SHA |

**Silver Liner Sheath**

- Designed to be worn underneath any liner, without compromising suspension
- Antimicrobial, Anti-odor, All Natural
- Thermodynamic: Cooler in summer or warmer in winter

| Description | Length | Product No. | |
|-------------|---------|---------------|----------------|
| | | Regular | with 3S Hole |
| BK XX-Short | 6"-8" | XSHXXS | XSHXXS3 |
| BK X-Short | 8"-10" | XSHXS | XSHXS3 |
| BK Short | 10"-12" | XSHS | XSHS3 |
| BK Medium | 14"-18" | XSHM | XSHM3 |
| BK Long | 20"-24" | XSHL | XSHL3 |
| AK Medium | 14"-18" | XSHA | XSHA3 |



Knit-Rite - Silver Sheath

- Sheaths control friction, adjust thickness, provide strength and aid skin care
- Knit-Rite Sheaths have transfers on both sheaths of the pair to help patients position the sheath to maximize comfort
- The Silver Sheath with X-Static® is anti-odor, it also wicks moisture away from the skin and aids in thermal consistency
- PROPER CARE: Sheaths may be hand or machine washed, use soap or mild detergent in warm water. (NO BLEACH.) Rinse thoroughly, air dry overnight or tumble dry low
- Sold by the Box (3 pair/box)



| Description | Fits Sizes | Length | Product No. |
|-----------------|--------------|---------|------------------|
| Narrow X-Short | #B, #0, #1 | 8"-10" | 2HPXNAXS |
| Narrow Short | #B, #0, #1 | 10"-14" | 2HPXNASH |
| Regular X-Short | #1, #2 | 8"-10" | 2HPXRGXS |
| Regular Short | #1, #2 | 10"-14" | 2HPXRGSH |
| Regular Medium | #1, #2 | 16"-20" | 2HPXRGMD |
| Wide Short | #3 and above | 10"-14" | 2HPXWD SH |

Hole in Toe Prosthetic Sheath

- Sheaths control friction, adjust thickness and aid skin care
- Can be worn with prosthetic socks or alone
- 100% stretch nylon
- The qualities of stretch nylon enable fitting of somewhat greater widths or lengths
- Have transfers on both sheaths of the pair to help patients position sheath to give greatest comfort
- Sold by the Box (3 pair/box)
- Available colors: White and Dusk



| Sock Size | Fits Sizes | Length | Product No. | |
|---------------|------------|--------|-------------------|-------------------|
| | | | White | Dusk |
| Short/Narrow | #B, #0 | 12" | 2H-PWNH-SH | 2H-PBNH-SH |
| Short/Regular | #1, #2 | 12" | 2H-PWRH-SH | 2H-PBRH-SH |
| Short/Wide | #3 & above | 12" | 2H-PWWH-SH | 2H-PBWH-SH |

Super Stretch Sheath

- Super Stretch Sheaths have Lycra added to the nylon to accommodate an even greater range and more conforming fit, enhance flexion, movement and improve suspension
- The surface texture provides improved suspension
- Sold by the Box (3 pair/box)

Sheaths



| Sock Size | Fits Sizes | Length | Product No. | |
|----------------|------------|----------|-----------------|-----------------|
| | | | White | Dusk |
| Short/Regular | #1, #2 | 10"-14"+ | 2HPSWRSH | 2HPSBRSH |
| Medium/Regular | #1, #2 | 16"-20"+ | 2HPSWRMD | 2HPSBRMD |

Prosthetic Sheath

- Sheaths control friction and aid skin care
- Can be worn with prosthetic socks or alone
- 4 widths: child, narrow, regular, and wide-sheath corresponds to prosthetic sock sizes
- The qualities of stretch nylon enable fitting of somewhat greater widths or lengths
- Have transfers on both sheaths of the pair to help patients position sheath to give greatest comfort
- 100% stretch nylon; Machine wash/dry
- Sold by the Box (3 pair/box)
- Available colors: White and Dusk




| Sock Size | Fits Sizes | Length | Product No. | |
|-----------------|-----------------------|---------|-----------------|-----------------|
| | | | White | Dusk |
| Child Short | #A, Lycra spandex top | 6"-8" | 2HPWCHSH | 2HPBCHSH |
| Child 12" | #A, non-elastic top | 12" | 2HPWCH12 | 2HPBCH12 |
| Narrow Short | #B, #0, #1 | 10"-14" | 2HPWNASH | 2HPBNASH |
| Narrow Medium | #B, #0, #1 | 16"-20" | 2HPWNAMD | 2HPBNAMD |
| Narrow Long | #B, #0, #1 | 22"-28" | 2HPWNALG | 2HPBNALG |
| Regular X-Short | #1, #2 | 8"-10" | 2HPWRGX5 | 2HPBRGX5 |
| Regular Short | #1, #2 | 10"-14" | 2HPWRGSH | 2HPBRGSH |
| Regular Medium | #1, #2 | 16"-20" | 2HPWRGMD | 2HPBRGMD |
| Regular Long | #1, #2 | 22"-28" | 2HPWRGLG | 2HPBRGLG |
| Wide Short | #3 & above | 10"-14" | 2HPWWDSH | 2HPBWDSH |
| Wide Medium | #3 & above | 16"-20" | 2HPWWDMD | 2HPBWDMD |
| Wide Long | #3 & above | 22"-28" | 2HPWWDLG | 2HPBWDLG |

SPS Hole in Toe Sheath

| Description | Product No. |
|-------------|-----------------|
| A/K | SPS-AK-2 |
| B/K | SPS-BK-2 |

A/K Brim Sheath

- Designed for above knee patients, for all socket styles
- Separates and protects the soft tissue in the top brim of the socket; used to bridge the gap between the top edge of any roll-on liner and the groin
- Brim sheath uses stretchy, slick nylon material to fit shape appropriately
- Includes gripper elastic to secure to outer socket



| Sock Size | Product No. |
|-----------|-----------------|
| Small | 2HBSAKSM |
| Medium | 2HBSAKMD |
| Large | 2HBSAKLG |

Silipos - Silosheath Original

- 2 ply, 2 mm gel
- Minimizes skin breakdown
- Reduces friction, shearing and abrasion on residual limb
- Cushions and protects
- Softens, soothes and improves the appearance of skin
- Product is washable and reusable

Double Cushion Silosheath™

- 2 ply, 3 mm gel, extra gel at distal end
- Minimizes skin breakdown
- Reduces friction, shearing and abrasion on residual limb
- Cushions and protects
- Softens, soothes and improves the appearance of skin
- Product is washable and reusable

Silosheath Extra Life™

- 3 ply, 2 mm gel, increased durability
- Minimizes skin breakdown
- Reduces friction, shearing and abrasion on residual limb
- Cushions and protects
- Softens, soothes and improves the appearance of skin
- Product is washable and reusable




| Size | Description | Gel Length | Approx. Length | Approx. Distal Circ. | Approx. Prox. Circ. | Product No. | | |
|------|--------------------|----------------|-----------------------|-------------------------|-------------------------|--------------|----------------|--------------|
| | | | | | | Original | Double Cushion | Extra Life |
| 1 | Child | 8" (20 cm) | 10"-13" (33 cm) | 6½"-9" (17-23 cm) | 6½"-10" (17-25 cm) | 92815 | 13415 | 13815 |
| 2 | Short/ Narrow | 10" (25 cm) | 15"-17" (38-43 cm) | 7"-10½" (18-27 cm) | 7"-12" (18-30.5 cm) | 92825 | 13425 | 13825 |
| 3 | Short/ Regular | 10" (25 cm) | 15"-17" (38-43 cm) | 8"-12½" (20-32 cm) | 8"-14" (20-35.5 cm) | 92835 | 13435 | 13835 |
| 4 | Medium/ Narrow | 13" (33 cm) | 17"-19" (43-48 cm) | 8"-12½" (20-32 cm) | 8"-14" (20-35.5 cm) | 92845 | 13445 | 13845 |
| 5 | Medium/ Regular | 10" (25 cm) | 17"-19" (43-48 cm) | 9"-14" (23-35.5 cm) | 9"-16" (23-40.5 cm) | 92855 | 13455 | 13855 |
| 6 | Large/ Narrow | 13" (33 cm) | 18"-21" (46-53 cm) | 9"-14" (23-35.5 cm) | 9"-16" (23-40.5 cm) | 92865 | 13465 | 13865 |
| 7 | Large/ Regular | 13" (33 cm) | 18"-21" (46-53 cm) | 11"-16" (28-41 cm) | 11"-18" (28-46 cm) | 92875 | 13475 | 13875 |
| 8 | Large/ Wide | 13" (33 cm) | 18"-21" (46-53 cm) | 12"-19" (30.5-48 cm) | 12"-21" (30.5-53 cm) | 92885 | 13485 | 13885 |

Truform - Sheath

- Effectively minimized friction at coupling
- 100% stretch nylon yarn, reinforced distal end, turned welt construction at top
- Sold 3 pair/box
- BK proximal circ: 12"-20", distal circ, 7"-12"
- AK proximal circ: 14"-18", distal circ, 9"-14"

| Size | Length | Product No. | |
|----------|---------|-----------------|-----------------|
| | | White | Dusk |
| BK Small | 10"-14" | 0830FBKS | 0830WBKS |
| BK Med | 16"-20" | 0830FBKM | 0830WBKM |
| BK Large | 22"-26" | 0830FBKL | 0830WBKL |
| AK Med | 14"-18" | 0830FAKM | 0830WAKM |

**Bell-Horn - Two-Way Stretch Seamless Elastic**

- Mid weight, nylon covered rubber

| Size | Fits Top | Length | Material | Color | Product No. |
|--------|----------|--------|----------|-------|--------------|
| Small | 16"-21" | 12" | Nylon | Beige | 185 S |
| Medium | 18"-23" | 14" | Nylon | Beige | 185 M |
| Large | 20"-25" | 16" | Nylon | Beige | 185 L |

**Comfort Products - BK Shriner**


- Tapered compression from the closed distal end to the open proximal end from 25 mm Hg to 20 mm Hg
- Pocket type toe with special side stitching for additional comfort
- Cotton comfort grip top makes donning easier, especially helpful for geriatric patients



| Length/ Size | Color Top Seam | Circumference | | | Product No. | | |
|-----------------|-------------------|----------------|-----------------|--------------|---------------|---------------|--|
| | | Bottom | Top | 9" | 12" | 14" | |
| Small | Gold | 6"-7 1/2" | 10"-12 1/4" | SS9S | SS12S | SS14S | |
| Medium | Red | 7 1/2"-8 1/2" | 12 1/4"-14 1/2" | SS9M | SS12M | SS14M | |
| Large | Blue | 8 1/2"-9 3/4" | 14 1/2"-17 1/4" | SS9L | SS12L | SS14L | |
| X-Large | Green | 9 3/4"-11 1/2" | 17 1/4"-20 1/2" | SS9XL | SS12XL | SS14XL | |

X-Silver BK Shriner

- X-Static Silver provides all natural antimicrobial protection
- Tapered compression from the closed end to the open proximal end
- Easy grip top eases donning
- New gel band at the top of shriner



| Length/ Size | Color Top Seam | Circumference | | | Product No. | | |
|-----------------|-------------------|---------------|-----------|---------|-------------|----------|--|
| | | Bottom | Top | 9" | 12" | 14" | |
| Small | Gold | 6"-7½" | 10"-12½" | SSX9S | SSX12S | SSX14S | |
| Medium | Red | 7½"-8½" | 12¼"-14½" | SSX9M | SSX12M | SSX14M | |
| Large | Blue | 8½"-9¾" | 14½"-17½" | SSX9L | SSX12L | SSX14L | |
| X-Large | Green | 9¾"-11½" | 17¼"-20½" | SSX9XL | SSX12XL | SSX14XL | |
| XX-Large | Black | 11½"-13½" | 20½"-24½" | SSX9XXL | SSX12XXL | SSX14XXL | |

X-Silver AK Shriner

- Toe: Pocket type with special stitching minimizes seam thickness and possible irritation
- Nylon covered body plus fine gauge fabric yields a softer feel
- Wide one-way stretch elastic top belt has adjustable Velcro closure for comfort and maximum support
- Waist measurement 34"; Contact SPS to special order 38" or 42" waist measurements



| Length/ Size | Color Top Seam | Circumference | | | Product No. | | |
|-----------------|-------------------|---------------|------|----------|-------------|-----------|--|
| | | Bottom | Top | 8" | 10" | 12" | |
| Small | Gold | 13½" | 15¾" | SSX8SAK | SSX10SAK | SSX12SAK | |
| Medium | Red | 14¾" | 17½" | SSX8MAK | SSX10MAK | SSX12MAK | |
| Med/Large | Black | 16" | 19" | SSX8MLAK | SSX10MLAK | SSX12MLAK | |
| Large | Blue | 17½" | 20½" | SSX8LAK | SSX10LAK | SSX12LAK | |
| X-Large | Green | 19½" | 22" | SSX8XLAK | SSX10XLAK | SSX12XLAK | |

Juzo - Dynamic (Varin) Soft-In Shriner

- Soft, smooth fabric for sensitive skin
- Maintains shape of stump
- Stay-put band has silicone border
- Flat, non-irritating distal seam
- Available in BK, AK with hip attachment and AK with body part
- To Order, refer to Sizing Charts on next page:

Amputation Level: CF (BK); FG (AK)

Lengths available: Small (S), Medium (M) or Large (L); refer to chart below for sizing parameters

Hip Attachment options: HA to include attachment or leave blank to omit

Circumference sizes: I, II, III, IV, V or VI per sizing chart on next page



| Description | Product No. |
|----------------------------------|-------------|
| Compression Class 1 (20-30 mmHg) | 3511???? |
| Compression Class 1 (30-40 mmHg) | 3512???? |

Silver Varin Soft-In

- Knitted with 17% X-Static silver threads
- Anti-microbial, anti-odor, cool wearing and hypoallergenic
- May reduce phantom pain
- Amputation Level: CF (BK); FG (AK)
- Lengths available: Small (S), Medium (M) or Large (L); refer to chart below for sizing parameters
- Hip Attachment options: HA to include attachment or blank to omit
- Circumference sizes: I, II, III, IV, V or VI per sizing chart below

| Description | Product No. |
|-----------------------------------|-----------------|
| Compression Class 1 (20-30 mm Hg) | 3561???? |
| Compression Class 1 (30-40 mm Hg) | 3562???? |

| Level of Body | Short (S) | Medium (M) | Long (L) | Extra Long (XL) |
|-----------------|--------------|---------------|-------------|--------------------|
| CF (Below Knee) | 12" 30 cm | 15" 38 cm | 18" 46 cm | — |
| FG (Above Knee) | 8" 20 cm | 10" 25 cm | 12" 30 cm | 14" 35 cm |

**Ordering Example**

Model **3511**
 Amputation Level **FG**
 Length **L**
 Hip Attachment **HA**
 Circumference **V**

3511FGLHA V

*Silicone border option is available

20-30 and 30-40 mmHg
Size Chart for Juzo® Above Knee Shriner

| | I | II | III | IV | V | VI |
|----------|------------------|------------------|------------------|------------------|------------------|------------------|
| G | 44 cm 17 1/4" | 48 cm 19" | 52 cm 20 1/2" | 56 cm 22" | 60 cm 23 3/4" | 64 cm 25 1/4" |
| F | 41 cm 16 1/4" | 44 cm 17 1/4" | 47 cm 18 1/2" | 50 cm 19 3/4" | 53 cm 20 3/4" | 56 cm 22" |
| E | 31 cm 12 1/4" | 34 cm 13 1/2" | 37 cm 14 1/2" | 40 cm 15 3/4" | 43 cm 17" | 46 cm 18" |

AK Lengths
 Extra Long - 14" / 35 cm
 Long - 12" / 30 cm
 Medium - 10" / 25 cm
 Short - 8" / 20 cm

20-30 and 30-40 mmHg
Size Chart for Juzo® Below Knee Shriner


| | 0 | I | II | III | IV | V | VI |
|----------|------------------|------------------|------------------|------------------|------------------|------------------|------------------|
| F | 38 cm 15 1/4" | 41 cm 16 1/4" | 44 cm 17 1/4" | 47 cm 18 1/2" | 50 cm 19 3/4" | 53 cm 20 3/4" | 56 cm 22" |
| E | 28 cm 11 1/4" | 31 cm 12 1/4" | 34 cm 13 1/2" | 37 cm 14 1/2" | 40 cm 15 3/4" | 43 cm 17" | 46 cm 18" |
| C | 26 cm 10 1/2" | 29 cm 11 1/2" | 32 cm 12 1/2" | 35 cm 13 3/4" | 38 cm 15" | 41 cm 16 1/4" | 44 cm 17 1/2" |

BK Lengths
 Long - 18" / 46 cm
 Medium - 15" / 38 cm
 Short - 12" / 30 cm


Knit-Rite - Compressogrip Elastic Shrinkers

- Tubular elastic compression material for control or prevention of edema following amputation surgery
- Speeds fabrication, fitting and delivery of a prosthesis by uncomplicating application of an elastic bandage
- Approximate single layer pressure at 50% stretch is 20-25 mmHg
- Doubled layers will double compression at both 50% and 100% stretch
- Applications of more than two layers should be monitored as compression may exceed limits of desired therapeutic pressure
- 92% polyester, 8% spandex; latex free; washable


| Size | Length | | Fits Circumference | | Product No. | |
|------|--------|----|--------------------|-------|-----------------|-----------------|
| | in. | cm | in. | cm | White | Beige |
| 2 | 18 | 45 | 5-7 | 13-18 | 2H0S0218 | 2H0B0218 |
| 2 | 24 | 60 | 5-7 | 13-18 | 2H0S0224 | 2H0B0224 |
| 2.5 | 18 | 45 | 7-10 | 18-25 | 2H0S0S18 | 2H0B0S18 |
| 2.5 | 24 | 60 | 7-10 | 18-25 | 2H0S0S24 | 2H0B0S24 |
| 3 | 18 | 45 | 9-12 | 23-30 | 2H0S0318 | 2H0B0318 |
| 3 | 24 | 60 | 9-12 | 23-30 | 2H0S0324 | 2H0B0324 |
| 3 | 30 | 75 | 9-12 | 23-30 | 2H0S0330 | 2H0B0330 |
| 3.5 | 18 | 45 | 10-13 | 25-33 | 2H0SS318 | 2H0BS318 |
| 3.5 | 24 | 60 | 10-13 | 25-33 | 2H0SS324 | 2H0BS324 |
| 4 | 18 | 45 | 12-16 | 30-40 | 2H0S0418 | 2H0B0418 |
| 4 | 24 | 60 | 12-16 | 30-40 | 2H0S0424 | 2H0B0424 |
| 4 | 30 | 75 | 12-16 | 30-40 | 2H0S0430 | 2H0B0430 |
| 5 | 18 | 45 | 15-20 | 38-50 | 2H0S0518 | 2H0B0518 |
| 5 | 24 | 60 | 15-20 | 38-50 | 2H0S0524 | 2H0B0524 |
| 5 | 30 | 75 | 15-20 | 38-50 | 2H0S0530 | 2H0B0530 |
| 6 | 18 | 45 | 18-24 | 45-60 | 2H0S0618 | 2H0B0618 |
| 6 | 24 | 60 | 18-24 | 45-60 | 2H0S0624 | 2H0B0624 |
| 6 | 30 | 75 | 18-24 | 45-60 | 2H0S0630 | 2H0B0630 |



APPLY THE FIRST LAYER



POSITION THE RING SNUGLY



APPLY THE SECOND LAYER TO A POINT OF DESIRED MAXIMUM COMPRESSION.

To Select Size

- Match "Fits Circumference" to:
- B/K circumference 2" below medial tibia tubercle
 - A/K circumference 2" from distal end
 - Allow for Swelling

Nylon Bone Ring

- Use with compressogrip tubular compression stockinette
- Pack of 100

| Bone Ring Size | Shrinker Size | Product No. |
|----------------|---------------|-----------------|
| 5/8" | 2", 2 1/2" | 5B705505 |
| 3/4" | 2", 2 1/2" | 5B705506 |
| 7/8" | 3", 4" | 5B705507 |
| 1" | 5" | 5B705510 |
| 1 1/8" | 6" | 5B705511 |

Compressogrip B/K Shrinker with X-Static®


- X-static silver fiber inhibits bacteria for more hygienic post-op environment
- Patients report relief from phantom pain
- Compressogrip shrinkers are effective, easy to apply; practical and affordable to inventory for immediate delivery
- Color: Gray; latex free, washable

| Size | Fits Circumference | Length | Product No. |
|------|--------------------|-------------|-----------------|
| #3 | 9"-12" (23-30 cm) | 24" (60 cm) | 2H0X0324 |
| #4 | 12"-16" (30-40 cm) | 24" (60 cm) | 2H0X0424 |
| #5 | 15"-20" (38-50 cm) | 24" (60 cm) | 2H0X0524 |
| #6 | 18"-24" (45-60 cm) | 24" (60 cm) | 2H0X0624 |
| #3.5 | 10"-13" (25-33 cm) | 30" (75 cm) | 2H0XS330 |
| #4 | 12"-16" (30-40 cm) | 30" (75 cm) | 2H0X0430 |
| #5 | 15"-20" (38-50 cm) | 30" (75 cm) | 2H0X0530 |
| #6 | 18"-24" (45-60 cm) | 30" (75 cm) | 2H0X0630 |

**Tapered Compressogrip B/K Elastic Shrinker**

- Additional width at the top accommodates fuller thighs and reduces rolling and bunching behind the knee
- Used for control or prevention of edema, for limb volumetric maintenance, to minimize scarring and to give compressive support for pain relief
- Latex free and washable
- White/Beige - 92% polyester, 8% spandex; Gray - 87.4% polyester, 4.6% X-STATIC®, 8% spandex
- When determining size, allow for swelling

| Size | Fits Circumference | | | Product No. | | |
|------|--------------------|-------|-------------|-----------------|-----------------|-----------------|
| | in. | cm. | Length | White | Beige | Gray |
| #4 | 12"-16" | 30-40 | 24" (60 cm) | 2H0ST424 | 2H0BT424 | 2H0XT424 |
| #4 | 12"-16" | 30-40 | 30" (75 cm) | 2H0ST430 | 2H0BT430 | 2H0XT430 |
| #5 | 15"-20" | 38-50 | 24" (60 cm) | 2H0ST524 | 2H0BT524 | 2H0XT524 |
| #5 | 15"-20" | 38-50 | 30" (75 cm) | 2H0ST530 | 2H0BT530 | 2H0XT530 |
| #6 | 18"-24" | 45-60 | 24" (60 cm) | 2H0ST624 | 2H0BT624 | 2H0XT624 |
| #6 | 18"-24" | 45-60 | 30" (75 cm) | 2H0ST630 | 2H0BT630 | 2H0XT630 |




| Flat Width | Top |
|--------------|---|
| 4" (10.2 cm) | 5" flat width; tapers last 4" of length |
| 5" (12.7 cm) | 6.5" flat width; tapers last 5" of length |
| 6" (15.2 cm) | 7.5" flat width; tapers last 6" of length |

To Select Size

1. Match "Fits Circumference" to B/K circumference 2" below medial tibia tubercle
2. Compressogrip fabric gives 20-25 mm Hg compression at 50% stretch; at 100% stretch the compression is almost doubled
3. Double layers will almost double compression at both 50% and 100% stretch
4. Application of more than 2 layers may exceed the limits of desired therapeutic pressure

B/K Shrinker

- Made from Rayon/Lycra core yarn
- Softer, more comfortable, easily donned and may be machine washed and dried, warm temp. no bleach
- Medium compression (gray top) exerts 25-30 mm Hg pressure at 50% stretch for tender stumps, nighttime wear and limb size measurement
- Heavy compression (green top) exerts 55-60 mm Hg pressure at 50% stretch for reduction of edema and shaping of limb
- Replace ?? with length 10", 12", 14"
- Double taper has different measurements
- Corespun yarn results in softer fabric, better fit and permits stretch over the circumference range



Double Taper Illustration

| Pressure Weight | Width | Flat Measurement | | | | Fits Circumference | | | | Product No. | |
|-----------------|--------|------------------|-----|--------------|--------|--------------------|---------|--------------|---------|---------------|-------------------|
| | | Regular Taper | | Double Taper | | Regular Taper | | Double Taper | | Regular Taper | Double Taper |
| | | Toe | Top | Toe | Mid | Toe | Top | Toe | Top | | |
| Heavy | Narrow | 3" | 5" | 3" | 3 1/2" | 5" | 9"-11" | 13"-14" | 9"-11" | 15"-17" | 2HPHRN?? 2HPHDN?? |
| Heavy | Medium | 4" | 6" | 4" | 4 1/2" | 7" | 12"-14" | 15"-17" | 12"-14" | 18"-20" | 2HPHRM?? 2HPHDM?? |
| Heavy | Wide | 5" | 7" | 5" | 5 1/2" | 8" | 15"-17" | 18"-20" | 15"-17" | 21"-24" | 2HPHRW?? 2PHHDW?? |
| Medium | Narrow | 3" | 5" | 3" | 3 1/2" | 5" | 9"-11" | 13"-14" | 9"-11" | 15"-17" | 2HPMRN?? 2HPMDN?? |
| Medium | Medium | 4" | 6" | 4" | 4 1/2" | 7" | 12"-14" | 15"-17" | 12"-14" | 18"-20" | 2HPMRM?? 2HPMDM?? |
| Medium | Wide | 5" | 7" | 5" | 5 1/2" | 8" | 15"-17" | 18"-20" | 15"-17" | 21"-24" | 2HPMRW?? 2HPMDW?? |

Compressogrip A/K Shrinker

- Clinically effective
- Universal for left or right application
- Multi-directional stretch for improved fit and comfort
- Latex free; easy to use and apply

| Size | Fits Circ. | Flat Width | Product No. |
|----------------------|--------------------|------------|-------------|
| X-Small | 12"-17" (30-43 cm) | 5" | 2H-0S05-XS |
| Small | 15"-20" (38-50 cm) | 6" | 2H-0S06-SM |
| Medium | 18"-24" (45-60 cm) | 7 1/2" | 2H-0SS7-MD |
| Large | 24"-32" (60-80 cm) | 9" | 2H-0S09-LG |
| Medium (with gusset) | 18"-24" (45-60 cm) | 7 1/2" | 2HGSS7MD |
| Large (with gusset) | 24"-32" (60-80 cm) | 9" | 2HGS09LG |




Compressogrip A/K Shriner with X-Static

- X-static silver fiber inhibits bacteria for more hygienic post-op environment
- Clinically effective; latex free, easy to use and apply
- Universally effective for left and right application
- Multi-directional stretch for improved fit and comfort
- Latex free; easy to use and apply

| Size | Fits Circ. | Flat Width | Product No. |
|---------|--------------------|------------|-------------------|
| X-Small | 12"-17" (30-43 cm) | 5" | 2H-0X05-XS |
| Small | 15"-20" (38-50 cm) | 6" | 2H-0X06-SM |
| Medium | 18"-24" (45-60 cm) | 7½" | 2H-0XS7-MD |
| Large | 24"-32" (60-80 cm) | 9" | 2H-0X09-LG |

**AK Shriner**


- Easy to adjust, hook-and-pile, elongated side opening waist band, adaptable for left or right amputation
- Optional auxiliary support for proximal edge at groin can be attached at any point on waist band or easily detached and removed (item 6W-7322-3)
- Tapered design for gradient pressure; heavy (60-65 mm Hg) or medium (30-35 mm Hg) pressure available
- Superior 3-dimensional support and shaping of distal end; garment length from groin to distal end should be one to two inches shorter than actual measurement so proper snugging of distal end will not result in roll at groin
- Replace ?? with length 7", 10", 13" (groin to distal)



| Flat Measurement | | Fits Circumference | | Pressure Product No. | |
|------------------|-----|--------------------|-----------|----------------------|----------------|
| Toe | Top | Toe | Top | Heavy | Medium |
| 4 | 6 | 12"-14½" | 15"-17" | 2PH64?? | 2PM64?? |
| 4 | 7 | 12"-14½" | 17½"-20½" | 2PH74?? | 2PM74?? |
| 4 | 8 | 12"-14½" | 21"-23" | 2PH84?? | 2PM84?? |
| 5 | 7 | 15"-17½" | 17½"-20½" | 2PH75?? | 2PM75?? |
| 5 | 8 | 15"-17½" | 21"-23" | 2PH85?? | 2PM85?? |
| 5 | 9 | 15"-17½" | 23½"-25½" | 2PH95?? | 2PM95?? |
| 6 | 8 | 18"-22" | 21"-23" | 2PH86?? | 2PM86?? |
| 6 | 9 | 18"-22" | 23½"-25½" | 2PH96?? | 2PM96?? |
| 6 | 10 | 18"-22" | 26"-29" | 2PH06?? | 2PM06?? |

Compressogrip Tubular Elastic Bandage

- Edema reduction and control
- Casting compression and soft tissue control
- More consistent compression than elastic wrap
- Note: Single layer supportive compression (20-25 mm Hg), double layer provides therapeutic compression (35-50 mm Hg) at 50% stretch
- All sizes provided in 11 yard (10 meter) rolls
- 92% Polyester, 8% Spandex, Latex free
- Machine wash and dry, warm, no bleach
- Fabric will relax 20% in length after first washing, autoclaveable



| Size | Product No. | | | | |
|------|----------------------------|-----------------|-----------------|-----------------|-------------------------|
| | Circumference in. cm | White Boxed | Beige Boxed | Black Boxed | Natural Plastic Pkg. |
| 2 | 5-7 12.5-17.5 | 6S762802 | 6S763102 | — | 6S76CT02 |
| 2½ | 7-10 18-25 | 6S7628S2 | 6S7631S2 | — | 6S76CTS2 |
| 2¾ | 8-11 20-28 | 6S7628T2 | 6S7631T2 | — | 6S76CTT2 |
| 3 | 9-12 23-30 | 6S762803 | 6S763103 | 6S763303 | 6S76CT03 |
| 3½ | 10-13 25-33 | 6S7628S3 | 6S7631S3 | — | 6S76CTS3 |
| 4 | 12-16 30.5-40.5 | 6S762804 | 6S763104 | 6S763304 | 6S76CT04 |
| 5 | 15-20 38-51 | 6S762805 | 6S763105 | 6S763305 | 6S76CT05 |
| 6 | 18-24 46-61 | 6S762806 | 6S763106 | 6S763306 | 6S76CT06 |
| 8 | 24-32 61-81 | 6S762808 | 6S763108 | 6S763308 | 6S76CT08 |
| 10 | 30-40 76-102 | 6S762810 | 6S763110 | — | 6S76CT10 |
| 12 | 36-48 91-122 | — | — | — | 6S76CT12 |
| 14 | 42-56 107-142 | — | — | — | — |

Post-Op Dressing Retention Sock

- Utilizes a very soft, compressive fit, with a belt suspension to cover and retain postoperative dressings for transfemoral amputees
- Very soft Lycra and CoolMax body for low compression
- Universal right or left belt system to maintain secure suspension

| Size | Fits | Length | Product No. |
|-----------------|---------|--------|-----------------|
| 4" Post Op Sock | 10"-12" | 28" | 1RP10428 |
| 6" Post Op Sock | 14"-18" | 32" | 1RP10632 |
| 8" Post Op Sock | 18"-24" | 32" | 1RP10832 |



Rx Textiles - Tensitube® Blue

- A latex-free tubular elastic shrinker that helps reduce and control hypertrophic scarring and edema on residual limbs
- A ring eases application
- Shrinker packs are primarily used for post-op residual limb shrinker
- Easily applied and offers controlled compression which accelerates the prosthetic fit
- Perfect for casting and remains effective after repeated washings



| Tensitube Shrinker Packs | | | |
|---------------------------------|--------------------|----------------|----------------|
| Size | Product No. | | |
| | White | Beige | Silver |
| 2" x 18" | 0100100 | 0140100 | 0150100 |
| 2" x 24" | 0100200 | 0140200 | 0150200 |
| 3" x 18" | 0100700 | 0140700 | 0150700 |
| 3" x 24" | 0100800 | 0140800 | 0150800 |
| 3" x 30" | 0100900 | 0140900 | 0150900 |
| 4" x 18" | 0101000 | 0141000 | 0151000 |
| 4" x 24" | 0101100 | 0141100 | 0151100 |
| 4" x 30" | 0101200 | 0141200 | 0151200 |
| 5" x 18" | 0101300 | 0141300 | 0151300 |
| 5" x 24" | 0101400 | 0141400 | 0151400 |
| 5" x 30" | 0101500 | 0141500 | 0151500 |
| 6" x 18" | 0101600 | 0141600 | 0151600 |
| 6" x 24" | 0101700 | 0141700 | 0151700 |
| 6" x 30" | 0101800 | 0141800 | 0151800 |

Supreme Medical - Elastic Tubular Support Bandages

- 11 yard rolls
- Sold by the box (1 roll per box)
- 12 boxes (rolls) per case

| Description | Size | Product No. |
|--------------------|-------------|--------------------|
| D | 3" | MOL1437 |
| E | 3½" | MOL1434 |
| F | 4" | MOL1438 |
| G | 4½" | MOL1439 |
| Large | 15-17" | MOL1442 |



Truform - Medium Weight BK Stump Shrinker

- Two-way stretch medium weight beige spandex construction
- Lengths: 8", 10", 12", and 14"
- Determine length, measure from distal end to knee joint, plus 6"
- Replace ?? with length in inches



| Size | Knee Measurement | Product No. |
|----------|--|---------------------------|
| X-Small | fits to 12 ¹ / ₄ " | 31 cm 0839BK??-XS |
| Small | fits to 13 ¹ / ₂ " | 34.5 cm 0839BK??-S |
| Medium | fits to 14 ¹ / ₂ " | 37 cm 0839BK??-M |
| Large | fits to 15 ³ / ₄ " | 40 cm 0839BK??-L |
| X-Large | fits to 17" | 43 cm 0839BK??-XL |
| XX-Large | fits to 18" | 46 cm 0839BK??-XXL |

Lightweight BK Stump Shrinker

- Two-way stretch nylon covered with real rubber core beige elastic
- Length: 14"
- Determine length, measure from distal end to knee joint, plus 6"



| Size | Knee Measurement | Product No. |
|--------|------------------|----------------------------|
| Small | fits 8-12" | 20.5-30.5 cm 0832-S |
| Medium | fits 12-16" | 30.5-40.5 cm 0832-M |
| Large | fits 16-20" | 40.5-51 cm 0832-L |

BK Stump Shrinker

- One-way stretch surgical elastic with padding at distal end
- Lengths: 6", 8", 10", 12", 14"
- Determine length, measure from distal end to knee joint, plus 6"
- Replace ?? with length in inches



| Size | Calf Circumference | Product No. |
|---------|--------------------|-----------------------------|
| X-Small | 9-11" | 23-28 cm 0833BK??-XS |
| Small | 11-13" | 28-33 cm 0833BK??-S |
| Medium | 13-15" | 33-38 cm 0833BK??-M |
| Large | 15-17" | 38-43 cm 0833BK??-L |
| X-Large | 17-19" | 43-48 cm 0833BK??-XL |

AK Stump Shrinker

- Convenient one-piece construction
- Surgical weight elastic, one-way stretch, contoured at groin with padded draw-string bottom
- Universally adjustable elastic waistband with hook-and-pile fastening
- Lengths: 6", 8", 10", 12" and 14"
- Replace ?? with length



| Size | Thigh Circ. | Product No. |
|---------|-------------|--------------------|
| Small | 13-16" | 0831AK??-S |
| Medium | 16-19" | 0831AK??-M |
| Large | 19-22" | 0831AK??-L |
| X-Large | 22-25" | 0831AK??-XL |

AK Stump Shrinker with Expansion Gusset

- Single piece one-way stretch surgical elastic construction
- Universally adjustable elastic waistband with hook-and-pile fastening
- Slightly tapered expansion gusset, 4 1/2" wide by 3 1/4" deep, at the inner proximal edge where the stump meets the trunk of the body
- Gusset prevents irritation over the inner thigh near the groin area
- Lengths: 6", 8", 10", 12" and 14"
- Replace ?? with length



| Size | Thigh Circ. | Product No. |
|---------|-------------|--------------------|
| Small | 13-16" | 0837AK??-S |
| Medium | 16-19" | 0837AK??-M |
| Large | 19-22" | 0837AK??-L |
| X-Large | 22-25" | 0837AK??-XL |

Medium Weight AK Stump Shrinker

- Single piece two-way stretch medium weight beige spandex elastic with hook-and-pile fastening
- Wide contour cut around groin for comfort
- Easy to apply, effectively delivers 35 mm Hg compression
- Lengths: 8", 10", 12" and 14"



| Size | Thigh Circ. | Product No. |
|----------|-------------|---------------------|
| Small | 13-16" | 0839AK??-S |
| Medium | 16-19" | 0839AK??-M |
| Large | 19-22" | 0839AK??-L |
| X-Large | 22-25" | 0839AK??-XL |
| XX-Large | 25-27" | 0839AK??-XXL |