


To Determine proper Aegis SoFlex liner sizing, two measurements are required:

1. Take circumference measurement 2" (5-6cm) proximal to distal end of residual limb.
2. Take length measurement from mid-patella tendon to distal end of residual limb.

For liner size and item #, see chart below.

Note: Aegis SoFlex liners come with full fabric cover, in either locking or cushion (non-locking) versions.



AEGIS SoFlex

1. Limb Cir.	2. Length	Locking	Cushion
17-19 cm	4"	SBK 14-4L	SBK 14-4C
	6"	SBK 14-6L	SBK 14-6C
	8"	SBK 14-8L	SBK 14-8C
19-21 cm	4"	SBK 16-4L	SBK 16-4C
	6"	SBK 16-6L	SBK 16-6C
	8"	SBK 16-8L	SBK 16-8C
21-23 cm	4"	SBK 18-4L	SBK 18-4C
	6"	SBK 18-6L	SBK 18-6C
	8"	SBK 18-8L	SBK 18-8C
23-25 cm	4"	SBK 20-4L	SBK 20-4C
	6"	SBK 20-6L	SBK 20-6C
	8"	SBK 20-8L	SBK 20-8C
	10"	SBK 20-10L	SBK 20-10C
25-27 cm	4"	SBK 22-4L	SBK 22-4C
	6"	SBK 22-6L	SBK 22-6C
	8"	SBK 22-8L	SBK 22-8C
	10"	SBK 22-10L	SBK 22-10C
27-29 cm	4"	SBK 24-4L	SBK 24-4C
	6"	SBK 24-6L	SBK 24-6C
	8"	SBK 24-8L	SBK 24-8C
	10"	SBK 24-10L	SBK 24-10C
30-32 cm	4"	SBK 26-4L	SBK 26-4C
	6"	SBK 26-6L	SBK 26-6C
	8"	SBK 26-8L	SBK 26-8C
	10"	SBK 26-10L	SBK 26-10C
32-34 cm	4"	SBK 28-4L	SBK 28-4C
	6"	SBK 28-6L	SBK 28-6C
	8"	SBK 28-8L	SBK 28-8C
	10"	SBK 28-10L	SBK 28-10C
35-37 cm	4"	SBK 30-4L	SBK 30-4C
	6"	SBK 30-6L	SBK 30-6C
	8"	SBK 30-8L	SBK 30-8C
	10"	SBK 30-10L	SBK 30-10C