

# SmartKnit AFO and KAFO Socks

## SMARTKNIT® SEAMLESS AFO AND KAFO INTERFACE SOCKS

Each patented SmartKnit® AFO and KAFO sock is made the same way that a caterpillar spins its cocoon—starting at the toe and circling to the top, resulting in a seam-free, wrinkle-free sock that fits like a second layer of skin, provides exceptional softness, and wicks perspiration. Socks feature a “heel-less” design to ensure a perfect fit every time. Each corespun sock is knit with CoolMax®, Polyester, or X-STATIC®—The Silver Fiber® with Lycra® spandex for superior moisture wicking and anti-microbial properties to inhibit odor in the sock. Machine wash and dry, medium temperature, no bleach. Patent No. 6,158,253.



|                                                  | SIZE          | COLOR  | WIDTH / LENGTH           | SOLD |
|--------------------------------------------------|---------------|--------|--------------------------|------|
| <b>SmartKnit® Child Size AFO Interface Socks</b> |               |        |                          |      |
| <b>1SP1-AF-NS</b>                                | Infant Small  | White  | 1.5" x 7.5" (4 x 19cm)   | Pair |
| <b>1SP1-AF-NF</b>                                | Infant        | White  | 2" x 9" (5 x 23cm)       | Pair |
| <b>1SP1-AF-CS</b>                                | Child Small   | White  | 2.25" x 11" (5.5 x 28cm) | Pair |
| <b>1SP1-AF-CR</b>                                | Child Regular | White  | 2.50" x 14" (6.5 x 36cm) | Pair |
| <b>1SP1-AF-CX</b>                                | Child X-Long  | White  | 2.75" x 20" (7 x 51cm)   | Pair |
| <b>1S21-AF-NF</b>                                | Infant        | Pink   | 2" x 9" (5 x 23cm)       | Pair |
| <b>1S21-AF-CS</b>                                | Child Small   | Pink   | 2.25" x 11" (5.5 x 28cm) | Pair |
| <b>1S21-AF-CR</b>                                | Child Regular | Pink   | 2.50" x 14" (6.5 x 36cm) | Pair |
| <b>1S31-AF-NF</b>                                | Infant        | Purple | 2" x 9" (5 x 23cm)       | Pair |
| <b>1S31-AF-CS</b>                                | Child Small   | Purple | 2.25" x 11" (5.5 x 28cm) | Pair |
| <b>1S31-AF-CR</b>                                | Child Regular | Purple | 2.50" x 14" (6.5 x 36cm) | Pair |

|                                                                                     |               |          |                          |      |
|-------------------------------------------------------------------------------------|---------------|----------|--------------------------|------|
| <b>SmartKnit® Child Size AFO Interface Socks with X-STATIC® – The Silver Fiber®</b> |               |          |                          |      |
| <b>1NX1-AF-NF</b>                                                                   | Infant        | Navy     | 2" x 9" (5 x 23cm)       | Pair |
| <b>1NX1-AF-CS</b>                                                                   | Child Small   | Navy     | 2.25" x 11" (5.5 x 28cm) | Pair |
| <b>1NX1-AF-CR</b>                                                                   | Child Regular | Navy     | 2.50" x 14" (6.5 x 36cm) | Pair |
| <b>1SX1-AF-NF</b>                                                                   | Infant        | Grey     | 2" x 9" (5 x 23cm)       | Pair |
| <b>1SX1-AF-CS</b>                                                                   | Child Small   | Grey     | 2.25" x 11" (5.5 x 28cm) | Pair |
| <b>1SX1-AF-CR</b>                                                                   | Child Regular | Grey     | 2.50" x 14" (6.5 x 36cm) | Pair |
| <b>1SX1-AF-CX</b>                                                                   | Child X-Long  | Grey     | 2.75" x 20" (7 x 51cm)   | Pair |
| <b>1S51-AF-CS</b>                                                                   | Child Small   | Charcoal | 2.25" x 11" (5.5 x 28cm) | Pair |
| <b>1S51-AF-CR</b>                                                                   | Child Regular | Charcoal | 2.50" x 14" (6.5 x 36cm) | Pair |
| <b>1BX1-AF-NF</b>                                                                   | Infant        | Black    | 2" x 9" (5 x 23cm)       | Pair |
| <b>1BX1-AF-CS</b>                                                                   | Child Small   | Black    | 2.25" x 11" (5.5 x 28cm) | Pair |
| <b>1BX1-AF-CR</b>                                                                   | Child Regular | Black    | 2.50" x 14" (6.5 x 36cm) | Pair |

|                                                   |            |        |                          |      |
|---------------------------------------------------|------------|--------|--------------------------|------|
| <b>SmartKnit® Child Size KAFO Interface Socks</b> |            |        |                          |      |
| <b>1SP1-AF-CL</b>                                 | Child Long | White  | 2.50" x 18" (6.5 x 46cm) | Each |
| <b>1S21-AF-CL</b>                                 | Child Long | Pink   | 2.50" x 18" (6.5 x 46cm) | Each |
| <b>1S31-AF-CL</b>                                 | Child Long | Purple | 2.50" x 18" (6.5 x 46cm) | Each |

|                                                                                      |            |       |                          |      |
|--------------------------------------------------------------------------------------|------------|-------|--------------------------|------|
| <b>SmartKnit® Child Size KAFO Interface Socks with X-STATIC® – The Silver Fiber®</b> |            |       |                          |      |
| <b>1NX1-AF-CL</b>                                                                    | Child Long | Navy  | 2.50" x 18" (6.5 x 46cm) | Each |
| <b>1SX1-AF-CL</b>                                                                    | Child Long | Grey  | 2.50" x 18" (6.5 x 46cm) | Each |
| <b>1BX1-AF-CL</b>                                                                    | Child Long | Black | 2.50" x 18" (6.5 x 46cm) | Each |



KAFO's longer lengths provide full leg protection.



SmartKnit® AFO and KAFO Socks available in new, fun colors



No heel assures custom fit