

Code L1832/1833*

GenuTrain® S Pro

For improved lateral stability.

Indications

- mild to moderate instability
- moderate to severe knee OA
- post injury management

Pain Relief and Stabilization

- **Medical-strength compression** – utilizes graduated medical compression to reduce edema, accelerate pain reduction and enhance mobility.
- **Intermittent massage** – Anatomically-designed medical silicone inserts, provides permanent massage through soft tissue stimulation while active, further enhancing pain relief and promoting faster rehabilitation.
- **Enhanced proprioception** – GenuTrain S Pro stimulates proprioception, providing improved muscle and joint control.

Wearing Comfort

3D anatomically contoured knit, assures optimal support and stability across the knee joint delivering appropriate weight bearing stress relief.

Patient Compliance

Compliance is totally dependent on patient comfort throughout the healing process. Bauerfeind uses high-tech breathable, stretchable and moisture-dissipating textiles with an even softer knit for extra sensitive areas such as behind the knee.

Bauerfeind products are made in Germany by a team dedicated to quality, clinical effectiveness and detailed craftsmanship.

*Suggested codes are based on approved coding from PDAC as of 01/18. It is the responsibility of the provider to determine the appropriate billing code, as well as, whether the use of a product complies with medical necessity and other documentation requirements of the payer.

Now
available as
**CUSTOM-
MADE**



GenuTrain® S Pro



Knit

3D anatomically contoured fit providing medical-grade compression

Strap

adjustable tension straps provide enhanced knee stability and proper anatomical fit

Patella Pad

Provides targeted compression and pain relief

Hinge

easy and quick joint control adjustment

adjustment to 30°, 60° and 90° flexion, as well 20° extension


Sidebars

Mirror anatomic shape of the leg and can be adapted to patient needs (thermo-moldable)

Now available in **CUSTOM MADE SIZES**



also available in black



A 14 cm 5½ inch
B 12 cm 4¾ inch

Size	1	2	3	4	5	6	7	Order number
Circumfer. 4¾" below knee	11 - 12¼	12¼ - 13¾	13¾ - 14¾	14¾ - 15¾	15¾ - 17	17 - 18½	18½ - 19¼	110413500*0+0°
Circumfer. 12 cm below knee	28-31	31-34	34-37	37-40	40-43	43-46	46-49	
Circumfer. 5½" above knee	15 - 16⅞	16⅞ - 17¾	17¾ - 18½	18½ - 19¾	19¾ - 20⅞	20⅞ - 22	22 - 23¼	
Circumfer. 14 cm above knee	38-41	41-44	44-47	47-50	50-53	53-56	56-59	



*1 = nature, *7 = black, *8 = titan, +6 = right, +7 = left, ° = size

BAUERFEIND USA, INC. 1230 Peachtree St. NE • Suite 1000 • Atlanta, GA 30309
P 1 800 423 3405 • P (404) 201 7800 • F (404) 201 7839
E-mail info@bauerfeindusa.com